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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication
of this manual, or part of it, without the written agreement of
GIBERTINI ELETTRONICA Srl*

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WARRANTY
WARNING
DECLARATION OF CONFORMITY

INDEX

1. WARNINGS.....	18
2. TECHNICAL CHARACTERISTICS.....	19
3. LOCATION.....	19
4. INSTALLATION.....	20
5. OPERATING CONTROLS OF THE BALANCE.....	21
6. OPERATING INSTRUCTIONS.....	22
7. KEYBOARD OPERATION.....	22
8. TARE.....	22
9. MAIN MENU AND OPERATING FUNCTIONS.....	23
9.1. UNIT 2.....	24
9.1.1. STANDARD BALANCE FUNCTION.....	24
9.1.2. PIECE COUNTING FUNCTION.....	24
9.1.3. PERCENTAGE FUNCTION.....	25
9.2. SETUP.....	26
9.2.1. SERIAL.....	27
9.2.1.1. BAUD.....	27
9.2.1.2. BYTE FOR.....	27
9.2.1.3. PROTOC.....	27
9.2.1.4. TRANS OPT.....	27
9.2.2. CALIBR.....	27
9.2.2.1. CAL OFF.....	27
9.2.2.2. CAL ON.....	27
9.2.2.3. CAL AUTO.....	27
9.2.3. A-ZERO.....	28
9.2.3.1. ZERO OFF.....	28
9.2.3.2. ZERO ON.....	28
9.2.4. FILTER.....	28
9.2.4.1. SLO.....	28
9.2.4.2. FAST.....	28
9.2.5. SET DEF.....	28
9.2.5.1. STORE.....	28
9.2.5.2. END.....	28
9.3. CAL E.....	28
9.4. CAL INT.....	29
9.2.5.1. STORE.....	29
9.2.5.2. END.....	29
9.5. TIMED AUTOMATIC CALIBRATION.....	29
10. QUICK RESTORE OF DEFAULT PARAMETERS.....	29
11. SERIAL OUTPUT.....	30
11.1. SERIAL OUTPUT CONNECTION.....	30
11.2. DATA TRANSMISSION FORMAT.....	30
11.3. REMOTE CONTROL COMMANDS.....	31

NOTE

**Read the Instructions carefully before operating the Balance.
This is an high accuracy Instrument and must be handled with care.**

1. WARNINGS

The power supply pack of the Crystal balance complies with Electrical Protection Class II.

Before connecting the power, check that the operating voltage, indicated on the Power supply pack, matches the line voltage.

CRYSTAL:

MAY ONLY BE USED WITH THE POWER SUPPLY PACK SUPPLIED.

MUST ONLY BE OPERATED IN DRY, CLEAN AND VIBRATION-FREE ENVIRONMENT.

MUST NOT BE OPERATED IN AREAS WHERE THERE ARE EXPLOSION RISKS.

NOTE:

In case of improper setting of the balance parameters during setup, the balance may power-up in an incorrect configuration.

The user can restore the factory configuration to avoid the malfunction as follows:

- ◆ Disconnect the supply of the balance
- ◆ Hold the On/Off key pressed
- ◆ Reconnect the supply
- ◆ Release the On/Off key when the display shows 'Set Factory'

Now the factory default settings are restored.

ALL ELECTRONIC BALANCES, WHICH ARE SENSITIVE TO CHANGES OF GRAVITY, MUST BE CHECKED AND CALIBRATED IN THE PLACE WHERE THEY ARE TO BE USED (CE 2014/30/UE).

ALL OUR INSTRUMENTS MEET INTERNATIONAL EMC-EMI-RFI STANDARDS ACCORDING TO CE 2014/31/UE DIRECTIVE ON ELECTROMAGNETIC COMPATIBILITY.

2. TECHNICAL CHARACTERISTICS

	MODELS					
	CRYSTAL 100	CRYSTAL 100 CAL IN	CRYSTAL 200 SMI	CRYSTAL 200	CRYSTAL 300	CRYSTAL 500
Capacity	110 g		210 g		310 g	510 g
Resolution	0,1mg					1mg
Linearity	± 0,2 mg					± 2 mg
External Calibration Weight:	50 g or multiples		100g or multiples			
Stabilization time (fast/slow)	6 / 10 s					5 / 8 s
Operating temperature	15 ÷ 25 °C					15 ÷ 30 °C
Pan diameter	85 mm					110 mm
Operating Voltage	100 / 240 Vac -15/+10% (range 85 to 264 V) 50/60Hz					
Power consumption	200 mA					
Data output	RS232 I/O					
Total dimensions	215 x 350 x 310 mm (w x d x h)					
Weighing chamber dimensions	210 x 155 x 220 mm (w x d x h)					
Weight	7 Kg (power supply included)					
Internal Calibration	NO	YES	NO	YES	YES	YES

3. LOCATION

Choose an appropriate location taking into account these criteria:

Vibrations and air drafts in the ambient must be avoided.

The supporting stand must be rigid.

The mains line must be stable. It is not advisable to use the same electric line of the fluorescent lamps or other electric apparatus which can be switched on and off randomly.

The balance must be kept far from magnetic fields.

Environmental conditions with wide temperature changes must be avoided.

Ambient temperature above or below temperature in the table may affect the operation of the balance and its accuracy.

In very dry environments electrostatic charges may build up, especially if glass or porcelain objects are used and even more with plastic (plastic should be avoided in any case).

WARNING !

Both the balance pan and the sample to be weighed must be very clean and kept at the same temperature.

4. INSTALLATION

Instructions to remove the locking device of the balance

- ◆ Unpack the balance and all its accessories.
- ◆ Unscrew the knob on the locking bracket as indicated by the arrow (anti clockwise).
- ◆ Retain the locking bracket in case the balance has to be returned to the supplier. This bracket can be refastened to the balance by replacing the knob.


FIG. 1

Level the balance using the spirit level on the rear and rotating the levelling feet under the front part of the basement (Fig 3).

Insert the drip cup in its support, log the protection ring on top of the drip cup (NO for models 500), insert the pan.

Insert the power supply cable located on the rear panel (Fig 2).

Insert the power supply in the main line plug checking the supply voltage printed on the power supply case (100 / 240 Vac) - (Fig 2).

FIG. 2 - Rear panel details

5. OPERATING CONTROLS OF THE BALANCE


FIG. 3 – General view

- 1- Main Display
- 2- On/Off Pushbutton
- 3- Tare Pushbutton
- 4- Mode Pushbutton
- 5- Print Pushbutton
- 6- Weighing Pan
- 7- Levelling feet

6. OPERATING INSTRUCTIONS

Wait for about 3 hours after power-up before performing any high accuracy weighing operation.

It is advisable to leave the balance permanently on (connected to the AC line).

It is possible to switch off only the display with the ON/OFF key, this allows a quick recovery when the balance is switched on.

Crystal balances have ranges with capacity of 110 g, 210 g, 310 g, 510 g with up to 0.1 mg resolution.

Pressing the ON/OFF key, after the display test, the indication 0,0000 g. appears.

If the instrument has not been used for more than an hour, perform a weighing operation with at least 10 g before high precision measurements.

7. KEYBOARD OPERATION

The keyboard has four control keys: PRINT, MODE, TARE, ON-OFF.

Every key has a specific function, while keys MODE, TARE and ON-OFF also activate a second function.

If the keys are activated for less than 1 second the first function is activated, when pressed for more the second function will be active.

The table below shows the operation capabilities.

Normal use

KEY NAME	KEY ACTIVATED FOR LESS THAN 1 SECOND	KEY ACTIVATED FOR MORE THAN 1 SECOND
On/Off	Display illumination on/off	Balance on/off
Mode	See additional functions	See additional functions
Tare	Tare function	Tare function
Print	Print function	No function

Main & configuration menu

KEY NAME	KEY ACTIVATED FOR LESS THAN 1 SECOND	KEY ACTIVATED FOR MORE THAN 1 SECOND
On/Off	Confirm selection	No function
Mode	Scroll available items	No function
Tare	Exit current selection	No function
Print	No function	No function

8. TARE

After a container is placed on the pan, push the TARE key to zero the reading.

The balance has a tare capacity up to its total measuring range.

Place the sample to be weighed in the container.

The display shows the sample net weight.

It is possible to use the tare key again and again after adding several samples, up to the maximum capacity of the balance, indicated by a bargraph on the left side of the LCD display.

9. MAIN MENU AND OPERATING FUNCTIONS

The balance can be operated in the following configurations:

Balance with gram, carat, ounce, pound measuring unit
 Piece counting balance
 Percentage computing balance

The selection is performed as follows:


Hold MODE key pressed until the display shows "UNIT 2", release the key; pressing key MODE (short time) it is possible to select the functions "UNIT 2", "SETUP", "CAL E", "CAL INT", wrapping around at the end of the list. For the models without internal calibration the function "CAL INT" will not appear.

The ON/OFF key activates the selected submenu, while the TARE key exits the main menu.

When the submenu is entered, the keys operate as above: pressing ON/OFF the selection is confirmed, while pressing TARE the selection is left unchanged.

NOTE: If no key is activated for 20 seconds, the menu level is incremented up to the end, without altering any parameter and the balance resumes its normal operation.

The Main menu flow is the following:

UNIT 2

submenu defining the additional operating functions of the balance.

SETUP

function used to configure the balance


CAL E

activation of the external calibration (manual calibration with user supplied weight)

CAL Int

activation of the internal calibration (automatic calibration with built-in weight) if present.

9.1. UNIT 2

With this function the balance can be programmed to operate with different measuring units. CRYSTAL can weight also in ounce (OZ), pounds (LB) and carats (CT). The balance can also be programmed for operation as:

STANDARD BALANCE
 PIECE-COUNTING BALANCE
 PERCENTAGE BALANCE

Select the function and confirm with the ON-OFF key.
 Select the desired unit or operating mode, using the MODE key.
 Confirm the unit/function selected with the ON-OFF key.
 Now the balance is operating with the new unit or function.

9.1.1. STANDARD BALANCE FUNCTION

When not operating as piece-counting or percentage balance the keys act in normal mode. The key MODE (activated for less than 1 second) switches between gram and the second unit. If a measuring unit different from gram (main operating unit) is selected, the number of decimal places is changed and the resolution of the balance (total number of divisions) is reduced, so that the last digit always exhibits a 1 count increment. When a second unit is selected the balance automatically stores the information so that, even if the power is removed, the selection is retained and, when the balance is switched on, the selected unit will be active.

9.1.2. PIECE COUNTING FUNCTION

The function is selected as above with the menu item PCS. A new reference may be acquired pressing the MODE key for more than 1 second. The display will show 10 PCS blinking: this is the first sample acquisition quantity. Pressing the MODE key for less than 1 second, other sample quantities, 20 PCS and 50 PCS, are selected. Now the reference quantity may be weighted, pressing the ON-OFF key, when the stable symbol lights on the display, the new reference value is acquired and the balance starts to display the piece count.

Note:

If the weight of the sample is too small, acquisition is interrupted and the message "PCS ERROR" appears until the TARE key is activated.

9.1.3. PERCENTAGE FUNCTION

The function is selected as above with the menu item %.

A new reference may be acquired pressing the MODE key for more than 1 second.

The display will show – 100 – blinking.

Now the reference quantity may be weighted, pressing the ON-OFF key, when the stable symbol lights on the display, the the new reference value is acquired and the balance starts to display the weight as a percentage of the reference acquired.

Notes:


If the weight of the sample is too small, acquisition is inhibited and the message “PCS ERROR” appears until the TARE key is activated.

The balance can measure percentages up to 500 % of the reference if the full scale is not exceeded.

Example:

A reference weight of 180 g is acquired as reference and the displays shows 100.00 %, then a weight of 35 g is placed on the pan and the display now shows 19,44 %, if the weight is increased up to 95 g, the display will show 52,77 %.

9.2. SETUP

9.2.1. SERIAL

This menu controls the serial output control parameters.

9.2.1.1. BAUD

Defines the transmission baudrate of the serial output from 1200 to 9600 baud.

9.2.1.2. BYTE FOR

Defines the format of the byte transmitted and may be set as:

7 E 1	7 data bit	even parity	1 stop bit
7 O 1	7 data bit	odd parity	1 stop bit
7 N 2	7 data bit	no parity	2 stop bit
7 E 2	7 data bit	even parity	2 stop bit
7 O 2	7 data bit	odd parity	2 stop bit
8 N 1	8 data bit	no parity	1 stop bit
8 E 1	8 data bit	even parity	1 stop bit
8 O 1	8 data bit	odd parity	1 stop bit
8 N 2	8 data bit	no parity	2 stop bit

9.2.1.3. PROTOC

Is the protocol of the serial data output and may be:

EUROPE
CRYSTAL
PRINTER

For a detailed description of the different protocols refer to the section relative to the data transmission.

9.2.1.4. TRANS OPT

Selects the continuous or on-demand data transmission.

When the continuous mode is selected, the balance transmits a data at every conversion.

May be used to drive peripherals than can't issue a print command to the balance.

This function can be stored by STORE function (see 9.2.5)

9.2.2. CALIBR

The Crystal balance may be calibrated with an external mass or with the built-in calibration mass, the type of calibration is a function of the model and can be selected with this menu.

If the balance has a built-in calibration mass the following items are displayed:

- 9.2.2.1. CAL OFF** → Both external and internal calibration are disabled
- 9.2.2.2. CAL ON** → Both external and internal (if installed) calibration are enabled
- 9.2.2.3. CAL AUTO** → With this selection the calibration is enabled and a timer activates the internal calibration every 3 hour. The automatic calibration is delayed if, during last 10 minutes, the balance has been operated or serial command has been issued.

9.2.3. A-ZERO

If the Auto-Zero function is active and when no weights are on the pan, the automatic zero tracking correction is active: the balance remains at a net value equal to zero. The function is useful to cancel small drift of the zero reading.

9.2.3.1. ZERO OFF

The zero tracking function is inactive.

9.2.3.2. ZERO ON

The function is active and small changes in zero reading are corrected.

9.2.4. FILTER

The setting of the filter menu controls the dynamic response of the balance. Selections are.

9.2.4.1. SLO slow response applicable to disturbed environments (high accuracy weighing)

9.2.4.2. FAST fast response applicable to quickly changing weights (dosing)

9.2.5. SET DEF

All setting, changed by the setup menu, are valid as soon as they are changed. However if the balance is powered off they are lost and when powered up the balance will use different parameters. To maintain a custom configuration that will be active at power-on the function SET DEF was implemented.

9.2.5.1. STORE → current settings are stored as user default

9.2.5.2. END → user default setting is left unchanged

9.3. CAL E

Once the external calibration function is started the display shows "LOAD" and the reference weight acquisition is active: a reference weight of E2 class is required for calibration, with value function of the range, the value of the weight is automatically detected and checked by the balance. When the value of the reference weight is acquired, the display shows the "UNLOAD" message: at this time the weight must be removed from the pan and the balance enters the zero acquisition phase. At the end of the cycle the message "CALIBRATED" appears if the calibration is valid and the new calibration factor is stored, if anything fails, the message "CAL ERROR" appears and the calibration factor of the balance remains unchanged.

Notes:

Balances with ranges higher than 100 g use a calibration weight of 100 g and multiples.

Balances with ranges lower than 120 g use a calibration weight of 50 g and multiples.

9.4. CAL INT

If the balance is equipped with the internal calibration mass the message "CAL INT" appears as a choice of the main menu and the internal automatic calibration is available. Once this function is activated, the automatic calibration process begins.

The display shows "CAL AUTO" during the calibration process while the functions of reference and zero weight acquisition are performed automatically by the balance.

During the calibration the process is monitored for proper execution.

At the end of the cycle the message "CALIBRATED" or "CAL ERROR" appear as a result of the calibration process as above.

9.2.5.1. STORE → current settings are stored as user default

9.2.5.2. END → user default setting is left unchanged

9.5. TIMED AUTOMATIC CALIBRATION

If the balance is equipped with the internal calibration mass, the timed calibration periodically starts an internal automatic calibration.

In order to operate, 3 hours must be elapsed from the last calibration, the balance must be inactive for at least 10 minute before the calibration time, no serial command was send to the balance.

The display shows "CAL AUTO" during the calibration process while the functions of reference and zero weight are performed automatically by the balance.

At the end of the cycle the message " CALIBRATED " or "CAL ERROR" appear as a result of the calibration process as above.

In this way the operator is able to determine the status of the last timed calibration.

To restore the normal operation the key TARE must be activated.

10. QUICK RESTORE OF DEFAULT PARAMETERS

If configuration parameters are modified, so that a wrong setup is stored, a quick restore function of the factory parameters is available.

Remove the power supply plug from the rear of the balance

Press and hold pressed the ON/OFF key

Reinsert the power supply plug

Wait until the display shows "Crystal"

Release the ON/OFF key

The display will show "SET FACTORY"

Factory default parameters are restored

Note:


When the balance configuration was modified and its operation is correct, the use of the SET DEF command of the setup menu will permanently store the new configuration parameters in memory so that they are not lost on power down.

11. SERIAL OUTPUT

The serial output of the balance uses the EIA RS232 standard
Balance and computer/ peripheral settings must be the same.

11.1. SERIAL OUTPUT CONNECTION

The connection of the balance is:

11.2. DATA TRANSMISSION FORMAT

The format of the data transmitted by the balance is different as a function of the selected protocol by the menu option 9.2.1.3.
EURO and CRYSTAL protocols are suitable for computer interfacing, the PRINTER protocol may be used to drive external printers.

CRYSTAL Protocol

The answer to any type of command is formatted with a measure field, a unit field and a status flag field as follows:

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	B	F1	F2	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

D0 ...D9	Measure	with sign & decimal point, right justified
B	Blank	Space between measure and unit
U	Unit	Measuring unit
F1, F2	Status	Flags indicating the status of the balance
CR	Carriage Return	ASCII code <CR>
LF	Line Feed	ASCII code <LF>

The status flags values of F1 are:

D	Transmitted data is valid
O	Balance in over-range
U	Balance in under-range
T	Tare in progress
C	Calibration in progress
Z	Zero acquire
I	Init system

The status flags values of F2 are:

S	Measure is stable
I	Measure is not stable
E	Measure is in error
A	Automatic calibration acknowledged
B	Calibration in progress

EURO protocol

The format of the euro protocol is:

D0	D1	D2	D3	D4	D5	D6	D7	D8	CR	LF
1	2	3	4	5	6	7	8	9	10	11

D0 ... D8	Measure	with sign & decimal point
CR	Carriage Return	ASCII code <CR>
LF	Line Feed	ASCII code <LF>

PRINTER protocol

Is activated by the "PRINT" pushbutton. The output string uses 10 characters for the measure including sign and decimal point, the value is right justified and filled with spaces from the left. The format of the output is:

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Do to D9	Measure	with sign & decimal point, right justified
B	Blank	Space between measure and unit
U	Unit	Measuring unit
CR	Carriage Return	ASCII code <CR>
LF	Line Feed	ASCII code <LF>

Note :

Data will be printed if and only if the weight on pan is stable.

11.3. REMOTE CONTROL COMMANDS

The balance responds to the following commands:

Command	Function
B	Send data on display
C	Perform calibration
F	Stop continuous transmission of data on display
I	Start continuous transmission of data on display
T	Tare command

Every remote control command must end with <CR>.

Line-feed characters are considered don't care symbols

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