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This warranty certificate entitles either "free-of-charge" replacement or repair of the components recognized faulty by the manufacturer, provided the instrument is returned "freight-pre-paid" to the distributor.

The following are excluded:

- ALL GLASS-WARE

- Damages caused during the shipment
- Faulty power connections
- Accidental unforeseen faults
- Damages caused by improper use and/or not in conformity with the enclosed instructions
- Product liability is revoked if anyone tampers with the instrument or its additional parts anyway

**FOR RESELLER
ONLY**

Purchaser's name
Purchaser's address
Zip City Country.....
Date of purchase..... Seller's name.....

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WARRANTY
 WARNING
 DECLARATION OF CONFORMITY

INDEX

1. WARNINGS	16
2. TECHNICAL CHARACTERISTICS	17
3. LOCATION	17
4. INSTALLATION	18
5. OPERATING CONTROLS OF THE BALANCE	19
6. OPERATING INSTRUCTIONS	20
7. KEYBOARD OPERATION	20
8. TARE	20
9. MAIN MENU AND OPERATING FUNCTIONS	21
9.1. UNIT 2	21
9.1.1. STANDARD BALANCE FUNCTION	22
9.1.2. PIECE COUNTING FUNCTION	22
9.1.3. PERCENTAGE FUNCTION	22
9.2. SETUP.....	23
9.2.1. SERIAL	24
9.2.2. CALIBR.....	24
9.2.3. A-ZERO	25
9.2.4. FILTER	25
9.2.5. SET DEF	25
9.3. CAL E	25
10. QUICK RESTORE OF DEFAULT PARAMETERS	26
11. SERIAL OUTPUT	26
11.1. SERIAL OUTPUT CONNECTION.....	26
11.2. DATA TRANSMISSION FORMAT	26
11.2.1. CRYSTAL protocol	26
11.2.2. EURO protocol	27
11.2.3. PRINTER protocol.....	27
11.3. SERIAL COMMANDS	28

NOTE

**Read the Instructions carefully before operating the Balance.
This is an high accuracy Instrument and must be handled with care.**

1. WARNINGS

**The power supply pack of the Crystal 1000 complies with Electrical Protection Class II.
Before connecting the power, check that the operating voltage, indicated on the Power supply pack, matches the line voltage.**

CRYSTAL 1000:

MAY ONLY BE USED WITH THE POWER SUPPLY PACK SUPPLIED.

MUST ONLY BE OPERATED IN DRY, CLEAN AND VIBRATION-FREE ENVIRONMENT.

MUST NOT BE OPERATED IN AREAS WHERE THERE ARE EXPLOSION RISKS.

NOTE:

In case of improper setting of the balance parameters during setup, the balance may power-up in an incorrect configuration.

The user can restore the factory configuration to avoid the malfunction as follows:

- ◆ Disconnect the supply of the balance
- ◆ Hold the On/Off key pressed
- ◆ Reconnect the supply
- ◆ Release the On/Off key when the display shows 'Set Factory'

Now the factory default settings are restored.

2. TECHNICAL CHARACTERISTICS

	CRYSTAL 1000
Capacity & Resolution	1010 g / 1 mg
Linearity	± 2 mg
Calibration Weight	500 g or multiples
Stabilisation time (fast/slow)	2 / 4 S
Operating temperature	18 - 28 °C (10 - 30 °C at reduced accuracy)
Pan diameter	Ø 110 mm
Operating Voltage	100 to 240 V -15%/+10% (range 85 to 264 V) 50 / 60Hz
Power consumption	200 mA
Data output	EIA RS232 I/O
Total dimensions	215 x 350 x 310 mm (w x d x h)
Weighing chamber dimensions	210 x 155 x 220 mm (w x d x h)
Weight	7 Kg (power supply included)
Internal Calibration	NO

3. LOCATION

Choose an appropriate location taking into account these criteria:

Vibrations and air drafts in the ambient must be avoided.

The supporting stand must be rigid.

The mains line must be stable. It is not advisable to use the same electric line of the fluorescent lamps or other electric apparatus which can be switched on and off randomly.

The balance must be kept far from magnetic fields.

Environmental conditions with wide temperature changes must be avoided.

Ambient temperature above or below temperature in the table may affect the operation of the balance and its accuracy.

In very dry environments electrostatic charges may build up, especially if glass or porcelain objects are used and even more with plastic (plastic should be avoided in any case).


WARNING !

Both the balance pan and the sample to be weighed must be very clean and kept at the same temperature.

4. INSTALLATION

Instructions to remove the locking device of the balance

- ◆ Unpack the balance and all its accessories.
- ◆ Unscrew the knob on the locking bracket as indicated by the arrow (anti clockwise).
- ◆ Retain the locking bracket in case the balance has to be returned to the supplier. This bracket can be refastened to the balance by reolacing the knob.


Level the balance using the spirit level on the rear and rotating the levelling feet under the front part of the basement (Fig 3).

Insert the pan.

Insert the power supply cable in the plug located on the rear panel (Fig 2).

Insert the power supply in the main line plug checking the supply voltage printed on the power supply case (100 / 240 Vac) - (Fig 2).

FIG. 2 - Rear panel details

5. OPERATING CONTROLS OF THE BALANCE


FIG. 3 – General view

- 1- Main Display
- 2- On/Off Pushbutton
- 3- Tare Pushbutton
- 4- Mode Pushbutton
- 5- Print Pushbutton
- 6- Weighing Pan
- 7- Levelling feet

6. OPERATING INSTRUCTIONS

Wait for about 3 hours after power-up before performing any high accuracy weighing operation. It is advisable to leave the balance permanently on (connected to the AC line). It is possible to switch off only the display with the ON/OFF key, this allows a quick recovery when the balance is switched on.

7. KEYBOARD OPERATION

The keyboard has four control keys: PRINT, MODE, TARE, ON-OFF. Every key has a specific function, while keys MODE, TARE and ON-OFF also activate a second function. If the keys are activated for less than 1 second the first function is activated, when pressed for more the second function will be active. The table below shows the operation capabilities.

Normal use

KEY NAME	KEY ACTIVATED FOR LESS THAN 1 SECOND	KEY ACTIVATED FOR MORE THAN 1 SECOND
On/Off	Display illumination on/off	Balance on/off
Mode	See additional functions	See additional functions
Tare	Tare function	Tare function
Print	Print function	No function

Main & configuration menu

KEY NAME	KEY ACTIVATED FOR LESS THAN 1 SECOND	KEY ACTIVATED FOR MORE THAN 1 SECOND
On/Off	Confirm selection	No function
Mode	Scroll available items	No function
Tare	Exit current selection	No function
Print	No function	No function

8. TARE

After a container is placed on the pan, push the TARE key to zero the reading. The balance has a tare capacity up to its total measuring range. Place the sample to be weighed in the container. The display shows the sample net weight. It is possible to use the tare key again and again after adding several samples, up to the maximum capacity of the balance, indicated by a bargraph on the left side of the LCD display.

9. MAIN MENU AND OPERATING FUNCTIONS

The balance can be operated in the following configurations:

Balance with gram, carat, ounce, pound measuring unit
 Piece counting balance
 Percentage computing balance

The selection is performed as follows:


Hold MODE key pressed until the display shows "UNIT 2", release the key; pressing key MODE (short time) it is possible to select the functions "SETUP" and "CAL E".

The ON/OFF key activates the selected submenu, while the TARE key exits the main menu.

When the submenu is entered, the keys operate as above: pressing ON/OFF the selection is confirmed, while pressing TARE the selection is left unchanged.

NOTE: If no key is activated for 20 seconds, the menu level is incremented up to the end, without altering any parameter and the balance resumes its normal operation.

The Main menu flow is the following:

UNIT 2

submenu defining the additional operating functions of the balance.


SETUP

function used to configure the balance

CAL E

activation of the external calibration (manual calibration with user supplied weight)

9.1. UNIT 2

With this function the balance can be programmed to operate with different measuring units. CRYSTAL 1000 can weight also in ounce (OZ), pounds (LB) and carats (CT).

The balance can also be programmed for operation as:

STANDARD BALANCE
PIECE-COUNTING BALANCE
PERCENTAGE BALANCE

Select the function and confirm with the ON-OFF key.
Select the desired unit or operating mode, using the MODE key.
Confirm the unit/function selected with the ON-OFF key.
Now the balance is operating with the new unit or function.

9.1.1. STANDARD BALANCE FUNCTION

When not operating as piece-counting or percentage balance the keys act in normal mode.
The key MODE (activated for less than 1 second) switches between gram and the second unit.
If a measuring unit different from gram (main operating unit) is selected, the number of decimal places is changed and the resolution of the balance (total number of divisions) is reduced, so that the last digit always exhibits a 1 count increment.
When a second unit is selected the balance automatically stores the information so that, even if the power is removed, the selection is retained and, when the balance is switched on, the selected unit will be active.

9.1.2. PIECE COUNTING FUNCTION

The function is selected as above with the menu item PCS.
A new reference may be acquired pressing the MODE key for more than 1 second.
The display will show 5 PCS blinking: this is the first sample acquisition quantity.
Pressing the MODE key for less than 1 second, other sample quantities, 10 PCS, 20 PCS, 50 PCS and 100 PCS, are selected.
Now the reference quantity may be weighted, pressing the ON-OFF key, when the stable symbol lights on the display, the new reference value is acquired and the balance starts to display the piece count.

Note:

If the weight of the sample is too small, acquisition is interrupted and the message "PCS ERROR" appears until the TARE key is activated.

9.1.3. PERCENTAGE FUNCTION

The function is selected as above with the menu item %.
A new reference may be acquired pressing the MODE key for more than 1 second.
The display will show – 100 – blinking.
Now the reference quantity may be weighted, pressing the ON-OFF key, when the stable symbol lights on the display, the new reference value is acquired and the balance starts to display the weight as a percentage of the reference acquired.

Notes:

If the weight of the sample is too small, acquisition is inhibited and the message "PCS ERROR" appears until the TARE key is activated.
The balance can measure percentages up to 500 % of the reference if the full scale is not exceeded.


Example:

If a reference weight of 180 g is acquired the display shows 100.00 %. Now, if a weight of 35 g is placed on the pan the display shows 19,44 %. If the weight is increased up to 95 g, the display will show 52,77 %.
Now, if the weight is increased more than 900 g the display will show the message OVER PERC.
This message means that whole percentage function capacity (500 %) has been exceeded.

Pressing the key MODE, you can switch from main unit to percentage unit to know the weight value loaded on the pan.

To go out from OVER PERC condition, remove the weight from the pan.

9.2. SETUP

9.2.1. SERIAL

This menu controls the serial output control parameters.

9.2.1.1. BAUD

Defines the transmission baudrate of the serial output from 1200 to 9600 baud.

9.2.1.2. BYTE FOR

Defines the format of the byte transmitted and may be set as:

DATA FORMATS	DATA BIT	PARITY CHECK TYPE	BIT STOP
7 E 1	7	EVEN	1
7 O 1	7	ODD	1
7 N 2	7	NO PARITY	2
7 E 2	7	EVEN	2
7 O 2	7	ODD	2
8 N 1	8	NO PARITY	1
8 E 1	8	EVEN	1
8 O 1	8	ODD	1
8 N 2	8	NO PARITY	2

9.2.1.3. PROTOC

Is the protocol of the serial data output and may be:

EUROPE
CRYSTAL
PRINTER

For a detailed description of the different protocols refer to the section relative to the data transmission.

9.2.1.4. TRANS OPT

Selects three different communication mode:

- continue** When this mode is selected, the balance transmits a data at every conversion. May be used to drive peripherals than can't issue a print command to the balance.
- request** The balance transmits a data one more time when a valid command is issued by serial port.
- data stable** The balance transmits a data one more time when a valid command is issued by serial port if and only if the weight on the pan is stable. If the non-stable condition persist then, after 15 seconds, whatever there will be a transmission.

This parameter could be stored by STORE function (see 9.2.5)

9.2.2. CALIBR

The calibration function can be selected with this menu.

- 9.2.2.1. CAL OFF** → Both external and internal calibration are disabled
- 9.2.2.2. CAL ON** → Both external and internal (if installed) calibration are enabled

9.2.3. A-ZERO

If the Auto-Zero function is active and when no weights are on the pan, the automatic zero tracking correction is active: the balance remains at a net value equal to zero. The function is useful to cancel small drift of the zero reading.

9.2.3.1. ZERO OFF

The zero tracking function is inactive.

9.2.3.2. ZERO ON

The function is active and small changes in zero reading are corrected.

9.2.4. FILTER

The setting of the filter menu controls the dynamic response of the balance. Selections are.

9.2.4.1. SLO slow response applicable to disturbed environments (high accuracy weighing)

9.2.4.2. AVG average value for both high and low accuracy

9.2.4.3. FAST fast response applicable to quickly changing weights (dosing)

9.2.5. SET DEF

All setting, changed by the setup menu, are valid as soon as they are changed.

However if the balance is powered off they are lost and when powered up the balance will use different parameters.

To maintain a custom configuration that will be active at power-on the function SET DEF was implemented.

9.2.5.1. STORE → current settings are stored as user default

9.2.5.2. END → user default setting is left unchanged

9.3. CAL E

Once the external calibration function is started the display shows "LOAD" and the reference weight acquisition is active: a reference weight is required for calibration, with value function of the range, the value of the weight is automatically detected and checked by the balance.

When the value of the reference weight is acquired, the display shows the "UNLOAD" message: at this time the weight must be removed from the pan and the balance enters the zero acquisition phase.

At the end of the cycle the message "CALIBRATED" appears if the calibration is valid and the new calibration factor is stored, if anything fails, the message "CAL ERROR" appears and the calibration factor of the balance remains unchanged.

NOTE - The calibration process not starts if a weight of more than 2 gram is on the pan at calibration start time.

10. QUICK RESTORE OF DEFAULT PARAMETERS

If configuration parameters are modified, so that a wrong setup is stored, a quick restore function of the factory parameters is available.

Remove the power supply plug from the rear of the balance
 Press and hold pressed the ON/OFF key
 Reinsert the power supply plug
 Wait until the display shows "Crystal"
 Release the ON/OFF key
 The display will show "SET FACTORY"
 Factory default parameters are restored

Note:

When the balance configuration was modified and its operation is correct, the use of the SET DEF command of the setup menu will permanently store the new configuration parameters in memory so that they are not lost on power down.

11. SERIAL OUTPUT

CRYSTAL serial output is compatible with EIA RS232 standard.
 Every answer from balance is sent every 100 ms.

11.1. SERIAL OUTPUT CONNECTION

BALANCE	DB9 CONNECTOR (9 pole)	DB25 CONNECTOR (25 pole)
RS 232 Out	Pin 2	Pin 3
RS 232 In	Pin 3	Pin 2
GROUND	Pin 5	Pin 7

11.2. DATA TRANSMISSION FORMAT

The format of the data transmitted by the balance is different as a function of the selected protocol by the menu option 9.2.1.3.

EURO and CRYSTAL protocols are suitable for computer interfacing; PRINTER protocol may be used to drive external printers.

11.2.1. CRYSTAL protocol

The answer to any type of command is formatted with a data field, a unit field and a status field as follows:
 a measure field
 a unit field
 a status field

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	B	F1	F2	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

The measure field (1..10) is the reading on the display.
 It is right justified with decimal point and measure sign.
 Unit field shows the current unit weight
 Status field shows the current balance status
 B is the blank character.
 Field N°18 and 19 are:

CR carriage returns ASCII code.

LF line feed ASCII code.

Status information of the F1 character

D	Transmitted data is valid
O	Balance in over range
U	Balance in under range
T	Tare
C	Calibration process in progress
Z	Zero acquisition
I	Init test

Status information of the F2 character

S	Measure is stable
I	Measure is not stable
E	Error
A	Automatic calibration (when F1 = C)
L	Load weight(manual calibration) (when F1 = C)
U	UnLoad weight (manual calibration) (when F1 = C)
D	Calibration done (when F1 = C)
B	Calibration busy(already started)
O	Calibration OFF
P	Percentage

11.2.2. EURO protocol

The answer to any type of command is formatted as follow:

D0	D1	D2	D3	D4	D5	D6	D7	D8	CR	LF
1	2	3	4	5	6	7	8	9	10	11

The measure field (1..9) is the display reading.
It is right justified with decimal point and the measure sign.

Field N°10 and 11 are:
CR carriage returns ASCII code.
LF line feed ASCII code.

11.2.3. PRINTER protocol

Is activated by the "PRINT" pushbutton. The output string uses 10 characters for the measure including sign and decimal point, the value is right justified and filled with spaces from the left. The format of the output is:

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

The measure field (1..10) is the display reading.
It is right justified and has the decimal point and the measure sign.
Unit field shows the current unit weight
B is the blank character.
Field N°15 and 16 are:
CR carriage returns ASCII code.
LF line feed ASCII code.

Note: Data will be printed if and only if the weight on pan is stable.

11.3. SERIAL COMMANDS

The balance responds to the following commands only if a CRYSTAL or EURO protocol is selected:

SERIAL REMOTE CONTROLS TABLE

COMMAND	FUNCTION
B	Send data on display
C	Start calibration
I	Start continuous transmission of data on display
F	Stop continuous transmission of data on display
T	Tare command

Notes:

Every remote control command must end with <CR>.
Line-feed characters are ignored.

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11. SERIAL OUTPUT

CRYSTAL serial output is compatible with EIA RS232 standard.

Every answer from balance is sent every 100 ms.

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11.2.1. CRYSTAL protocol

The answer to any type of command is formatted with a data field, a unit field and a status field as follows:

a measure field
 a unit field
 a status field

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	B	F1	F2	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

The measure field (1..10) is the reading on the display.

It is right justified with decimal point and measure sign.

Unit field shows the current unit weight

Status field shows the current balance status

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The measure field (1..10) is the display reading.
It is right justified and has the decimal point and the measure sign.

Unit field shows the current unit weight

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Notes:

Every remote control command must end with <CR>. Line-feed characters are ignored.

customer care

info@gibertini.com

Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

**GIBERTINI
ELETTRONICA s.r.l.**

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434
Fax (+39) 02 3541438
www.gibertini.com
sales@gibertini.com

info@gibertini.com
Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434
Fax (+39) 02 3541438
www.gibertini.com
sales@gibertini.com