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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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

ENGLISH

NOTE

Read the Instructions carefully before operating the Microbalance
This is an high accuracy Instrument and must be handled with care.

1. WARNINGS

The power supply pack of the MicroCrystal microbalance complies with Electrical Protection Class II.

Before connecting the power, check that the operating voltage, indicated on the Power supply pack, matches the line voltage.

MICROCRYSTAL MICROBALANCE:

- MAY ONLY BE USED WITH THE POWER SUPPLY PACK SUPPLIED
- MUST ONLY BE OPERATED IN DRY, CLEAN AND VIBRATION-FREE ENVIRONMENT
- MUST NOT BE OPERATED IN AREAS WHERE THERE ARE EXPLOSION RISKS.

NOTE

In case of improper setting of the microbalance parameters during set-up, the microbalance may power-up in an incorrect configuration.

The user can restore the factory configuration to avoid the malfunction as follows:

- Disconnect the supply of the microbalance
- Hold the On/Off key pressed
- Reconnect the supply
- Release the On/Off key when the display shows 'Set Factory'

Now the factory default settings are restored.

NOTE

THE GLASS FOR PAN MUST BE CLEANED ONLY WITH CLOTH “**ANTISTATICUM**” OR SIMILAR TO ELIMINATE STATIC ELECTRICITY.

2. TECHNICAL CHARACTERISTICS

FUNCTIONS

- LCD display with small size digits, for easier and more immediate reading
- Waterproof and acid resistant membrane keyboard. Easy to use with the ON/OFF, TARE and MODE keys
- Indication of reached stable weight
- Bar-graph indicator of dosage and remaining capacity of the microbalance
- Parameters configurable by menu: reading in g (grams), lb (pound), oz (ounce), ct (carats), pcs (pieces), % (percentage)

- Capacity: 1000 mg
- Readability: 0.001 mg
- Repeatability: ± 0.001 mg
- Linearity: ± 0.002 mg
- Response time: 6 seconds
- Data output: RS232
- Pan diameter: 20 mm standard and for filters max 60 mm
(on request: bigger diameter)
- Automatic calibration with external mass (supplied with standard equipment)
- Functioning temperature: 20°C ÷ 30°C (corrected with a micro-processor)
- Power supply: 100/240 VAC
- Absorption: 200 mA
- Dimensions (WxDxH): 215 x 385 x 230 mm
- Dimensions (WxDxH): 240 x 180 x 120 mm
- Net weight: 5.2 kg

STANDARD EQUIPMENT

Mass in class E1 for the end of range calibration; grab/forceps for the mass, antistatic cloth

ON REQUEST

USB data output

Ion generator for the elimination of the electrostatic charges

3. LOCATION

Choose an appropriate location taking into account these criteria:


- Vibrations and air drafts in the ambient must be avoided.
- The supporting stand must be rigid.
- The mains line must be stable. It is not advisable to use the same electric line of the fluorescent lamps or other electric apparatus which can be switched on and off randomly.
- The microbalance must be kept far from magnetic fields.
- Environmental conditions with wide temperature changes must be avoided. Ambient temperature above 28 °C or below 18 °C may affect the operation of the microbalance and its accuracy.
- Avoid very dry environmental conditions, with very low Relative Humidity. In these conditions, electrostatic charges build up easily on materials such as porcelain and glass, and even more on plastic materials.

4. INSTALLATION

Install the microbalance as follows:

- Unpack the microbalance and all its accessories.
- Insert the weighing pan.
- Level the microbalance using the spirit level on the rear of the microbalance and rotating the levelling feet under the front part of the microbalance basement.
- Insert the power supply cable in the connector located on the rear of the microbalance
- Insert the power supply in the mains line plug, checking that the voltage it's the same printed on the power supply case.


At first start up, after the installation, wait for at least 3 hours before performing any weighing. The microbalance is thus allowed to stabilise mechanically and thermally, from possible stresses due to transportation or other factors.

5. OPERATING CONTROLS OF THE MICROBALANCE

Control panel (with test display visualization)

1. Bargraph
2. Display
3. Key PRINT
4. Key MODE
5. Key TARE
6. Key ON-OFF

6. OPERATING INSTRUCTIONS

Wait at least 1 hour before performing any weighing at every daily start up.

It is advisable to leave the microbalance always ON (connected to the mains), in order to have the microbalance always thermally stabilised and ready for use. In this way it is possible to avoid the thermal stabilisation time.

In this way, to limit the electrical consumption, it is possible to turn off the display light, by pressing for more than one second the ON/OFF key. The display will show the message "OFF". To reactivate the light, press again the ON/OFF key.

If the instrument is not used for more than one hour, it is advisable to perform a weighing with at least 100 milligrams before performing other high precision weighings.

ATTENTION

Both the microbalance pan and the sample to be weighed must be very clean and kept at the same temperature.

In presence of particularly dry climate, or in environments with low Relative Humidity, avoid using materials that cumulate easily electrostatic charged such as porcelain, glass and in particular plastic materials (these should be avoided at any rate). In case these materials must be used, clean accurately the parts with suitable antistatic products.

7. KEYBOARD OPERATION

The keyboard has four control keys: PRINT, MODE, TARE, ON-OFF.

Every key has a specific function, while keys MODE and ON-OFF also activate a second function.

If the keys are activated for less than 1 second the first function is activated, when pressed for more the second function will be active.

The table below shows the operation capabilities.

Normal use

KEY NAME	KEY ACTIVATED FOR LESS THAN 1 SECOND	KEY ACTIVATED FOR MORE THAN 1 SECOND
On/Off	Display illumination on/off	Microbalance on/off
Mode	See additional functions	See additional functions
Tare	Tare function	No function
Print	Print function	No function

Main & configuration menu

KEY NAME	KEY ACTIVATED FOR LESS THAN 1 SECOND	KEY ACTIVATED FOR MORE THAN 1 SECOND
On/Off	Confirm selection	No function
Mode	Scroll available items	No function
Tare	Exit current selection	No function
Print	No function	No function

8. TARE

After a container is placed on the pan, push the TARE key to zero the reading.

The microbalance has a tare capacity up to its total measuring range.

Place the sample to be weighed in the container.

The display shows the sample net weight.

It is possible to use the tare key again and again after adding several samples, up to the maximum capacity of the microbalance, indicated by a bargraph on the left side of the LCD display.

9. MAIN MENU AND OPERATING FUNCTIONS

The microbalance can be operated in the following configurations:

The selection is performed as follows:


Hold MODE key pressed until the display shows "SET UP", release the key; pressing key MODE (short time) it is possible to select the functions "CAL E", wrapping around at the end of the list.

The ON/OFF key activates the selected submenu, while the TARE key exits the main menu.

When the submenu is entered, the keys operate as above: pressing ON/OFF the selection is confirmed, while pressing TARE the selection is left unchanged.

NOTE: If no key is activated for 20 seconds, the menu level is incremented up to the end, without altering any parameter and the microbalance resumes its normal operation.

The Main menu flow is the following:


SETUP

function used to configure the microbalance

CAL E

activation of the external calibration (manual calibration with user supplied weight)

9.1 SETUP

9.1.1. SERIAL

This menu controls the serial output control parameters.

9.1.1.1.1 BAUD

Defines the transmission baud rate of the serial output from 1200 to 9600 baud.

9.1.1.2 BYTE FOR

Defines the format of the byte transmitted and may be set as:

7 E 1	7 data bit	even parity	1 stop bit
7 O 1	7 data bit	odd parity	1 stop bit
7 N 2	7 data bit	no parity	2 stop bit
7 E 2	7 data bit	even parity	2 stop bit
7 O 2	7 data bit	odd parity	2 stop bit
8 N 1	8 data bit	no parity	1 stop bit
8 E 1	8 data bit	even parity	1 stop bit
8 O 1	8 data bit	odd parity	1 stop bit
8 N 2	8 data bit	no parity	2 stop bit

9.1.1.3 PROTOC

Is the protocol of the serial data output and may be:

CRYSTAL PRINTER

For a detailed description of the different protocols refer to the section relative to the data transmission.

9.1.1.4 TRANS OPT

Selects the continuous or on-demand data transmission.

When the continuous mode is selected, the microbalance transmits a data at every conversion.

May be used to drive peripherals that can't issue a print command to the microbalance.

This function can be stored by STORE function (see 9.2.5)

9.1.1.5 STABLE OPT

Data stable transmission:

- STABLE ON**
- With the EUROPE protocol the balance transmits only once, as soon as a valid command is sent, but only if the data of the balance is stable.
 - With the CRYSTAL protocol the balance transmits like the EUROPE protocol but, in presence of an unstable data, it will carry out however a transmission after 15 second. This is in order to not leave the system blocked.
- STABLE OFF** The balance transmits in whichever condition.

9.1.2. CALIBR (see pag. 24)

The MicroCrystal microbalance must be calibrated with an external mass. By this menu it is possible to inhibit the process calibration.

9.1.2.1 CAL OFF → Calibration disabled

9.1.2.2 CAL ON → Calibration enabled

9.1.3. A-ZERO

If the Auto-Zero function is active and when no weights are on the pan, the automatic zero tracking correction is active: the microbalance remains at a net value equal to zero. The function is useful to cancel small drift of the zero reading.

9.1.3.1 ZERO OFF

The zero tracking function is inactive.

9.1.3.2 ZERO ON

The function is active and small changes in zero reading are corrected.

9.1.4. FILTER

The setting of the filter menu controls the dynamic response of the microbalance. Selections are:

9.1.4.1 SLO

slow response applicable to disturbed environments (high accuracy weighing)

9.1.4.2 AVG

average value for both high and low accuracy

9.1.4.3 FAST

fast response applicable to quickly changing weights (dosing)

9.1.5. IDLE SET

9.1.5.1 IDLE ON

Enables the automatic switch-off of the display when the balance is unutilised for more than 10 minutes.

9.1.5.2 IDLE OFF

Disables the automatic switch-off of the display

9.1.5. SET DEF

All setting, changed by the set-up menu, are valid as soon as they are changed.

However if the microbalance is powered off they are lost and when powered up the microbalance will use different parameters.

To maintain a custom configuration that will be active at power-on the function SET DEF was implemented.

9.1.5.1 END → user default setting is left unchanged

9.1.5.2 STORE → current settings are stored as user default

9.2. CALIBRATION WITH MASS OF 1g E1 (calibration with external mass)

To activate this function operate as follows:

- press the MODE for 3 seconds – UNIT 2 appears – press again and SETUP appears – press again and CALIBRATION.
- To confirm press ON/OFF
- The glass protection is UP and display shows load
- Put the mass on the PAN – glass protection down
- The microbalance executes the calibration and the written UNLOAD appears.
- The glass protection in UP. Remove the mass from the pan.
- The glass protection down
- If all the phases are respected and the weights are correct, the calibration is accepted and the display will show for some seconds the message "CALIBRATED".
- If, on the contrary, there was any error, the display will show the message "CAL ERROR" and the previous calibration value will be maintained.

NOTE

The calibration process is interrupted if a weight of more than 2 mg is left on the pan.

10. QUICK RESTORE OF DEFAULT PARAMETERS

If configuration parameters are modified, so that a wrong set-up is stored, a quick restore function of the factory parameters is available.

- Remove the power supply plug from the rear of the microbalance
- Keep pressed the ON/OFF key
- Reinsert the power supply plug
- Wait until the display shows "Crystal"
- Release the ON/OFF key
- The display will show "SET FACTORY"
- Factory default parameters are restored

NOTE

When the microbalance configuration was modified and its operation is correct, the use of the SET DEF command of the set-up menu will permanently store the new configuration parameters in memory so that they are not lost on power down.

11. SERIAL OUTPUT

The serial output of the microbalance uses the EIA RS232 standard Microbalance and computer/ peripheral settings must be the same.

11.1. SERIAL OUTPUT CONNECTION

The connection of the microbalance is:

Micromicrobalance	DB9		DB25
RS232 out	2	→	3
RS232 in	3	←	2
GND	5	—	7

11.2. DATA TRANSMISSION FORMAT

The format of the data transmitted by the microbalance is different as a function of the selected protocol by the menu option 9.2.1.3.

EURO and CRYSTAL protocols are suitable for computer interfacing, the PRINTER protocol may be used to drive external printers.

CRYSTAL Protocol

The answer to any type of command is formatted with a measure field, a unit field and a status flag field as follows:

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	B	F1	F2	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

D0 ...D9	Measure	with sign & decimal point, right justified
B	Blank	Space between measure and unit
U	Unit	Measuring unit
F1, F2	Status	Flags indicating the status of the microbalance
CR	Carriage Return	ASCII code <CR>
LF	Line Feed	ASCII code <LF>

The status flags values of F1 are:

D	Transmitted data is valid
O	Microbalance in over-range
U	Microbalance in under-range
T	Tare in progress
C	Calibration in progress
Z	Zero acquire
I	Init system

The status flags values of F2 are:

S	Measure is stable
I	Measure is not stable
E	Measure is in error
A	Automatic calibration acknowledged
B	Calibration in progress

PRINTER protocol

Is activated by the "PRINT" pushbutton. The output string uses 10 characters for the measure including sign and decimal point, the value is right justified and filled with spaces from the left. The format of the output is:

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

Do to D9	Measure	with sign & decimal point, right justified
B	Blank	Space between measure and unit
U	Unit	Measuring unit
CR	Carriage Return	ASCII code <CR>
LF	Line Feed	ASCII code <LF>

NOTE

Data will be printed if and only if the weight on pan is stable.

11.3. REMOTE CONTROL COMMANDS

The microbalance responds to the following commands:

Command	Function
B	Send data on display
C	Perform calibration
F	Stop continuous transmission of data on display
I	Start continuous transmission of data on display
T	Tare command

- Every remote control command must end with <CR>.
- Line-feed characters are considered don't care symbols

customer care

info@gibertini.com
Tel. (++39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI
ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (++39) 02 3541434
Fax (++39) 02 3541438
www.gibertini.com
sales@gibertini.com