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**INSTRUMENTS FOR USE IN
 DANGEROUS ENVIRONMENTS**
 CESA 03 ATEX 308 Certification
 ICIM quality system certification EC N. 427
 Notified body 0425

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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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1.1.1 Scale

The system is constituted by:

1. Power supply, classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.
2. Balance classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.

1.1.2 Platform

The system is constituted by:

1. Power supply, classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.
2. Weighing Unit classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.
3. Display Unit classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.

1.1.3 OEM version

The system is constituted by:

1. Power supply, classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.
2. Weighing Unit classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.
3. Display Unit classified as category 2G apparatus suitable for use in Zone 1,2 Gas Group IIB.

1.1.4 Installation requirements

The mains voltage must be less than 250 Vac.

The power supply must be protected by a differential interrupting device able to disconnect the system with a differential current sensitivity of no more than 30 mA.

The power supply must be protected by a fault limiting device able to limit the maximum short circuit current to no more than 70 A for a duration of less than 50 msec.

The ambient operating temperature must be limited to the range 0°C to +50°C for safe operation.

The classification T5 indicates that the maximum surface temperature of the system is 100 °C.

Before powering up the system the control procedure listed below must be performed.

For detailed information refer to the standard EN 60079-14: Electrical apparatus for explosive gas atmospheres -- Part 14: Electrical installations in hazardous areas (other than mines).

1.2. MAINTENANCE

The normal maintenance procedure of Eex series must be performed following the steps:

- 1- Visual inspection of the installation
- 2- Verification of tightening of screws of the power supply and sealing of the filling opening
- 3- Verification of the integrity of cable entries and cable glands on the power supply
- 4- Verification of the integrity of cables
- 5- Verification of tightening of the fiber optic feed-troughs
- 6- Verification of the quality of the ground connection

-- WARNING --

IN ANY CASE THE POWER SUPPLY ENCLOSURE MUST NOT BE OPENED.

IF CLEANING OF THE SYSTEM IS REQUIRED, IT MUST BE PERFORMED IN A WAY TO AVOID ELECTROSTATIC DISCHARGE THAT CAN ACTIVATE AN EXPLOSIVE ATMOSPHERE.

FOR ANY REPAIR CONCERNING THE INSTRUMENT, THAT CAN AFFECT IN ANY WAY THE SAFETY OF THE EQUIPMENT, THE USER MUST RETURN IT TO THE PRODUCER. NO AUTHORIZATION IS GIVEN TO PERFORM ANY TYPE OF SERVICING ON SAFETY RELATED COMPONENTS OR PARTS.

ANY OPTION, LIKE THE RS232 CONVERTER, MUST BE INSTALLED IN SAFE AREA.

DETAILED INFORMATION ON MAINTENANCE CAN BE FOUND IN REFER TO THE STANDARD EN 60079-17: INSPECTION AND MAINTENANCE OF ELECTRICAL INSTALLATIONS IN HAZARDOUS AREAS (OTHER THAN MINES).


1.3. INSTRUCTION TO AVOID MACHANICAL RISKS

Don't exceed the maximum value of **SAFE OVERLOAD**, indicated for the specific type of instrument.

Don't subject the instrument to an **HIGH AND CONTINUOUS LEVEL OF VIBRATION**, so that the reading is not adequately stable.

BOTH SITUATIONS CAN PROVOKE IGNITION OF THE POTENTIALLY EXPLOSIVE ATMOSPHERE PRESENT IN THE AMBIENCE.

1.4. INSTALLATION PLAN

INSTALLATION PLAN

Note

1. If the power supply is connected to the line in Hazardous area Zone 1 or Zone 2, proper protection of the connection must be provided.
2. The protective device must comply with the requirements of installation listed in the manual.
3. The power supply must be installed in fixed position
4. Connection cables must be protected against damage and are not bent at sharp angles

2. HOW TO INSTALL

NOTE

BALANCES AND SCALES ARE PERFECTLY TESTED AND CALIBRATED IN OUR LABORATORIES BUT, LIKE ALL WEIGHING INSTRUMENTS, THEY ARE INFLUENCED BY THE GRAVITY FORCE THAT CHANGES ACCORDING TO THE GEOGRAPHIC ZONE AND THE ALTITUDE. THEREFORE WE SUGGEST TO CALIBRATE THEM IN THE INSTALLATION PLACE.

2.1. EUROPE / EURO-PAINTEX MODEL

Place the instrument on a horizontal, rigid, vibration-free flat, far from heat sources and draughts.
 Remove the locking device of the balance (fig. A) unscrewing the knob on the locking bracket (12) as indicated by the arrow (anti clockwise). Retain the locking bracket in case the balance has to be returned to the supplier.
 Fit the pan
 Connect the power supply cable to the balance (Figg. C-D)
 Connect the cable to the differential magneto-thermic circuit breaker
 Switch on the instrument. Before using the instrument allow 30 minutes (Warm - up)
 Calibrate the instrument with a calibration mass (See tables at Point 3.4).


FIG. A – Locking bracket

2.2. TMB MODELS UNTIL 300 kg

ATTENTION – The industrial scales from 25 to 300 kg don't need unlocking operations as they don't have locking devices.

Place the instrument on a horizontal, rigid, vibration-free flat, far from heat sources and draughts.
 Connect the power supply cable to TMB and TMB-display cable (Fig. E)
 Connect the cable to the differential magneto-thermic circuit breaker.
 Switch on the instrument. Before using the instrument allow 30 minutes (Warm - up).
 Calibrate the instrument with the calibration mass (See tables at Point 3.4).

2.3. TMB MODELS FROM 600 TO 3000 kg

Put the platform in a destination side and then:

Remove the screws and then the plate.

Place the platform and level it by levelling feet "1". Check that all feet touch the bench surface.

Unscrew the yellow bolts "2" till the pan is unlocked.

Unscrew the screw "3" and remove the locking bracket "4" that locks the cell "A".

Connect the reading unit cable "5" and the power supply cable "6" to the connector of the weighing unit case "A".

Unblock the platform lever screwing completely the screw "7B". **Don't move the screw "7A"**

Reassemble and fix the plate of the platform.

Connect the cable "5" to the reading unit.

Connect the main cable to the differential magneto-thermic circuit breaker

Switch on the instrument. Before using the instrument allow 30 minutes (Warm - up).

Calibrate the instrument with the calibration mass (See tables at Point 3.4).

Note – The locking devices can vary from model to model.


FIG. B – Unlocking and connections of platform 600÷3000 kg

3. BALANCE DESCRIPTION

3.1. EUROPE MODEL


FIG. C – EUROPE Ex balance: assembly view – control panel

- 1- Bar-graph / percent used capacity
- 2- Display
- 3- TARE key
- 4- MODE key
- 5- ON/OFF key - Start counting

3.2. EURO-PAINTEX MODEL



FIG. D - EURO-PAINTEX Ex Balance: assembly view-control panel

- 1- Bar-graph / percent used capacity
- 2- Display
- 3- TARE key
- 4- MODE key
- 5- ON/OFF key - Start coating

3.3. TMB MODELS



FIG. E - Assembly view – Control panel

- 1- Bar-graph / percent used capacity
- 2- Display
- 3- TARE key
- 4- MODE key
- 5- ON/OFF key - Start counting


3.4. CALIBRATION MASSES TABLE

TABLE 1 - EUROPE EX balance features				
MODELS	CAPACITY	READABILITY	PAN mm	CALIBRATED MASSE
EUROPE 500 Ex	510g	0.001 g	Ø 110	200 g
EUROPE 1700 Ex	1700g	0.01 g	Ø 150	1 kg
EUROPE 6000 Ex	6300 g	0.1 g	Ø 190	1 kg or multiples
EURO-PAINTEX 8.5	8500 g	0,1 g	Ø 190	1 kg or multiples

TABLE 2 - EUROPE EX platform features				
MODELS	CAPACITY	READABILITY	PAN mm	CALIBRATED MASS
TMB 25 Ex	25 kg	0,5 g	400 x 500	1 kg and multiples up to max capacity
TMB 35 AR Ex	35 kg	up to 6,5 kg = 0,1 g from 6,5 to 35 kg = 1 g	400 x 400	
TMB 50 Ex	50 kg	1 g	400 x 500	
TMB 150 Ex	150 kg	5 g	600 x 700	10 kg and multiples up to max capacity
TMB 300 Ex	300 kg	10 g	700 x 900	
TMB 600 Ex	600 kg	20 g	1250 x 1250	50 kg and multiples up to max capacity
TMB 1500 Ex	1500 kg	50 g	1250 x 1250	100 kg and multiples up to max capacity
TMB 1500 A Ex	1500 kg	50 g	1250 x 1500	
TMB 3000 Ex	3000 kg	100 g	1500 x 1500	
TMB 26 Ex CE	26 kg	0,5 g	400 x 500	1 kg and multiples up to max capacity
TMB 52 Ex CE	52 kg	1 g (0,1 barred)	400 x 500	

4. HOW TO INSTALL THE PLATFORM INTO THE CAVITY (pan at floor level)

- 1- Unscrew the screws "A" and remove the pan "E".
- 2- Bring the platform frame "B" carefully near the four corners (see the figure). Pay attention: DON'T REMOVE the two bolt "F".
- 3- Put the platform "B" into the cavity "I" and the frame "L". In the meantime insert the cable "C2" and "D1" through the conduit. Connect the cable "C2" of the power supply to platform.
- 4- Centre the platform "B" using the screw "H" mounted on the four sides of the platform. Fasten the nut after the pan has been centred in the middle of the floor cavity frame "L".
- 5- Referring to the onboard bubble level, carefully level the platform moving the four adjustable feet "G". Before fastening the nut of the four feet, help yourself with the surveyor's level in order to get it perfectly levelled with the frame "L".
- 6- Unblock the two bolts "F" and check that the up-frame must be free. Fix the pan "E" with the screw "A".
- 7- Connect the cable "D1" to readout "D" and the main cable to the differential magneto-thermic circuit breaker, following specific instructions.
- 8- Switch on the instrument. Before using the instrument allow 30 minutes (Warm - up).
- 9- CALIBRATE THE INSTRUMENT WITH THE CALIBRATION MASS (See tables at Point 3.4).

5. OPERATING INSTRUCTIONS

5.1. GENERAL CONSIDERATIONS

Place the balance on a rigid, horizontal and vibration-free support, away from air flow and heat sources, or wherever it should be ready to use.

Connect the power supply. The balance is now ON and the display shows a message formed by a letter followed by numbers. This indicates the instrument program version. This is an useful information for our after-sales service.

After about 2 seconds, the display shows:

HELLO

The instrument then performs a display test as segments start flashing.

-8888888

After automatic zero setting

0.0

the balance is ready to use.

Before starting operations, allow the instrument a 30 minute warm up time. In fact the electronic components of the instrument need a thermal stabilization time throughout weighing values are not accurate.

Always leave the instrument ON during the working day. By using the ON/OFF key it is possible to turn off the display. In this way the instrument is always ready to use, without the warm up period.

Besides the numeric values, the display can show two kinds of message: a steady message regarding a normal situation or a flashing message indicating an abnormal situation or requesting intervention by the user.

When an object is put on the pan, the display shows its mass value. If this value is higher than the capacity of the balance, the display shows a flashing message.

-HI-

WARNING

**Never put objects on the pan with a mass 50% higher than the capacity of the balance.
Even if the mass value is lower than the balance capacity, the object always must be put on the pan carefully.**

If for any reason, the pan is not correctly positioned, the display shows a flashing message:

-LO-

5.2. TARE KEY

With this key, it is possible to zero all values, positive or negative, up to the capacity of the balance.

This key is particularly useful when weighing powders or liquids or anything else requires a container.

Put the empty container on the pan and press the TARE key.

Remove it from the balance and fill it. During this operation, the display shows a negative value equivalent to the container mass.

When the filled container is replaced on the pan, the display will show the net weight.

5.3. BAR GRAPH

The capacity percentage display (Bar-graph) (6) is on the left side of the panel. Each line represents 10% of the capacity. This information points out the percent used capacity when the pan is loaded and the balance has just been zeroed.

5.4. CALIBRATION (MODE KEY)

Calibration is not necessary for the instrument operation. You can periodically test the balance with a calibrated mass supplied on request.

Go on as follow:

- Press the MODE key until the display shows the “CAL” message and then press the ON/OFF key.

- The display shows:

ACAL

- until the instrument stability is reached. Then the flashing message appears on the display:

LOAD

- Put the calibration mass on the pan.
- While the instrument is storing the mass value, the display shows the message

CALC

- and the mass value.

The calibration mass value can be the one indicated in “Balance features” or a multiple of it.

If the calibration mass is not correct, the display shows the flashing message

ErrC

The same message is shown when a mass is not put on the pan within 10 seconds from the “LOAD” message. In both situations press the TARE key. The instrument retains the last calibration value.

5.5 DISPLAY MESSAGES

Fixed messages inform the operator that elaborations are being performed. No action is needed.

Flashing messages indicate an anomalous condition or the request of operator intervention.

MESSAGES	MEANING
Fx-x	Firmware release indication. The letters "x" are actually appropriate numbers.
rx-x	As Fx-x but referred to a firmware located in ROM (masked).
----	The reading is being set to zero. Pressing key TARE activates the message.
Er18	Unstable reading while storing the tare. Press TARE to cancel. The previous tare value is retained. If this message appears at the switching ON of the instrument, it cannot be cancelled. Unplug the instrument from the mains for a few seconds and plug it in again: if the error reappears, the instrument is faulty or badly installed (e.g. on an unstable support).
Er19	Unstable reading during a calibration operation. Press key TARE to cancel. The previous calibration value is retained.
Er23	Pan not completely unloaded at the switching ON of the balance. Press key TARE to cancel. The weight present on the pan is taken as initial tare. It can be excluded (see P3 parameter).
ErrC	Calibration error. Wrong calibrating weight or wrong procedure. Press key TARE to cancel. The previous calibration value is retained.
ErrP	Pre-tare error. The sum of the existing tare value and the given value of the pretare exceeds the full scale capacity. Press TARE to clear.
-HI-	Overload condition. Remove the excess weight from the pan.
-LO-	The weight is less than the minimum allowed. (Pan in the wrong position or not present). Remove the cause to cancel.
PrE	Indication for the beginning of the introduction of the pretare via small keyboard (It stays on for about 2 sec.).
ACAL	Start of the calibration procedure. It stays on till the existing weight is taken as a tare.
LOAD	Request to place the calibration weight on the pan.
CALC	During the time needed to take the calibration weight and to calculate the calibration factor. At the end of this operation the calibrating weight value is displayed or the message ErrC is displayed (see above).
Px=yy	Function for the modification of the functional parameters of the instrument. The letters "x" & "y" represent respectively the parameter number and the current value of the parameter (see Table of point 5). The access to the function is obtained by pressing the key MODE repeatedly till the message Par is displayed and by pressing the key ON/OFF to confirm. Press the key MODE to pass to the following parameter. Press the key ON/OFF to modify the parameter value. Press the key TARE to go back to the operative mode.
EEPI	Initialization cycle of the non volatile internal memory (EEPROM). It may appear but only the first time that the electronic card is switched ON. If this message appears during normal use, this is an indication of a malfunction.
EEPEX	Reading error of the "x" parameter. It is an indication of a malfunction of the EEPROM. The message stays on for about ten seconds. The instrument will operate with the parameter standard values of "Table " (see point 5).
Pnor	All the parameter initial values have been restored. It is activated when the instrument has been unplugged from the mains and plugged in again after a few seconds, while pressing the key ON/OFF. (It stays on for about one second).
bIL	Return to the balance operative mode. It is activated by pressing the key MODE and confirmed with the key ON/OFF .
CAL	Calibration procedure. It is activated by pressing the key MODE and confirmed with the key ON/OFF .
PAr	Menu for the parameter selection. It is activated by pressing the key MODE and confirmed with the key ON/OFF .
CPZ	Piece counting operation. It is activated by pressing the key MODE and confirmed by pressing the key ON/OFF .
P10 ~ P100	Sample of 10~100 pieces in the CPZ mode. It is flashing when it indicates that the operator must place the piece sample on the pan and press the key ON/OFF to confirm.
Error	Only in the CPZ mode. The sample weighs less than one division of the balance scale. Press key TARE to cancel.
Pcalc	Only in the CPZ mode. The instrument is taking the weight of the sample and calculating the weight/pieces conversion factor.

RESTART (ON REQUEST)

How to active the “restart” tare memory function

Press the “MODE” key till the display shows ZERO.
 Press the ON/OFF key till the display shows **On** (Enable) or **OFF** (Disable).
 Confirm your choice by pressing the “TARE” key.
 The function ON or OFF will be stored.

6. COMMUNICATION TO HOST COMPUTER

6.1. GENERAL CONSIDERATIONS

The communication can be performed in a mono-directional way i.e. from the instrument to the host computer, or in a bi-directional way i.e. to and from host computer.

In the mono-directional mode (from instrument to host computer) it is possible to select the data transmission in the **continuous** mode or in the **on request mode**. In the first case the following data are sent to the output: each weight figure as soon as available, the messages related to the tare and calibration function when they are being activated, messages of any anomalous conditions such as over range or others.

When the transmission is bi-directional, a weight transmission command or a function command (e.g. Tare or Calibration) may be sent from the host computer. The instrument will send the weight current value if the command is a weight request, or an appropriate replay message after a function command. The **continuous** transmission mode, described above, may be activated or deactivated with an appropriate command from the host computer.

6.2. FORMATS

The selectable communication parameters are as follows:

- Baud rate from 75 to 9600 bit/sec. non split (the transmission and reception rates must be the same).
- Parity EVEN, ODD, or NONE, word length 7 or 8 bits, 1 or 2 stop bit.
- Continuous option or on request.

See details about the available options in tables **Parameter 0** and **Parameter 1** (point 8).

The data at the output of the instrument is always a sequence of 8 **ASCII** coded characters followed by CR (Carriage Return = **0D** hex) and LF (Line Feed = **0A** hex) for a total of **10 characters**.

The format of the output weight data is as follows:

stffffff + CR + LF	
s =	sign "-" or code space
t =	most significant figure of the weight data or "*" character indicating an unsteady weight value when "stability indicator RS232" is active (see P3 parameter)
f =	number or code space in place of each non significant zero and the decimal point in the appropriate position.

The message format, on the other hand, is as follows:

mmmmmmmm + CR + LF

m= character of the message or code space.

The command from host computer is a sequence of one or several **ASCII** coded **capital** characters, followed by CR. The maximum number of characters accepted before the command is executed or rejected is ten, even if the character CR has not yet been sent.

The first letter of the sequence identifies the command; if it is not one of those listed below, the command is void and no return message is sent.

If the command syntax is correct and the possible parameters contained in the command (see "**M**", "**P**") are within the limits of the foreseen values, the instrument sends the message **ACK**; otherwise it sends the message **NACK**.

6.3. COMMAND OF TICKET PRINT

The ticket print gets started by the serial command L (see RS232 command table) as it follows:

```
TARA.: 19.81
NETTO.: 36.78
P.M.U.: 0.0272
Pezzi.: 1352
```

if the selected language is Italian (see P3 parameter), while if it is selected in English, it will be:

```
TARE.: 19.81
NET...: 36.78
W.A.U.: 0.0272
PCs...: 1352
```

The first two lines are always printed while the other two get printed only if a sampling operation has been done and CPZ option selected by the four-key keyboard.

TARE is the memorized weight value shown on the display before any TARE command: two consecutive tare commands should be performed for zerosetting.

W.A.U. is the average unit weight value of the pieces used in the sampling operation and it is calculated with one or two decimal figures more than those available to the balance.

Print commands

The print gets started in two different ways:

1. by serial command suitably introduced ("**L**" command)
2. by pedal or by the "**print**" key on the keyboard (when available in our instruments)
 - The first option is possible with balance connected to computer or to one of our printers supplied with **EE16K** option (code **30.10.066** and **30.10.067**). In this case any information required by the Customer (logo, address etc.) can be printed on the ticket head.
 - The second option simply allows the use of any printer with a serial interface. Also in this case it is possible to print on the ticket head any Customer's information through data stored in the balance rather than in the printer. Infact, 80 locations of the EEPROM are dedicated to this purpose so that whatever the information the Customer requires, it can be written and printed before numerical values.

There is the possibility to obtain via serial interface all values of **Tare**, **Net Weight**, **Average Unit Weight** and **Pieces** number by adding to the "**L**" command respectively "**T**", "**N**", "**M**" and "**P**".

6.4. SERIAL COMMANDS FOR THE USE OF BALANCE AS A TERMINAL

By "D" and "Rdx" commands the message transmitted from the balance can be displayed instead of the weight value or even the balance can be switched off, though keeping its functions unchanged. Infact even if the display is switched to show the external coming message, the weight can be obtained using the proper command.

"RKx" command allows the inhibition of the keyboard functioning and the transmission via serial line of a code corresponding to the key pressed by the user.

So the balance can be used to obtain the weight, to process data in its proper application (for example for the calculation of the leftover stock in the warehose after goods discharging) and to transmit the result as a numerical message. Even though with a limited character-set, the transmission of alphanumerical messages can be used to guide the user during performances with complex operating sequences.

Transmission of messages to be displayed(command "D")

Messages of maximum six characters can be transmitted. A longer message is cut off on the right. They should be preceded by letter "D" or by signal "-" if a negative number is transmitted: in this case, "-" signal will light up on the display even if the message includes alphabetic characters.

It is possible to have one only character with decimal point lit up. If it is an EUROMATIC the decimal point must be necessarily referred to the first digit.

The displayable characters are those referring to "ASCII Table of codes recognized by EUROPE" and codes not admitted will be in blanks. The received message will be stored in memory and displayed only if external message display mode is on (see **Rdx** command).

7. TABLE OF THE SERIAL COMMANDS

7.1. TABLE OF AVAILABLE COMMANDS

Command	Function	Reply 1	Reply 2	Reply 3
B	Sending of a single data (2)	dato / data		
C	Calibration	ACAL	Err18 LOAD (4)	Err19 ERRORC ACK (4)
- <message>	Send message to display with "-" (see "D" command)	ACK NACK		
D <message>	Send message to display (see characters table)	ACK NACK		
F	Continuous transmission end	END		
I	Continuous transmission start	INIT		
M <number> (3)	Send pre-tare	ERRORP ACK NACK		
N	Instrument model request	EUxxxx		
PA=W	EEPROM writing of actual parameters value in EPROM	ACK NACK		
Px=?	Reading the actual value of parameter "x"	Px=y NACK		
Px=y	Set parameter "x" to the value "y"	ACK NACK		
RDx	Selection display mode (see table)	ACK NACK		
RKx	Selection display mode (see table)	ACK NACK		
T (1)	Tare	TARE Err18	ACK (4)	
L	Print CPZ ticket	(Vedi testo)		
LT	Send TARE value	dato / data		
LN	Send net weight value	dato / data		
LM	Send average unit weight	dato / data		
LP	Send pieces value	dato / data		

- (1) The reply messages are not sent if the command T is used to cancel error situations.
- (2) The command sequence "I" and "B" activates the "verbose" mode, i.e. it authorizes the transmission of all the messages associated to each command (see (4)). The command "F" deactivates the "verbose" mode.
- (3) The <pretare> value may be given in the format deemed to be the most convenient. In any case the total number of characters cannot be more than 9 and the integer part must be less than the end of scale value. The reply message "ErrP" means that the preset value for the pre-tare added to the actual tare value totals more than the end of scale value.
- (4) Messages sent only when the "verbose" mode is operative (See (2)).

7.2. ASCII TABLE OF RECOGNIZED CODES (EXADECIMAL VALUE)

	20	30	40	50	60	70
0	blank	0		P		
1		1	A			
2		2	B (1)		b	r
3		3	C	S	c	
4		4	D (1)		d	
5		5	E	U		u
6		6	F			
7		7				
8		8	H		h	° (2)
9		9	I	Y		
A			J	Z		
B						
C			L			
D	-	=				
E					n	
F			O		o	

- (1) Displayed with corresponding small characters.
- (2) "Degree" symbol: its standard ASCII code should be F8 instead of 78.

7.3. DISPLAYING INTERNAL/EXTERNAL AND FLASHING MESSAGE

Command	Display	Flashing normal	Message source
RD0	Off	-	-
RD1	On	Normal	Internal
RD2	On	Flashing	Internal
RD3	On	Normal	External
RD4	On	Flashing	External

RD0 command switches off the display but the balance functions are kept unchanged.

7.4. DISABLING KEYBOARD AND SEND KEYS TO SERIAL LINE

Command	Keyboard	Send in line
RK0	Disabled	Disabled
RK1	Enabled	Disabled
RK2	Disabled	Enabled

If mode "transmission on line" is active, when a key is pushed by the user, the answer will be:

Pressed key	Transmitted sequence
Tare	RKTARE
Mode	RKMODE
On-Off	RKONOFF
Print (pedal)	RKPRINT

8. WORKING PARAMETERS

Parameter 0 (baud rate, transmission type)		
Parameter value	Baud rate	Transmission type
0	9600	on request
1	4800	on request
2	2400	on request
3	1200	on request
4	600	on request
5	300	on request
6	150	on request
7	75	on request
8	9600	continuous
9	4800	continuous
10	2400	continuous
11	1200	continuous
12	600	continuous
13	300	continuous
14	150	continuous
15	75	continuous

Parameter 1 (parity, Bit/char, stop bit)				
Parameter value	Parity	Bit/char	Stop Bit	Sending of data (1)
0	Even	7	1	Current
1	Odd	7	1	Current
ND	No	7	1	Current
2	Even	7	2	Current
3	Odd	7	2	Current
4	No	7	2	Current
5	Even	8	1	Current
6	Odd	8	1	Current
7	No	8	1	Current
ND	Even	8	2	Current
ND	Odd	8	2	Current
8	No	8	2	Current
9	Even	7	1	Stable
10	Odd	7	1	Stable
ND	No	7	1	Stable
11	Even	7	2	Stable
12	Odd	7	2	Stable
13	No	7	2	Stable
ND	Even	8	1	Stable
ND	Odd	8	1	Stable
14	No	8	1	Stable
ND	Even	8	2	Stable
ND	Odd	8	2	Stable
15	No	8	2	Stable

N.B.	“ND” means unavailable configuration.
(1)	Mode “Current” enables data transmission as soon as it is available, while mode “Stable” enables data transmission only when when it is stable.

Parameter 2 (up to 65000 reading counts)		
Parameter value	Conversions/second	Seconds/conversion
0	14,36	-
1	7,18	-
2	3,59	-
3	1,80	0,55
4	0,90	1,11
5	0,45	2,22

Parameter 2 (more than 65000 reading counts)		
Parameter value	Conversions/second	Seconds/conversion
0	19,52	-
1	9,77	-
2	4,88	-
3	2,44	0,41
4	1,22	0,82
5	0,61	1,64

Parameter 3 (brightness, steadiness, language, Error 23)				
Parameter value	Display brightness	Steadiness indicator on RS232 (1)	Language (2)	Error 23 indication (3)
0	Normal	Disabled	Italian	Enable
1	½ intensity	Disabled	Italian	Enable
2	Normal	Enable	Italian	Enable
3	½ intensity	Enable	Italian	Enable
4	Normal	Disabled	English	Enable
5	½ intensity	Disabled	English	Enable
6	Normal	Enable	English	Enable
7	½ intensity	Enable	English	Enable
8	Normal	Disabled	Italian	Disabled
9	½ intensity	Disabled	Italian	Disabled
10	Normal	Enable	Italian	Disabled
11	½ intensity	Enable	Italian	Disabled
12	Normal	Disabled	English	Disabled
13	½ intensity	Disabled	English	Disabled
14	Normal	Enable	English	Disabled
15	½ intensity	Enable	English	Disabled

(1)	It consists in the introduction of stable data in the weight value transmitted on the serial interface: if datum is unstable, the most significant weight value is replaced by character “***”.
(2)	Languages of messages on CPZ ticket head.
(3)	Error occurring when the balance is switched on with a loaded pan.

Parameter <u>4</u> (display type)	
Parameter value	Display type
0	<u>Grams</u> (operative mode)
1	<u>Net linearized</u> counts
2	<u>Net not linearized</u> counts
3	<u>Millivolt</u> (on f.s. +/-2,5V)
4	<u>Gross linearized</u> counts
5	<u>Gross not linearized</u> counts


Parameter <u>5</u> (Average)	
Parameter value	Number of reading in the average
0	Disabled average
1	2
2	4
3	8
4	16
5	32

Parameter <u>6</u> (zero-tracking)	
Parameter value	Seconds
0	Disabled
1	1
2	2
3	3
4	4
5	5
6	6
7	8
8	10
9	12
10	14
11	16
12	18
13	20
14	25
15	30


9. INTERCONNECTIONS (ON REQUEST)

The interconnection is with optical fiber cable as the type shown in the figure F.

Pay attention to the rule of the connection at the end of the cable: the polarizing must be upper.


PAINTEX series

EUROPE series

READING UNIT
(for TMB and
OEM version)

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WARNING

- Pin 4 is connected internally to pin 5.
- Pin 8 is connected internally to pin 20.
- All signal must conform RS232-C specifications for proper operation.

INTERFACE MODULE
- UNTIL 60 meters -


Standard (Female)

TR60-DCE-F 9 pin

Pin 2 : Rx
Pin 3 : Tx
Pin 5 : GROUND

WARNING

- Pin 1 is connected internally to pin 4 e 6.
- Pin 7 is connected internally to pin 8.
- All signal must conform RS232-C specifications for proper operation.

Optical fiber interconnection cable
Length on request (Max 100 meters)

FIG. F - Serial connection in optical fiber

customer care

info@gibertini.com
Tel. (++39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (++39) 02 3541434
Fax (++39) 02 3541438
www.gibertini.com
sales@gibertini.com