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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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*All rights reserved.
Forbidden **any form** of print, duplication, reproduction and publication
of this manual, or part of it, without the written agreement of
GIBERTINI ELETTRONICA Srl*

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**ATTENTION !
READ CAREFULLY THIS MANUAL, WHICH IS INTEGRAL PART OF THE
INSTRUMENT, BEFORE INSTALLING THE BALANCE**

1. IMPORTANT WARNINGS

N.B. - These "WARNINGS" are integral part of the instrument

This instrument is to be used strictly for scientific purposes and/or for internal factory control only.

NOT FOR USE IN RETAIL TRADES. THIS INSTRUMENT IS NOT TO BE USED IN AREAS WITH DANGER OF EXPLOSION.

It is important to note the following points for correct operation of the instrument:

1. Read the technical instructions given in the operating manual carefully. This manual contains information about installation and use. If these instructions are not carried out precisely, the accuracy of results may be affected. If in doubt, please consult the supplier's engineers.
2. In common with any electrical equipment, installation should be carried out by a competent person. In particular:
 - ensure that the instrument is correctly earthed;
 - do not install the instrument in areas of high fire risk, for example, in the presence of inflammable gases and vapours;
 - do not touch the instrument with wet hands.
 - switch off and disconnect the power supply before removing any cover;
3. Install the instrument on a vibration free base and away from draughts or sources of heat and cold. Have the instrument checked by qualified staff with test weighing after installation or after repositioning.
4. Avoid connection to a power supply showing high voltage variations due to other loads. The instrument should not be connected to the same electrical circuit as other high power consuming equipment as large fluctuations may cause some inaccuracy.
5. The instrument should be checked every 6 months for correct functioning, calibration and accuracy of the internal calibration mass by a suitable qualified person.
6. A routine calibration check should be carried out each day before use, using a standard mass. This mass should be independently checked monthly.
7. Check the zero function before each weighing operation.
8. If in doubt about any function please call a qualified person. Avoid interference by non-qualified persons.
9. If the instrument is used for weighing food or drink, ensure that other substances are not weighed on the unit, which might cause contamination.
10. If the instrument is used for weighing chemicals and other similar substances ensure that it is cleaned and any spillage removed to avoid potentially dangerous reactions.
11. Always clean the instrument after use to avoid subsequent inaccuracies and/or damage to the equipment.
12. Gibertini Elettronica s.r.l. does not accept any legal liability for damages or other consequences due to mishandling of the instrument or its accessories.

ADDITIONAL PROCEDURES FOR WEIGHING SUBSTANCE WITH GREAT ACCURACY

Gibertini Instruments are manufactured to be reliable in use with high accuracy. However, reduced accuracy is always possible due to ambient conditions, incorrect installation or improper use. Where greater accuracy is required, the following procedures should be observed:

- use modern instruments (manufactured within three years) which are fully maintained;
- use certified calibration masses before each weighing operation;
- make at least two measurements;
- consider the use of two instruments installed together and make weighings on both units, otherwise use calibrated masses with a certificate of calibration from a recognised institute.

**ALL ELECTRONIC BALANCES, WHICH ARE SENSITIVE TO CHANGES OF GRAVITY, MUST BE CHECKED
AND CALIBRATED IN THE PLACE WHERE THEY ARE TO BE USED (2014/31/EU).**

**ALL OUR INSTRUMENTS MEET INTERNATIONAL EMC-EMI-RFI STANDARDS ACCORDING TO 2014/30/EU
DIRECTIVE ON ELECTROMAGNETIC COMPATIBILITY.**

2. TECHNICAL CHARACTERISTICS AND GENERAL DESCRIPTION

MODELLO	7500 DR
Capacity g	950 / 7500
Readability g	0,01 / 0.1
Linearity g	$\pm 0,05 / \pm 0,1$
Repeatability g	$\pm 0,05 / \pm 0,1$
Response time	1 sec.
Pan dimens.	$\varnothing 190$
Warm-up time	≥ 15 minutes
Calibration	With external mass (optional) 1000 g F1 class or multiples
Operatine temp.	$+10^{\circ}\text{C} / +40^{\circ}\text{C}$
Data output	USB 2.0 – WIFI 802/11
Operatine voltage	100 \pm 240 V - 50Hz with external power supply
Power consumption	< 60 mA
Standard equipment	Power supply - User manual - CD with driver

In the front of the balance (Fig. 1) there is the control panel, which is constituted by a LCD display and four command keys. The function of the keys is described in the correspondent section.


Fig. 1 – Control panel

In the rear (Fig. 2) the balance presents the connector for the power supply, the connector for the USB 2.0 interface and the ON-OFF power switch (only for "BP" version).


Fig. 2 – Rear view

3. INSTALLATION

For an optimum installation of the balance respect the following indications:

- **Place the balance on a rigid surface, away from vibrations, heat sources and air currents.**
- The operating voltage must comply with the declared Technical Characteristics.
- Avoid connecting the balance to an electric network subject to great tension swings.
- Connect the power supply, **supplied in equipment**, to the balance and then to a 100÷240V network socket.

- NOTE -

The balance can also be powered via the USB connection.

It is very important to use the power supply provided with the scale if WIFI transmission is activated.

After the initial test the display will show "**0.0 g**". Now the balance is in **operating mode**. At the first ignition after the installation, or if the balance was switched-off for a long time, wait for 15 minutes minimum (preheating) before using the instrument so that will be thermally stabilised.

- NOTE -

All electronic balances, which are sensitive to changes of gravity, must be checked and calibrated in the place where they are to be used.

To make a calibration a **class F1 certified mass of 1000g (or multiples)** is needed. Please follow the instructions described at Section 6.4

4. DISPLAY AND COMMAND KEYS

4.1 DISPLAY

When the balance is in operating mode, the display can show characters, symbols and icons, like those of the left side figure. Their meaning is described in the following.

Seven segments **alphanumeric characters**, for reading data and visualising the messages

Symbols	Meaning
	Weight/Data stable
%	Percentage
	Bar-graph: percent used capacity

Units	Meaning
<i>Pcs</i>	Pieces
<i>ct</i>	Carats (1 ct = 0,2 g)
<i>lb</i>	Pounds (1 lb = 453,59237 g)
<i>oz</i>	Ounces (1 oz = 28,34952 g)
<i>g</i>	Grams (Default)

4.2 KEYS FUNCTION

Everyone of the 4 command keys (Fig. 3a and 3b) has one or more functions assigned that vary in base of the mode in use.


Fig. 3a – Standard command keys


Fig. 3b – DR model command keys

When we are in **operating mode** the functions of the keys are the following:

Key	Function
RANGE	In the first 950g of the capacity, if pressed for less than one second, it changes the balance division.
MODE	Pressed for less than one second, it selects the second unit, if enabled. If pressed for more than one second, the balance get into a configuration mode (see Sec. 6).
ZERO-T	It zeroes the value on the display
ON-OFF	If pressed for less than one second, it switches off the display illumination. If pressed longer, it puts the balance in stand-by mode (see Sec. 5.4).

When we are in **configuration mode** instead (see Sec. 6) the functions change as follows:

Key	Function
RANGE	No function.
MODE	Pressed for less than one second it scrolls the possible selections.
ZERO-T	Press for less than one second to exit from the selection function.
ON-OFF	Pressed for less than one second it confirms the showed selection.

ATTENTION !

If the balance stops working because of incorrect operations on the keyboard, operate as follows:

- a- unplug the power supply unit from the mains socket and wait some seconds
- b- plug in again the power supply unit and simultaneously press the **ON-OFF** key. Release the key when the display shows "**EU-C**".

The message "**SEt Factory**" will confirm the reset of default parameters (factory settings).

5. OPERATING MODE

ATTENTION - When using the balance, avoid loading on the pan weights heavier than the balance capacity. At any rate, the objects to be weighed must be placed on the pan with care.

5.1 WEIGHING MODE

After possible warm-up, the balance is ready to weigh with precision. After placing the object or the objects to be weighed on the pan, the display will show their weight. By pressing the **ZERO-T** key, the weight is set to zero, and the display will show "**0.0**". This operation is useful for dosing or when it is necessary to perform the tare of the containers necessary to weigh bulk objects, powders, liquids, little animals, etc.

- NOTE -

If there is an interruption in the mains supply or the power supply is accidentally disconnected during the weighing operations, the balance stores in memory the data displayed at the moment of the interruption.

For example: if the display was showing "3950 g" at the moment of the power interruption, the display will show "3950 g" when the power supply is restored. If the pan is unloaded before restoring the power supply, when the balance turns on again the display will show "0.0 g" or the negative value of the possible tares set.

5.4 ENERGY SAVING OPTIONS

The available options for energy saving are:

- Option 1:** Pressing the **ON-OFF** key for less than one second, the backlight of the display switches off. The balance remains fully operating.
- Option 2:** Pressing the **ON-OFF** key for more than one second, the balance switches to "stand-by" mode and the display will show "**OFF**". Press again the **ON-OFF** key to switch-on the balance.

When the balance is not used for a very long time power it off by unplugging the power supply unit from the mains.

At the next power up, wait at least 15 minutes (warm-up) before using the instrument.

6. CONFIGURATION

6.1 "WIFI" MENU

Allows you to activate / deactivate the WIFI transmission (**Default UIFI OFF**).

When the display show "UIFI" press the "ON-OFF" key to confirm.

The display will show "UIFI OFF", press the "MODE" key to show the "UIFI ON" and confirm by pressing the "ON/OFF" key.

By activating the WIFI transmission (UIFI On) the USB serial transmission is automatically deactivated and the current IP address will be show on the scale display, the default address is "**192.168.4.1**".


The scale operates as an Access Point and the menu items "**UIFI ADDRESS**" and "**UIFI SSID**" become active.

UIFI	→ UIFI OFF	→ Turn off the WIFI transmission
	→ UIFI ON	→ Activate the WIFI transmission
	→ UIFI ADDRESS	→ Allows to view the set IP address.
	→ UIFI SSID	→ Allows to view the identification name of the WIFI network.

6.2 NETWORK CONFIGURATION

For the configuration of the network parameters it is not necessary to use any dedicated software, the whole procedure can be performed with a normal internet browser (Firefox, Chrome, Edge etc.) and from any device and / or operating system, simply connect to the configuration page with this procedure:

- Open the search for wireless networks, the factory configuration of the wireless network is unique for each instrument and is composed as follows:
 - **SSID** = EU-C_XXXXXXXXXX (where XXXXXXXXXXXX correspond to the MACADDRESS of the device)
 - **PASS** = EU-C-PASS
- in the example image, the network has SSID EU-C_3C71BF2F9CAC:

- Once connected to the network by accessing with a browser at 192.168.4.1, a summary page of the current configuration will be displayed:

- By accessing the setup page, it will be possible to configure the instrument in Access Point or Station mode.

6.3 "CALIBRATION" FUNCTION

The function CALIBRATION is the last option in the main menu of the configuration mode and is used to calibrate the balance.

To make a calibration a **class F1 certified mass of 1000g (or multiples)** is necessary. Operate as follow:

- When the display shows "**CALIBRATION**", press the **ON-OFF** key to confirm.
- The display will show "**LOAD**"; put on the pan the calibration mass.
- When the display shows "**UNLOAD**", remove the mass from the pan.
- Once the calibration process is completed, the display will show for some seconds "**CALIBRATED**" and then "**0.0**". The balance is calibrated.

If, for some reason, the calibration failed, the display will show "**CAL ERROR**". In this case, the previous calibration value will be maintained.

ATTENTION !

- Calibrate periodically the balance, especially when the balance is used intensively
- The balance must be calibrated with the pan absolutely unloaded.
- To obtain a precise calibration, the balance must be powered-on minimum 15 minutes earlier, and must be thermally stabilised at any rate.
- Use only class F1 masses, best if certified.

7. SERIAL TRANSMISSION

Serial transmission is a USB 2.0 connection.
The answer from the balance is sent every 200 ms.

7.1 USB CONNECTION

it is necessary to install the drivers supplied with the instrument and use the supplied cable.

7.2 SERIAL COMMANDS

In the "CRYSTAL" and "EURO" protocols the balance responds to the following commands:

COMMAND	FUNCTION
B	Send data on display
C	Start the calibration process
I	Start continuous transmission of data
F	Stop continuous transmission of data
T	Tare command

Notes : Every remote control command must end with "CR".
Line-feed characters are ignored.

8. MAINTENANCE

Always keep the pan clean. To clean the balance, do not use solvents or acid or corrosive substances. Use a soft cloth, made humid, if necessary, with water or non-aggressive cleansers.

Calibrate periodically the balance with certified masses, especially when the balance is used intensively, or after long inactive periods. Use only a class F1 certified mass of 1000g (or multiples).

9. DISPOSAL – INFORMATION FOR USERS

According to the 2002/95/CE, 2002/96/CE and 2003/108/CE Directives, concerning the reduction in the use of hazardous substances in electrical and electronic apparatus, as well as the disposal of waste materials

The symbol of a crossed box applied on the apparatus or on the packaging indicates that the product must be collected separately from other waste materials at the end of its useful life.

The separate waste collection of the apparatus which has reached the end of its useful life is organised and managed by the producer. The user who desires to get rid of present apparatus must therefore contact the seller and follow the given instructions.

Suitable separate waste collection for future sending of the disused apparatus for recycling, treatment and environmentally friendly disposal, contributes towards preventing any possible negative effects on the environment and on health and encourages the reuse and recycling of the materials the apparatus is made of.

Unauthorised disposal of the product by the user will lead to payment of the administrative sanctions in force in the country were it is put on the market.

customer care

info@gibertini.com

Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434

Fax (+39) 02 3541438

www.gibertini.com

sales@gibertini.com