ISTRUZIONI D'USO US
ER MANUAL MODE D'E
MPLOI ISTRUZIONI D'
USO USER MANUAL M
ODE D'EMPLOI ISTRU
ZIONI D'USO USER MA
NUAL MODE D'EMPLOI
ISTRUZIONI D'USO US
ER MANUAL MODE D'E

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INDEX

1.	IMPORTANT WARNINGS.....	20
2.	SPECIFICATIONS	21
1.1	General description	21
2.2	Comparison with the official methods	21
2.3	Electrical Requirement of DEE and of VADE-4.....	22
2.4	Water Requirements	22
2.5	Dimensions and weights.....	22
3.	UNPACKING.....	22
3.1	Standard accessory kit.....	22
3.2	Description	23
4.	INSTALLATION.....	24
5.	ALCOHOL DISTILLATION (ASV determination).....	25
5.1	Distillation.....	25
5.2	Cleaning.....	27
6.	THE DISTILLATION PROCEDURE FOR SPIRITS	28
7.	VOLATILE ACIDITY WITH VADE-4	30
7.1	Installation.....	30
7.2	Preliminary operations.....	31
7.3	Determination of volatile acidity	32
7.4	Determination of sulphur dioxide (SO ₂).....	32
7.5	Comparison with other acidometers	32
7.6	Distillation for determination of ASV (Alcoholic Strength by Volume).....	32
8.	MAINTENANCE	33
8.1	Distillation chamber	33
8.2	Electrodes	33
8.3	Column.....	33
8.4	Glass nozzle.....	33
9.	DISPOSAL – INFORMATION FOR USERS.....	34

*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

*All rights reserved.
Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

1. IMPORTANT WARNINGS

N.B. - THESE “WARNINGS” ARE INTEGRAL PART OF THE INSTRUMENT

THIS INSTRUMENT IS TO BE USED STRICTLY FOR SCIENTIFIC PURPOSES AND/OR FOR INTERNAL FACTORY CONTROL ONLY. NOT FOR USE IN RETAIL TRADES.

THIS INSTRUMENT IS NOT TO BE USED IN AREAS EXISTING DANGER OF EXPLOSION.

ALL OUR INSTRUMENTS MEET INTERNATIONAL EMC-EMI-RFI STANDARDS ACCORDING TO 2004/108 EEC DIRECTIVE ON ELECTROMAGNETIC COMPATIBILITY.

It is important to note the following points for correct operation of the instrument:

- 1- Read the technical instructions given in the operating manual carefully. This manual contains information about installation and use. If these instructions are not carried out precisely, the accuracy of results may be affected. If in doubt, please consult the supplier's engineers.
- 2- In common with any electrical equipment, installation should be carried out by a competent person. In particular:
 - ensure that the instrument is correctly earthed;
 - do not install the instrument in areas of high fire risk, for example, in the presence of inflammable gases and vapours;
 - switch off and disconnect the power supply before removing any cover;
 - do not touch the instrument with wet hands.
- 3- Install the instrument on a vibration free base and away from draughts or sources of heat. Check the instrument by a suitable qualified person after installation or after repositioning.
- 4- Avoid connection to a power supply showing high voltage variations due to other loads. The instrument should not be connected to the same electrical circuit as other high power consuming equipment as large fluctuations may cause some inaccuracy.
- 5- The instrument should be checked every 6 months for correct functioning, calibration and accuracy, also of the internal calibration mass, by a suitable qualified person.
- 6- Each day, before use, a routine calibration check should be carried out using a reference standard, possibly certified. The reference standard should be checked at least once a month.
- 7- Check the zero function before each operation.
- 8- If in doubt about any function please call a qualified person. Avoid interference by non-qualified persons.
- 9- If the instrument is used for weighing or analyse alimentary substances, do not use it with other substances that can cause contamination.
- 10- If the instrument is used for weighing or analyse chemicals and other similar substances, ensure that it is always carefully cleaned to avoid potentially dangerous reactions.
- 11- Always clean the instrument after use to avoid subsequent inaccuracies and/or damage to the equipment.
- 12- Gibertini Elettronica s.r.l. does not accept any legal liability for damages or other consequences due to mishandling of the instrument or its accessories.

ADDITIONAL PROCEDURES FOR WEIGHING SUBSTANCE WITH GREAT ACCURACY

Gibertini Instruments are manufactured to be reliable in use with high accuracy. However, reduced accuracy is always possible due to ambient conditions, incorrect installation or improper use. Where greater accuracy is required, the following procedures should be observed:

- use modern instruments (manufactured within three years) which are fully maintained;
- check the instrument calibration before any operation;
- make at least two measurements;
- consider the use of two instruments installed together and make measurements on both units, otherwise use reference standard certified from a recognised institute.

2. SPECIFICATIONS

1.1 General description

An automatic distillation unit designed for distillation of wines, spirits, perfumes and beers for determination of alcohol content and, with VADE-4, volatile acidity.

The measured quantity of sample, (100 or 200 ml for alcohol content, 20 ml for volatile acidity), is placed in the distillation chamber, reagents are added to neutralise the sample and to make a conductive solution and then distillation can begin.

The distillate is collected in a volumetric flask mounted on a balance arm which will stop the distillation when a preset volume has been collected.

For volatile acidity, the distillate is collected in a 300 ml conical flask.

Distillation Unit can be used for other parameters (ex : sorbique acid, total SO₂)

Further information on the book (Analysis methods of wines and alcoholic beverages – Gibertini current edition).

Distillation Unit is equipped with ecological apparatus that reduce consumption of water-cooling : water pass into the refrigerator only during the distillation.

DEE and VADE-4 are produced according to directives 2006/42 EEC and 85/375/EEC for the safety of the analysts.

The instruments are also in accordance with the directive 2004/108 EEC (electromagnetic compatibility).

2.2 Comparison with the official methods

Whether for ASV – Alcoholic Strength by Volume (DEE) or for volatile acidity (DEE with VADE-4), the instruments are in conformity with OIV methods as prescribed by Regulation n. 606/2009 EEC.

2.2.1 Alcoholic strength by volume (ASV)

For the wines, the current regulations permit the use of any model of the apparatus for distillation or extraction with a steam current of water on condition that it satisfies the following test:

- ❖ Distil 5 times in succession a 10 % vol water-alcohol mixture. After the last distillation, the distillate should present an Alcoholic Strength by Volume of almost 9,9 % vol, in other words, during each distillation it does not lose a quantity of alcohol higher than 0,02 % vol.

2.2.2 Volatile acidity

The ECC legislation does not described one unique type of distillation apparatus, but requires that any apparatus or technique which satisfies the following 3 requirements, constitutes an official international apparatus or technique.

1. The steam of water produced by the generator must be free from CO₂.
Distil 20 ml water collecting 250 ml of distillate, adding 0,1 ml di NaOH 0,1 N and 2 drops phenolphthalein 1 % in alcohol: the colour must remain stable pink for at least 10 seconds.
2. 99,5 % of acetic acid contained in a 0,1 M solution must be found in the distillate.
Distil 20 ml acetic acid 0,1 M.
Collect 250 ml distillate and titrate with at least 19,9 ml NaOH 0,1 N.
3. A maximum of 0,5 % of a 1 M solution of lactic acid may distil.
Distil 20 ml of a M solution of lactic acid, collect 250 ml of distillate and titrate with maximum 1,0 ml NaOH 0,1 N.

2.3 Electrical Requirement of DEE and of VADE-4

Power source : 220 V -10 % / +15 % , 50 Hz
Power consumption: 3300 W max (DEE)
2600 W (VADE-4)

The unit must be connected to a fused and switched mains supply which is efficiently earthed.

2.4 Water Requirements

Minimum pressure 2 bar
Flow rate 5 litres/minute
Temperature max 20°C

2.5 Dimensions and weights

DEE PV	L280 x P460 x H1000 mm	25 kg
VADE-4	L140 x P405 x H307 mm	10 kg

3. UNPACKING

If the unit has been shipped by carrier, carefully unpack the cartons and stand the distillation unit on a bench. Identify the box containing the reagents, bottles and various connectors.

A visual check should be made of all items for evidence of any damage which, if found, must be reported to the carrier and the suppliers within 3 days.

3.1 Standard accessory kit

3.1.1 Glassware

- 3 x 100 Volumetric flasks with plug
- 1 x 200 Volumetric flasks with plug
- 2 x 300 ml Conical flasks with plug
- 1 x 500 ml Dispenser
- 1 x 100 ml Dropper bottle
- 3 x Dropper plastic pipettes

3.1.2 Reagents

- Distilled water
- 1 x 500 ml Calcium oxide CaO suspension 12%
- 1 x 100 ml Anti-foam
- 1 x 500 ml Solution for liqueurs

3.1.3 Fitting

- 1 Length water tubing and clips
- 1 Length drainage tubing and clip

ATTENTION: The flasks are calibrated and interchangeable in mass, not in volume.

Remember that our flask explicitly do not indicate ml, but only 100 and 200.

We do not withdraw 100 ml, but a repeated quantity. So accuracy in ml is not necessary, but it is imperative that the volume of wine and distillate is the same at same temperature.

3.2 Description

- 1) Power cable 220 V + ground Fuse 16 A
- 2) Sprinkler for rinsing distillation chamber
- 3) Drain tap of distillation chamber
- 4) Electrodes
- 5) Hole for VADE-4 glass nozzle
- 6) Warning light for open drain tap
- 7) Sealing bung of distillation chamber
- 8) Hand knob for sealing bung
- 9) Plexiglas screen
- 10) START push button
- 11) START indicator
- 12) WATER indicator
- 13) VOLATILE analysis indicator
- 14) VOLATILE analysis push button
- 16) ON/OFF switch with indicator
- 18) Distillation chamber
- 19) Cylindrical calibrated masses
- 20) Distillate tube
- 21) Precision balance


Fig. 1

- 1) Power cable
- 22) Rear panel fastening
- 23) Rear panel
- 24) Pressure regulator-Flow regulator group
- 25) Water supply connector
- 26) Water discharge
- 27) VADE-4 connector


Fig. 2

4. INSTALLATION

Remove the screws (Fig. 2 - n. 22) of the rear panel of the distillation unit and remove any foam packing around the internal and external glassware.

Release the internal balance mechanism by fully unscrewing the red painted screw which locks the balance (Fig. 3).


Fig. 3

Connect the large bore corrugated drainage tube to the stainless steel water discharge pipe. Clamp the tube with a hose-clip. (Fig. 4)

Connect the mains water supply to the nipple on the Tee piece using the clear PVC. tubing and clamp with a hose-clip. (Fig. 4)


Fig. 4

Close the Distillation Unit with rear panel and connect the power cable to a 220 V -16A power supply protected by fuses.

Turn the ON/OFF switch on the Control panel and turn the water supply ON. The water indicator on the front panel (Fig. 5) will be illuminated. (Note: if the water indicator light is not on, the distillation cycle will not start).


Fig. 5

5. ALCOHOL DISTILLATION (ASV determination)

5.1 Distillation

Check that the drainage tap on the distillation chamber is closed and that the warning light is off (Fig. 6).


Fig. 6

Rinse out a clean 100 ml volumetric flask with the sample to be analysed and then accurately fill the flask with the wine or spirit to be tested. The meniscus must be on the mark (Fig. 7).

For beers and sparkling wines, a degassing stage must be carried out by filtration prior to filling the flask.


Fig. 7

Carefully pour the contents of the flask into the distillation chamber and then use the distilled water sprinkler to wash out the flask still overturned. Use approximately 20 ml of water keep rotating the flask (Fig.. 8).


Fig. 8

Place the flask on the balance arm with the cylindrical mass (100 ml) on the hexagonal pin (Fig. 9).

To avoid an loss of alcohol, it advised pour into the volumetric flask about 2 ml of cold distilled water.

MASS (100 ml)


Fig. 9

Add the required distillation reagents: for wine and beer distillation 5/10 ml of milk of lime (max 10 ml) and 5/10 drops of anti-foam.

The alkalisation of sample by the addition of suspension of milk of lime is necessary to avoid the distillation of the other volatile substances (acetic acid, sulphurous anhydride, aldehydes, etc.) which will interfere with the density.

In fact the desired goal is to obtain a distillate similar to a hydroalcoholic solution, that is to say a pure alcohol-water mixture.

To obtain a satisfactory distillation, optimum sample pH is about 7,5. To check if the addition of milk of lime is enough can be used litmus paper.

Anyway is visual a change of sample colour, this means that is change sample pH.

Lower the plexiglas screen and close the distillation chamber: insert the plug in the distillation chamber and screw in the knob clockwise until stop.

If this safety screen is not lowered, distillation will not start. (Fig. 10)

ATTENTION: it is necessary to screw in completely the plug knob to guarantee a perfect sealing of the distillation chamber.



Fig. 10

If the 200 ml volumetric flask is used, place on the hexagonal pin the correspondent mass (200 ml). Add a double quantity of milk of lime and anti-foam.

Press the "Start" switch (Fig. 1 - n. 10) whereupon the "Start indicator" (Fig. 1 - n. 11) should also illuminate. The distillation process will now start.

Fig. 11

The distillate (~ 80 ml) will collect in the 100 ml volumetric flask (or 160 ml in the 200 ml volumetric flask). The normal time of distillation is 3-4 minutes for 100 ml and 5-6 minutes for 200 ml.

The actual weight of sample to be collected can be altered by adjusting the position of the hexagonal pin on the balance arm.

For normal use 80 ml should be collected but it must never exceed 100 ml.

When the set quantity is collected the unit will turn off automatically.

Wait some minutes because the distillate temperature must be the same of the sample before distillation. Between sample temperature before distillation and distillate temperature after can be accepted maximum ± 2 °C.

Add distilled water exactly until the mark (Fig. 7), then seal with the rubber bung and shake the flask. This sample is then ready for analysis by densimetry.

5.2 Cleaning

Remove the sealing bung from the distillation chamber and raise the plexiglas safety screen.

Open the drain tap. When the residue of the sample has been drained, wash out the distillation chamber with the sprinkler and then, when all water has been drained out, close the drain tap (Fig. 12).

The unit is then ready for a new sample.


Fig. 12

6. THE DISTILLATION PROCEDURE FOR SPIRITS

Spirits present different characteristics from those of wines:

- Values of ASV up to 70%
- Variable composition

In fact spirits with very different composition exist: from cream whisky to fruit creams, from very sweet liqueurs such as vodka with fruit or limoncello to liqueurs with extracts of aromatic plants or roots, or even drinks with egg or coffee.

In addition the higher alcoholic strength result in:

- A very slow distillation in order to make sure of the complete extraction of alcohol.
- Collection of a greater volume of the distillate (around 95 ml).

It is difficult to find a common distillation method for all the different kinds of spirits because the various ingredients interfere with the distillation in different ways, sometimes blocking the passage of the alcohol. Below is listed the basic procedure for spirits which allows a good distillation to be carried out in a short time, without excessive foam, and resulting in a clear distillate; the specific methods applicable to different alcohol products are indicated in the subsequent paragraph.

<i>The distillation procedure</i>	<i>Observations</i>
1) The sample is poured up to the mark of the special flask ³	It is recommended that the temperature of the sample is 20°C. One can also operate at room temperature: from the sampling and the bringing back to volume after distillation a maximum deviation of ± 2°C is allowed.
2) Pour the contents into the bulb of the distillatory	
3) Wash the flask 3 times with 10 ml of distilled water and add washings to distillatory bulb	
4) Add the “ solution for liqueurs ” ⁴ : solution 1:1 of NaCl 1% and [KAl(SO ₄) ₂] 12%	The aim is to increase the conductivity of the sample and to suppress the aromatic substances.
5) Add a few ml of antifoam ²	Its addition avoids the formation of foam which would be an obstacle to the normal distillation of the alcohol.
6) Add about 2 ml cold distilled water to the same flask used for measurement of the sample. Collect ≈ 95 ml of distillate.	It is possible to regulate the quantity of the distillate by adjusting “the pivot of mass”, this balance system controls the automatic stop . (Fig. 9)
7) Add distilled water until near volume, stopper and bring to volume when the temperature reaches 20°C or the temperature of sample before distillation (see point 1 – observations)	The use of a thermostatic bath is advised
8) Stopper the flask and mix with care to avoid bubbles	

Regarding brandy, etc. the method is the same as that indicated for wines the exception of the addition of a small quantity milk of lime : because of the low acidity of such products, 1-2 ml is sufficient.

In the case of a higher than 50 % vol alcoholic strength it is advisable to dilute to 1:2.

³ The flask are calibrated and interchangeable **in mass**, not in volume.

Remember that our flask explicitly **do not indicate "ml"** but only 100-200.

We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary.

⁴ See “Example” table for the suggested quantity.

Example:

Fill a 50 ml volumetric flask (class “A”) and pour the sample into distillation chamber.

Wash internally the volumetric flask with distilled water and pour also the washing water into distillation chamber.

Add into the distillation chamber about 50 ml distilled water with the required distillation reagents (see “Examples” table).

The distillate must be collected into a 100 ml volumetric flask (class “A”).

Wait some minutes because the distillate temperature must be the same of the sample before distillation.

Add distilled water exactly until the mark.

The sample is ready for analysis by densimetry

Multiply the value % vol by 2

EXAMPLES

SAMPLE	INGREDIENTS	D.E.E PROCEDURE	%vol D.E.E.	Declared
Liqueur type Strega	Alcohol water, sugar, natural flavours, saffron	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	40.37	40
Lemon cream	Water, sugar, lemon infusion, milk, cream, alcohol, glucose syrup	100 ml +12 ml antifoam	22.05	22
Grappa of raspberries	grappa, raspberries juice, sugar, natural flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	32.02	32
Limoncello	Alcohol, sugar, infusion of lemon rind, natural flavours	100 ml + 3 ml antifoam +10 ml milk of lime	30.29	30
Whisky cream	Dairy cream, sugar, Whisky, water, natural flavours	100 ml +15 ml antifoam	14.71	15
Grappa	Grappa	100 ml +5 drop antifoam +3 ml solution for liqueurs	40.12	40
Pina Colada	sugar, alcohol, cream, rum, milk, flavours	100 ml + 9 ml antifoam	16.80	17
China Martini	Sugar, alcohol, flavours, colorants: caramel; acidifier: citric acid	100 ml +15 ml antifoam + 10 ml milk of lime	31.00	31
White Vermuth	Wine, sugar, alcohol, flavours	100 ml +1.5 ml antifoam + 5 ml milk of lime	16.06	16
Strong anisette	Alcohol, sugar, natural flavours	100 ml +1.5 ml antifoam + 5 ml solution for liqueurs	40.02	40
Liqueur with coffee	sugar, alcohol, coffee, natural flavours	100 ml + 12 ml antifoam	32.89	33
Cherry	Alcohol, sugar, egriot juice, natural flavours	100 ml +1.5 ml antifoam + 5 ml solution for liqueurs	30.03	30
Tequila	Distillate of agaves techilana	100 ml +1.5 ml antifoam + 2 ml solution for liqueurs	38.14	38

SAMPLE	INGREDIENTS	D.E.E PROCEDURE	%vol D.E.E.	Declared
Pastis	Alcohol, sugar, 2 ⁰ / ₀₀ anethole, infusion of vegetable substances, natural flavours	100 ml +1.5 ml antifoam + 5 ml solution for liqueurs	44.76	45
Banana Vodka	vodka, sugar, banana juice, flavours	100 ml +1.5 ml antifoam + 3 ml solution for liqueurs	25.17	25
Brandy	Alcohol hydrate, Brandy, natural flavours, colorant: caramel	100 ml + 3 ml antifoam + 2 ml solution for liqueurs	36.07	36
Liqueur with chocolate and mint	Milk, sugar, cream, glucose syrup, alcohol, milk protein, infusion of cacao, flavours, colorant: caramel	100 ml + 12 ml antifoam	16.06	16
Cachaca	Distillate of sugar-cane	100 ml + 3 ml antifoam + 3 ml solution for liqueurs	43.66	43
Amaretto (alcoholate for exclusive use in confectionery)	Alcohol, caramel, flavours	25 ml + 75 ml distilled water (collect 100 ml) + 3 ml antifoam + 3 ml solution for liqueurs	value % vol × 4 = 70.00	70
Sake	Distillate of rice	100 ml + 1 ml antifoam + 10 ml solution for liqueurs	15.85	16
Sambuca	Alcohol, sugars, anisette, natural flavours	100 ml + 5 ml antifoam + 8 ml solution for liqueurs	39.86	40
Cointreau	Alcohol, sugars, orange flavours	100 ml + 5 ml antifoam + 7 ml solution for liqueurs	39.81	40

The samples are only indicative, so they must not be considered exhaustive.

The differences found from the experimental data on some samples compared to the values indicated on the label fall within the tolerance expected from the current legislation (0,3 % vol on indicated ASV)

For other information to see the GIBERTINI book "Analysis methods of wines and alcoholic beverages" – current issue

7. VOLATILE ACIDITY WITH VADE-4

7.1 Installation

Place the VADE-4 at the left side of the DEE and insert the VADE-4 cable (Fig. 13 - n. 5) in the corresponding DEE connector.

Connect both instruments to two separate 220 V - 15 A sockets.

Insert the glass nozzle (Fig. 13 - n. 2) in the spout of the distillation chamber between the two electrodes. Screw firmly but not too tightly.

Insert the vent tube (Fig. 13 - n. 1) in a discharge or in a empty bottle (on the end of the distillation there is a steam refuge).

Insert the feeding tube "WATER LOADING" (Fig. 13 - n. 3) in a tank of bidistilled water **free from CO₂**.

Press the VOLATILE push button on the DEE panel.

- 1) Vent tube
- 2) Steam output
- 3) Water loading tube
- 4) Power cable
- 5) DEE cable connection
- 6) Reset buttons (on bottom)
- 7) Vaporizer warning light
- 8) Heater warning light
- 9) Switch OFF/ON
- 10) Fuse


Fig. 13

7.2 Preliminary operations

Switch ON the VADE by means of the ON/OFF switch (Fig. 13 - n. 9) and wait 30 seconds.

It is advisable to switch OFF the VADE-4 if this doesn't use for many hours.

Lower the plexiglas screen and close the distillation chamber with the sealing bung.

The distillation chamber must be empty and clean.

Press START and let the VADE-4 work for a few minutes in order to reach the steady state temperature.

This operation must be performed whenever more than 10 minutes elapse between one distillation and the following one in order to eliminate of the condensed water along the tubes.

WARNING! - The CO₂ concentration must be checked at least once a day and also whenever the water is changed or added to the feeding tank. Proceed as follows:

- place a conical flask on the DEE balance (without the additional mass). Pour into the distillation chamber 20 ml of bidistilled water, the same used to feed the VADE-4.
- Lower the plexiglas shield, close the distillation chamber and press the START push button.
- The distillation stops automatically as soon as about 250 ml of distilled liquid have been collected.
- Add to the distillate a few drops of phenolphthalein at 1 % dilution and 0.1 ml of NaOH 0.1 N.

The liquid must keep a stable pink colour for at least 10 sec.

If this condition is not met the water must be decarbonised and/or the column must be cleaned.

7.3 Determination of volatile acidity

Pour into the distillation chamber 20 ml wine free CO₂

It is necessary (for example if the wine has a fixed acidity less than 6 g/l) to add 1 ml tartaric acid 50 %. Check expiry date.

Lower the plexiglas shield, close the distillation chamber.

Insert a conical flask on the balance, without masse and press START push button .

The distillation stops automatically as soon as about 250 ml of distilled liquid have been collected. The time of distillation is about 6 minutes

At the end of the distillation, add to the distillate a few drops of phenolphthalein at 1 % and titrate with NaOH 0,1 N until the appearance of a pink colour stable for at least 10 seconds.

$$\text{Gross volatile acidity} = n \times 0,3 \quad (\text{were } n \text{ are the millilitres of NaOH used})$$

In order to obtain the correct value of volatile acidity, it is necessary to determine the amount of sulphur dioxide present in the distillate, which quantity must be subtracted from the gross volatile acidity.

7.4 Determination of sulphur dioxide (SO₂)

Add to already titrated distillate a few drops of sulphuric acid solution diluted 1:4 and about 5 ml of starch past 1%. The solution will become colourless.

Titrate the free sulphur dioxide with I₂ N/100 until the solution turns blue (millilitres used = n').

Now make the solution alkaline by adding 20 ml of saturated sodium borate solution. The solution will turn pink .

Titrate with I₂ N/100 until the appearance of the blue colour (millilitres used = n'').

Subtract from ml of NaOH 0,1 N (n) a quantity equivalent to $\frac{n'}{10} + \frac{n''}{20} = n'''$

$$\text{Correct volatile acidity} = (n - n''') \times 0,3 \quad (\text{g/l of acetic acid})$$

Using the automatic enological analyser QUICK, produced by Gibertini Elettronica, the titration of the distillate for the determination of volatile acidity, whether gross or correct, is automatic.

7.5 Comparison with other acidometers

In the event that the results obtained with this system of distillation are compared with other instruments it is absolutely essential to start from the same initial volume (20 ml) and to collect the same final volume (250 ml)

Otherwise the results are not comparable. It must also be kept in mind that the distillation with DEE and VADE-4 is more rapid than other and consequently more sulphur dioxide (SO₂) tends to pass into the distillate.

So there will be a larger quantity to subtract in order to obtain the correct volatile acidity.

7.6 Distillation for determination of ASV (Alcoholic Strength by Volume)

To pass to the distillation for the determination ASV, disconnect the "VOLATILE" key. Unscrew the glass nozzle from the distillation chamber and re-close with the black plug without to force.

WARNING! On the front panel of VADE-4 there are two warning lights with the danger signal (Fig. 13 - n.8 and 8).

When the VADE-4 is working the warning lights must be normally switch off.

In case of lighting it is necessary go on in this way :

Control the water of supply of the VADE-4.

The failing of water can cause a excessive overheating.

Push, with a point, the buttons placed under the instrument (Fig. 13 - n. 6).

Continue to use the VADE-4 pay attention that the warning lights don't switch on again.

If this problem repeat again, call the SAT (Service) Gibertini.

8. MAINTENANCE

8.1 Distillation chamber

After each distillation the distillation chamber must be washed with water to prevent contamination of the successive sample.

Periodically clean the distillation chamber introducing approximately 500 ml of distilled or demineralised water and approximately 20 ml of sodium chloride NaCl 1% solution. Remove the glass nozzle and start the distillation in alcohol (without start the function VOLATILE).

To remove deposits or incrustations that the samples can have left on the electrodes and/or on the glass can be done in the two following ways:

- 1- Pour into the distillation chamber a dilute solution of HCl 0,1N. Await some hour and eventually use a tube-brush to clean the chamber.
- 2- Introduce into the distillation chamber the solution B, supplied with the instrument, in sufficient quantity to cover the incrustations. Also in this case to await some hour and use eventually a tube-brush.

In any case rinse always accurately the distillation chamber before reuse the instrument.

WARNING! Don't press the START button when into the distillation chamber there is a cleaning solution. For more safety remove the power socket from the main.

8.2 Electrodes

The electrodes, with the use, stretch to be consumed and therefore must be replaced in order to maintain an effective distillation. It is difficult to predict the frequency of replacement because it will depend upon total usage and type of sample (approximately 2000 tests).

8.3 Column

During the normal use of the DEE, the column doesn't require a particular cleanness.

However, if the column is contaminated by acid residues (for example if the distillate in steam flow doesn't respect the requisites of the point 1 for the volatile acidity, as indicated in the page 3), it is necessary to clean it in this way :

- ❖ fill the distillation chamber with about 500 ml of distilled water and 20 m of sodium chloride NaCl 1 % solution. Start distillation in alcohol MODE (VOLATILE key disable).

It is not necessary to collect the distillate. Stop the distillation holding pressed toward the balance for some second .

8.4 Glass nozzle

If the glass nozzle is incrustated, it is possible to clean in this way :

- ❖ move back the black plug with hole and the gaskets to the rubber tube
- ❖ insert the glass part in a cylinder filled by cleaning solution for nozzle (solution B)
- ❖ leave plunged for some hour and successively wash the nozzle plunging it in water.

9. DISPOSAL – INFORMATION FOR USERS

According to the 2002/95/CE, 2002/96/CE and 2003/108/CE Directives, concerning the reduction in the use of hazardous substances in electrical and electronic apparatus, as well as the disposal of waste materials.

The symbol of a crossed box applied on the apparatus or on the packaging indicates that the product must be collected separately from other waste materials at the end of its useful life.

The separate waste collection of the apparatus which has reached the end of its useful life is organised and managed by the producer. The user who desires to get rid of present apparatus must therefore contact the seller and follow the given instructions.

Suitable separate waste collection for future sending of the disused apparatus for recycling, treatment and environmentally friendly disposal, contributes towards preventing any possible negative effects on the environment and on health and encourages the reuse and recycling of the materials the apparatus is made of.

Unauthorised disposal of the product by the user will lead to payment of the administrative sanctions in force in the country where it is put on the market.

HOW TO CLEAN THE DISTILLATION CHAMBER AND THE RECTIFICATION COLUMN

	CLEANING MODE	WHEN TO DO IT	PROCEDURE
1-	Standard cleaning	Periodically, after a series of VOLATILE distillations.	<p>a. Pour in the distillation chamber about 20ml of distilled water, and perform distillation in VOLATILE mode.</p> <p>b. Repeat once or twice.</p>
2-	Cleaning of rectification column 1st mode	When the column appears to be polluted by acid residuals (i.e. when the distillate of water requires more than 0.1ml of NaOH 0.1N to be titrated).	<p>a. Fill the distillation chamber with about 500 ml of distilled water + about 20 ml of NaCl 1% and perform a distillation in ALCOHOL mode.</p> <p>N.B.- <i>This type of washing performs also a mechanical action on the column, since the water goes up from the distillation chamber into the column.</i></p>
3-	Cleaning of rectification column 2nd mode	When the column filling is made black by the products of the distillations, or when acidity persists during volatile distillations (the distillate of water requires more than 0.1ml of NaOH 0.1N to be titrated).	<p>a. Pour in the distillation chamber about 1ml of a bleach (sodium hypochlorite), add about 500 ml of distilled water and perform a distillation in ALCOHOL mode.</p> <p>b. Repeat cleaning as indicated in point 2 (1st mode).</p> <p>c. Pour in the distillation chamber about 20 ml of distilled water, and perform distillation in VOLATILE mode.</p> <p>d. Repeat distillation of water until the water distillate requires less than 0.1ml of NaOH 0.1N to be titrated.</p>
4-	Cleaning of distillation chamber	When the distillation chamber appears dirty or there are residual inside.	<p>a. Pour in the distillation chamber about 100 ml of bleach, and let it rest for about 30 minutes. DON'T DISTILL !!!</p> <p>b. Drain and wash.</p> <p>c. Pour in the distillation chamber about 20 ml of distilled water, and perform distillation in VOLATILE mode.</p> <p>d. Repeat distillation of water until the water distillate requires less than 0.1ml of NaOH 0.1N to be titrated.</p>

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