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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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*All rights reserved.
Forbidden **any form** of print, duplication, reproduction and publication of
this instruction manual, or part of it, without the written agreement of
GIBERTINI ELETTRONICA Srl*

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
1. IMPORTANT SAFETY WARNINGS

CE  **II 3(3)G Ex nA[lc Gc] IIC T4 Gc**

This instrument is dedicated to the usage in environments with risk of explosion, classified as **ZONE 2** and it is built with **method of protection nA[lc Gc]** for **gas group IIC** and **temperature class T4 (135°C)**, according to what required by regulation 2014/34/EU.

- This manual is integral part of the balance. It lists all the rules to operate safely and provides important indications for usage and maintenance of the balance.
- The balance can be used only in the dangerous zone specified (**ZONE 2**). It must never be used in risk zone different than the one indicated.
- Safety of the balance is granted only if used according to the instruction of this manual, and according to what indicated by the label data on the keyboard panel.
- All the plaques and labels on the balance are integral part of the balance, and they must not be removed, hidden or modified in any way.
- For cleaning and maintenance operation, please refer to what specified on this manual.
- Exercise all the necessary precautionary measures to prevent creation or build up of electrostatic charges. Do not cover the balance with charged material, or with material that may build up electrostatic charges. **The operator and the service technicians must wear suitable working outfits and shoes too.**
- Every repair intervention must be performed with complying pieces, and it must be performed by qualified and authorized technicians.
- The RS232 serial interface of the balance has intrinsic safety. Possible peripherals connected to the balance must be equipped with protection devices at least equal to the balance (see technical data table).
- The responsible for the usage of the balance must grant that the product is used complying with the use it was built for and must also verify that the safety regulation in force in the country of usage are abided.
- The information written hereby are up-to-date to the day of printing. The manufacturer reserves the right to change them without notice. It is advisable to verify possible updates or revisions.

NOT OBSERVING THESE INSTRUCTIONS OR USAGE NOT COMPLYING TO WHAT WRITTEN IN THIS MANUAL LEADS TO MAKE THE WARRANTY VOID AND TO REVOCATE THE COMPLIANCE TO THE ATEX REGULATION.

2. IMPORTANT UTILIZATION WARNINGS

N.B. - These "WARNINGS" are integral part of the instrument

**This instrument is to be used strictly for scientific purposes and/or for internal factory control only.
NOT FOR USE IN RETAIL TRADES.**

It is important to note the following points for correct operation of the instrument:

1. Read the technical instructions given in the operating manual carefully. This manual contains information about installation and use. If these instructions are not carried out precisely, the accuracy of results may be affected. If in doubt, please consult the supplier's engineers.
2. For the connection to the electrical circuit and for the grounding observe scrupulously the regulation of the country of utilisation.
3. Before connecting the power plug to the mains, be sure that the power connector is connected to the balance, and that the blocking mechanism is locked.
4. Before initiating the balance, be sure that the instrument and the power supply are in perfect condition. No modifications are to be made on the instrument and/or on the power supply.
5. Any operation of repair or maintenance must be carried out from qualified and authorised staff of GIBERTINI ELETTRONICA s.r.l.
6. Take all the necessary precautionary measures for the installation and usage of any electrical instrument under tension. In particular:
 - do not touch the instrument with wet hands.
 - switch off and disconnect the power supply before removing any cover;
7. Install the instrument on a vibration free base and away from draughts or sources of heat and cold. Have the instrument checked by qualified staff with test weighing after installation or after repositioning.
8. Avoid connection to a power supply showing high voltage variations due to other loads. The instrument should not be connected to the same electrical circuit as other high power consuming equipment as large fluctuations may cause some inaccuracy.
9. The instrument should be checked every 6 months for correct functioning, calibration and accuracy of the internal calibration mass by a suitable qualified person.
10. A routine calibration check should be carried out each day before use, using a standard mass. This mass should be independently checked monthly.
11. Check the zero function before each weighing operation.
12. If in doubt about any function please call a qualified person. Avoid interference by non-qualified persons.
13. If the instrument is used for weighing food or drink, ensure that other substances are not weighed on the unit, which might cause contamination.
14. If the instrument is used for weighing chemicals and other similar substances ensure that it is cleaned and any spillage removed to avoid potentially dangerous reactions.
15. Always clean the instrument after use to avoid subsequent inaccuracies and/or damage to the equipment.
16. GIBERTINI ELETTRONICA s.r.l. does not accept any legal liability for damages or other consequences due to mishandling of the instrument or its accessories.

ADDITIONAL PROCEDURES FOR WEIGHING SUBSTANCE WITH GREAT ACCURACY

Gibertini Instruments are manufactured to be reliable in use with high accuracy. However, reduced accuracy is always possible due to ambient conditions, incorrect installation or improper use. Where greater accuracy is required, the following procedures should be observed:

- use modern instruments (manufactured within three years) which are fully maintained;
- use certified calibration masses before each weighing operation;
- make at least two measurements;
- consider the use of two instruments installed together and make weighings on both units, otherwise use calibrated masses with a certificate of calibration from a recognized institute.

ALL ELECTRONIC BALANCES, WHICH ARE SENSITIVE TO CHANGES OF GRAVITY, MUST BE CHECKED AND CALIBRATED IN THE PLACE WHERE THEY ARE TO BE USED (2014/31/EU DIRECTIVE).

ALL OUR INSTRUMENTS MEET INTERNATIONAL EMC-EMI-RFI STANDARDS ACCORDING TO 2014/30/EU DIRECTIVE ON ELECTROMAGNETIC COMPATIBILITY.

3. TECHNICAL CHARACTERISTICS AND GENERAL DESCRIPTION

3.1 TECHNICAL CHARACTERISTICS

BALANCE			
MODELS	EU-C7500PT	EU-C7500DR	
Capacity	g	7500	950 / 7500
Readability	g	0,1	0,01 / 0,1
Linearity	g	± 0,3	± 0,05 / ± 0,1
Repeatability	g	± 0,1	± 0,05 / ± 0,1
Response time		1 sec.	1 sec.
Pan dimensions	mm	Ø 190	
Warm-up time		≥ 15 minutes	
Calibration (external)		With 1000 g F1 class mass (optional) or multiples	
Operating temperature		+10°C / +40°C (see note)	
ATEX classification		II 3(3)G Ex nA[lc Gc] IIC T4 Gc	
Operating voltage		Vin max: 15 Vdc, lin Max: 80 mA	
Power consumption		Max 80 mA	
Dimensions and weight		210x355x120 mm - 4 kg ~	
RS232 I/O INTERFACE			
Connector		9 poles sub D female	
Connections		PIN OUT: Pin 2 = TX Pin 3 = RX Pin 5 = GND	
Safety parameters		Uo = 11 V; Io = 26 mA; Po = 75 mW Ui = 11 V; Ii = 26 mA Ci = 0; Li = 0	
EXTERNAL POWER SUPPLY			
Classification		Not classified – Use in safety zone	
Supply voltage		100÷240 Vac, 50/60 Hz, 180 mA	
Output voltage		Vout: 9 Vdc, Iout: 1200 mA	
Connection cable		3 metres	

Note - The operating temperature in the table (+10° C / +40° C) refers to the limits of the ambient temperature within which the correct metrological operation of the balance is guaranteed. Regarding security, the instrument is in compliance with the requirements of directive **2014/34/EU**.

3.2 GENERAL DESCRIPTION

The balance is supplied with power supply and instruction manual.

In the front of the balance there is the control panel (Fig. 1) which is constituted by a LCD display and four command keys. The function of the keys is described in the correspondent section.


Fig. 1 – Control panel



Fig. 2 – Rear view

- 1- Serial number
- 2- RS232 connection
- 3- AC power connection
(with safety retainer)
- 4- Ground screw

4. INSTALLATION

INSTALLATION SCHEME

- (1) Weighing unit
- (2) Power supply for safe zone
- (3) Peripherals connected to the balance must have the intrinsic safety protection "nL" in conformity with the installation requirement described in the present manual. The cable to connect the peripheral to the weighing unit must be 10 meters long maximum.

- ATTENTION -

- The connection cables must be protected from damaging and they do not have to be coiled or folded with sharp angles. If possible canalise them.
- The balance is not supplied with a power-on switch. To power-on the balance connect the power supply to the mains.

For an optimum installation of the balance respect the following indications:

- **Place the balance on a rigid surface, away from vibrations, heat sources and air currents, especially for the EU-C7500DR model.**
- The operating voltage must comply with the declared Technical Characteristics.
- Avoid connecting the balance to an electric network subject to great tension swings.
- Connect the power supply, **supplied in equipment**, to the balance and then to a 100÷240V network socket.

After the initial test the display will show "0.0 g". Now the balance is in **operating mode**. At the first ignition after the installation, or if the balance was switched-off for a long time, wait for 15 minutes minimum (preheating) before using the instrument so that will be thermally stabilized.

- NOTE -

All electronic balances, which are sensitive to changes of gravity, must be checked and calibrated in the place where they are to be used.

To make a calibration a **class F1 certified mass of 1000g (or multiples)** is needed. Please follow the instructions described at Section 7.3

5. DISPLAY AND COMMAND KEYS

5.1 DISPLAY

When the balance is in operating mode, the display can show characters, symbols and icons, like those of the left side figure. Their meaning is described in the following.

Seven segments **alphanumeric characters**, for reading data and visualizing the messages

Symbols	Meaning
	Bar-graph: percent used capacity

Units	Meaning
	Weight/Data stable
g	Grams (Default)

5.2 KEYS FUNCTION

Everyone of the 4 command keys (Fig. 3) has one or more functions assigned that vary in base of the mode in use.


Fig. 3 – Command keys

When we are in **operating mode** the functions of the keys are the following:

Key	Function
RANGE-PRINT	In the first 950g of the capacity, if pressed for less than one second, it changes the balance division (only DR model). If pressed for more than one second, it activates the PRINT function (if the printing protocol is active, see Sec. 7.2).
MODE	Pressed for less than one second, it selects the second unit, if enabled (see Sec. 7.1). If pressed for more than one second, the balance get into a configuration mode (see Sec. 7).
ZERO-T	It zeroes the value on the display
ON-OFF	If pressed for less than one second, it switches off the display illumination. If pressed longer, it puts the balance in stand-by mode (see Sec. 6.5).

When we are in **configuration mode** instead (see Sec. 7) the functions change as follows:

Key	Function
RANGE-PRINT	No function.
MODE	Pressed for less than one second it scrolls the possible selections.
ZERO-T	Press for less than one second to exit from the selection function.
ON-OFF	Pressed for less than one second it confirms the showed selection.

ATTENTION !

If the balance stops working because of incorrect operations on the keyboard, operate as follows:

- a- unplug the power supply unit from the mains socket and wait some seconds
- b- plug in again the power supply unit and simultaneously press the **ON-OFF** key.
Release the key when the display shows "**EU-C**".

The message "**SEt Factory**" will confirm the reset of default parameters (factory settings).

6. OPERATING MODE

ATTENTION - When using the balance, avoid loading on the pan weights heavier than the balance capacity. At any rate, the objects to be weighed must be placed on the pan with care.

6.1 WEIGHING MODE

After possible warm-up, the balance is ready to weigh with precision. After placing the object or the objects to be weighed on the pan, the display will show their weight.

By pressing the **ZERO-T** key, the weight is set to zero, and the display will show "**0.0**". This operation is useful for dosing or when it is necessary to perform the tare of the containers necessary to weigh bulk objects, powders, liquids, little animals, etc.

- NOTE -

If there is an interruption in the mains supply or the power supply is accidentally disconnected during the weighing operations, the balance stores in memory the data displayed at the moment of the interruption.

For example: if the display was showing "3950 g" at the moment of the power interruption, the display will show "3950 g" when the power supply is restored. If the pan is unloaded before restoring the power supply, when the balance turns on again the display will show "0.0 g" or the negative value of the possible tares set.

6.2 ENERGY SAVING OPTIONS

The available options for energy saving are:

Option 1: Pressing the **ON-OFF** key for less than one second, the backlight of the display switches off. The balance remains fully operating.

Option 2: Pressing the **ON-OFF** key for more than one second, the balance switches to "stand-by" mode and the display will show "**OFF**". Press again the **ON-OFF** key to switch-on the balance.

Unplug the power supply unit from the mains when the balance is not used for a very long time. In this case, at the next power up, wait at least 15 minutes (warm-up) before using the instrument.

6.3 CALIBRATION

If the necessary calibration mass is available (1000 grams, or multiples, in class F1), a balance calibration is possible. **After at least 15 minutes of thermal stabilisation**, the calibration can be performed as follows:

- Press and hold the **MODE** button until the display shows "**CALIBRATION**"
- Press the **ON-OFF** key to confirm.
- The display will show "**LOAD**": put on the pan the calibration mass.
- When the display shows "**UNLOAD**", remove the mass from the pan.
- Once the calibration process is completed, the display will show for some seconds "**CALIBRATED**" and then "**0.0**". The balance is calibrated.

If, for some reason, the calibration failed, the display will show "**CAL ERROR**". In this case, the previous calibration value will be maintained.

ATTENTION !

- Calibrate periodically the balance, especially when the balance is used intensively
- The balance must be calibrated whit the pan absolutely unloaded.
- To obtain a precise calibration, the balance must be powered-on minimum 15 minutes earlier, and must be thermally stabilized at any rate.
- Use only class F1 masses, best if certified.

- NOTE -

The calibration process does not start if a weight of more than 2 grams is on the pan at calibration start time.

7. SERIAL TRANSMISSION

Serial transmission is compatible with EIA RS232 standard. Every answer from the balance is sent every 200 ms.

ATTENTION : The RS232 balance interface has an intrinsic safety protection device. The peripherals connected to the balance must have minimum an intrinsic safety protection equal to the balance device.

7.1 SERIAL OUTPUT CONNECTION

<i>BALANCE</i>		<i>DB9 CONNECTOR (9 poles)</i>	<i>DB25 CONNECTOR (25 poles)</i>
RS 232 Out	Pin 2	Pin 2	Pin 3
RS 232 In	Pin 3	Pin 3	Pin 2
GROUND	Pin 5	Pin 5	Pin 7

7.2 SERIAL COMMANDS

The balance responds to the following commands:

COMMAND	FUNCTION
B	Send data on display
C	Start the calibration process
I	Start continuous transmission of data
F	Stop continuous transmission of data
T	Tare command

Notes : Every remote control command must end with "**CR**".
Line-feed characters are ignored.

8. MAINTENANCE

The frequency of the cleaning operation depends on the use of the balance. It is however good practice, regarding the safety, to avoid accumulation of powder, dirt and/or liquids, **especially in case of dangerous, inflammables substances or that can cause sparks**. Keep always clean the pan and the body of the balance. To clean these parts do not use solvents or acid or corrosive substances. Use a soft cloth, made humid, if necessary, with water or non-aggressive cleansers; best if antistatic products are used.

Calibrate periodically the balance with certified masses, especially when the balance is used intensively, or after long inactive periods. Use only a class F1 certified mass of 1000g (or multiples).

9. DISPOSAL – INFORMATION FOR USERS

According to the 2011/65/CE and 2003/108/CE Directives, concerning the reduction in the use of hazardous substances in electrical and electronic apparatus, as well as the disposal of waste materials

The symbol of a crossed box applied on the apparatus or on the packaging indicates that the product must be collected separately from other waste materials at the end of its useful life.

The separate waste collection of the apparatus which has reached the end of its useful life is organised and managed by the producer. The user who desires to get rid of present apparatus must therefore contact the seller and follow the given instructions.

Suitable separate waste collection for future sending of the disused apparatus for recycling, treatment and environmentally friendly disposal, contributes towards preventing any possible negative effects on the environment and on health and encourages the reuse and recycling of the materials the apparatus is made of.

Unauthorised disposal of the product by the user will lead to payment of the administrative sanctions in force in the country were it is put on the market.

customer care

info@gibertini.com
Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434
Fax (+39) 02 3541438
www.gibertini.com
sales@gibertini.com