

GSA Grain Size Analyzer


Perform soil particle size analysis automatically! Up to 6 analyses at the same time in 5 *6 hours.



Compliant with Official Methods

ASTM D422 UNI CEN ISO/TS 17892-4





GRANULOMETRIC ANALYSIS

The particle size analysis is performed using two techniques:

- **sieving**, for the coarse fraction
- sedimentation, for the fine fraction

It allows us to know the content of primary particles in the soil, divided by size.

The primary solid component of the soil is characterized by particles of various sizes: from diameters of the order of centimeters to millimeters and tenths of microns for the finest components.

GRANULOMETRIC CURVE


From the analysis of a soil's particle size curve, it is possible to derive its type and degree of particle distribution.

The analysis results are more clearly represented using charts in ordinary or logarithmic scales that show:

- the **diameter** of the particles, on the x-axis
- the percentage by weight that has a specified diameter, on the y-axis

The particle size curve is a percentage curve and indicates the % by weight of the fraction that exceeds a given diameter.

Time	Reading	Temperature	Diameter	Percentag d
24	1,0078	21.9	0,00745	12.035
27	1,0075	21.9	0.00703	11.228
30	1,0073	21,9	0,00667	10,609
36	1,007	21.8	0,00609	9.51
42	1,0067	21,8	0,00564	8,593
48	1,0065	21.9	0.00527	7.84
54	1.0063	21.8	0.00487	7.21
60	1,0061	21.6	0.00472	6,711
72	1,0058	21,8	0,00431	5,729
84	1,0057	21,7	0,00399	5,233
96	1,0055	21.7	0,00373	4,603
108	1.0063	21,7	0.00352	4.13
120	1,0053	21,7	0,00334	3,997
144	1,0051	21.7	0.00305	3,305
168	1,005	21.6	0,00283	2,865
192	1,0049	21.6	0.00264	2,472
216	1,0048	21,5	0,0025	2,127
240	1.0047	21.5	0,00237	1,916
288	1,0046	21.4	0.00217	1,509
336	1,0045	21,3	0,00201	1,2
384	1,0045	21.2	0,00168	0,965
432	1.0045	21.1	0.00177	0,852
460	1.0045	21	0,00168	0.744



HOW DOES IT WORK?

GSA measures the **density of the suspension** at predefined time intervals. The result is expressed in g/kg. The water-soil suspension is prepared inside a glass cylinder containing water, soil, and sodium hexametaphosphate to **facilitate the dispersion of particles**. The sample is **automatically shaken for 10 minutes** at a controlled and constant speed. The density is read using a float attached to the **hydrostatic balance**.

For the determination of classes or fractions, the software applies **Stokes' law**, which establishes a relationship between the **diameter** of suspended granules, the **viscosity** of the liquid, the **falling speed**, and the **specific weight** of the granules. Particles with a larger diameter (sand) settle first, while finer particles (silt and clay) remain in suspension for a long time.

The **results of the particle size analysis**, divided by sieve diameter, are entered into a table, and subsequently, for each particle size class (coarse sand, fine sand, coarse silt, fine silt, and clay), a **statistical analysis** of the data is performed.

Through the data provided by the statistics, it is possible to give an indication of the **homogeneity of the soil** sample used in the texture analyses.

Table of the particle size classes of interest, divided based on the diameter intervals of the particles*

Coarse sandCoarse siltClay> 100 μm50 – 20 μm< 2 μm</td>Fine sand
100 – 50 μmFine silt
20 – 2 μm* It is possible to modify and
customize the table of intervals

MAIN FEATURES

- **Detection** of soil from 0.1 mm to 0.001 mm
- Measurement of density variation by sedimentation
- Complete analyses in 5 to 6 hours
- 1, 3, or 6 samples **simultaneously**
- Density range read from 0.900 to 1.0500
- Automatic compensation for temperature variation and Stokes' law
- Repeatability better than 2%
- · Data saved in PDF or CSV file

GSA stops automatically at the end of the analysis, without the need for operator intervention





FASTER AND SAFER RESULTS

- Complete analyses in 5 *6 hours without the aid of the operator
- GSA was designed to simultaneously perform analyses on up to 6 samples
- It acquires data at predetermined time intervals. Based on the sedimentation rate of the soil particles, GSA allows for extremely accurate results and a very high repeatability of measurements, avoiding manual measurement errors.

SAMPLE PREPARATION

- 25 g of ground
- 62 ml of dispersing agent 40% sodium hexametaphosphate
- 438 ml of distilled water (adjust to the reference mark on the cylinder)

The soil is left in contact with the dispersing solution for a time ranging from 5 to 7 hours, to allow optimal dispersion of colloids.

PARTICLE SIZE ANALYZER – FIELDS OF APPLICATION

The classification criteria for soils and the acceptance of construction materials for roads, airports, dams, etc., are primarily based on particle size analysis.

The information obtained from this type of test, together with permeability tests, is used to **study filtration problems**. The sizing of filters is determined by the appropriate selection of particle size curves of the materials that constitute them.

GSA may be used by:

- Geotechnical Laboratories
- Testing and Certification Laboratories
- Regional Agencies for Environmental Protection
- Universities
- Technology of materials, structures, civil construction, industrial research, and Cultural Heritage
- Infrastructure and environment



STANDARD EQUIPMENT (6 POSITIONS)

- 7 floats
- 1 screen 15"
- 1 mouse
- 1 keyboard
- 6 anchors
- 1 anchor stick
- 7 x 500 ml cylinders
- 500 ml sodium hexametaphosphate
- 3 l of distilled water
- calibration mass in E2 of 100 g

TECHNICAL DATA AND PARAMETERS

- ground density
- force of gravity
- data acquisition time
- sample quantity
- float center of gravity
- power supply
- supply voltage
- output
- dimensions mm (LxWxH)
- weight

2,65 g/cm3

9,80 m/s2

1, 2, 4, 8, 16, 30, 60, 120, 240, 480, 1440 minutes

25 g

80 mm

100/240 VAC via external supply, 50 Hz

230 VAC - 1A

USB

1200 x 500 x 800 mm

40 – 65 - 96 kg



Compliant with Official Methods

ASTM D422 UNI CEN ISO/TS 17892-4







