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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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1. TECHNICAL FEATURES

TECHNICAL FEATURES

	MIA 2020	MIA 2020 BASIC
Power tension	100/240 Vac 50/60Hz	
Power consumption	10 VA	
Operating temperatura	10 ÷ 30 °C (suggested 15 ÷ 25 °C)	
Dimensions	210 x 370 x 380 mm (L x P x A)	
Net weight	8,5 kg	
Data output	RS232 I/O	
Standard equipment	Power supply 1 hydrorepellent floater interchangeable in mass and volume Double wall cylinder – Certified hydroalcoholic solution Accredia RMP	

HYDROSTATIC BALANCE

Relative density measurement range	0,5 ÷ 2,25	
Division	0,00001	0,0001
Repeatability	± 0,00002	± 0,0002
Internal mass autocalibration	YES	
Response time	6 s	

THERMOMETER

Type	Platinum probe Pt100 1/3 DIN	
Temperature range:	from 5° to 35°C	
Division	0,05 °C	
Repeatability	± 0,05 °C (in the range 10 ÷ 30 °C)	

ALCOHOLIC STRENGHR BY VOLUME (ASV)

ASV Range	0,050 ÷ 99,980 %vol	0,1 ÷ 99,9 %vol
Sensibility	0,005 %vol	0,1 %vol
Accuracy and repeatability	± 0,01 %vol (from 15 to 25 °C)	± 0,1 %vol (from 15 to 25 °C)

TOTAL DRY EXTRACT – only with dedicated software (not included)

Range	0,0 ÷ 505,8 g/l of sucrose	
Accuracy and repeatability	± 0,3 g/l	

MUSTS EVALUATION- only with dedicated software (not included)

Relative density $d_{20/20}$ range	1,04 ÷ 1,16	
------------------------------------	-------------	--

ACCESSORIES ON REQUEST:

- Floaters with certificate
- ACCREDIA RMP certified solutions 5 ÷ 60 %ASV

FLOATERS INTERCHANGEABLE IN MASS AND VOLUME

A particular care is devoted to the construction and calibration of the floaters.

They are all interchangeable in mass and volume to simplify their use; for example, it is possible to use two different floaters to determine the relative density of a must or a wine; in this way possible sample pollution is avoided.

In order to achieve what mentioned before, the suspension wire¹ and the Superficial treatment of the floater's glass were object of thorough studies to nullify the problems caused by the surface tension, by the wire's own volume, and by the electrostatic charges.

The floaters volume can, on request, be certified by competent organisations, such as the L.N.E.
¹Paris, the I.N.R.I.M.Turin, the PTB Germany, etc.

¹ Realised with a non magnetic and stainless material, with a diameter less or equal to 0,20 mm

2. INSTALLATION

The instrument was designed and built to resist and operate also in harsh conditions. However it is advisable to install it:

- on a rigid base without vibrations
- in an environment with stable temperature and without excessive ventilation
- using a dedicated electrical line

2.1 INSTALLATION PROCEDURE

- Level the balance acting on the front supports so that the bubble is in the centre of the circle on the level;
- unlock the balance by rotating clockwise the knob on the right side of the instrument until it stops (Fig. 1 e Fig. 3 n.11);


Fig. 1

- hang in the floater on the dedicated support (figure 5 n.2).
- connect the circular plug of the adapter to the connector on the instrument indicated in figure 2 with n. 3. Insert the plug in electrical network socket;

- 1 – Level bubble
2 - RS232 connector
3 – Power adapter socket


Fig. 2: Rear view

3. DESCRIPTION

The instrument is constituted by an hydrostatic balance and an elaborating module integrated in the balance itself.

Figure 3

- | | |
|-------------------------------|---|
| 6- Floater support | 11- Knob for blocking balance |
| 7- Cylinder and floater | 12- Knob of the cylinder centering device |
| 8- Control panel/Display | 13- Thermo probe |
| 9- Feet for level regulation | |
| 10- Cylinder centering device | |

3.1 CONTROL PANEL

The control panel has a display (Fig. 3 n. 8 and Fig 4) with 4 function keys which are associated with a different meaning according to the menu shown on the display and a Tare key (T).

Figure 4 - MIA 2020 Basic display

NOTE

**After some minutes when the instrument is not used, the display shuts off automatically.
By pressing any keys, the display turns on again.**

By pressing the "Switch off" key, it is possible to manually switch off the display at any time, excluding the backlight, without cutting off the power to the instrument.
By pressing any keys, the display turns on again.

4. POWER ON OF THE UNIT

Upon power on, the test routine is displayed and various information appears including the firmware release number.

After a few seconds, the last operating screen used appears

WARNING !

**If at power up the display shows "The instrument requires maintenance",
press a key to continue (see paragraph 8)**

5. MEASUREMENT PROCEDURES

To view the screen for the determination of density and alcoholic strength by volume ASV it is necessary to repeatedly press the "Unit" key

The instrument can operate at a temperature even different from 20°C (provided that it is in the range between 15 and 25°C) and allows you to:

to perform the relative density $d_{20/20}$ measurements of a distillate;
to visualise in few seconds and completely automatically the Alcoholic Strength by Volume of a distillate, in agreement with the official methods OIV

- through the dedicated software (not included) to determine the total dry extract of a wine or of a must for specific gravity in the range 0.99 - 1.16;
- through the dedicated software (not included) to evaluate the musts, by calculating the values of the Baumé, Brix, Oechsle and Plato degrees for $d_{20/20}$ in the range 1.04 - 1.16;


By pressing the "Unit" key in succession, the following alternative units of measurement can be activated

- g weight
- mg weight
- mN / m surface tension
- g / cm³ density

By further pressing the key the display shows the main screen

Using the metrological plate as an accessory it is possible, in the mode in grams, check the balance with masses calibrated in class E2 or perform a repeatability check of the measure.


The use of the metrological plate is also recommended during any external calibration operations (see Par. 6.2)

5.1 RELATIVE DENSITY AT 20°C OF A DISTILLATE AND ASV (ALCOHOLIC STRENGTH BY VOLUME)

It is done as follows:

- fill up the double coated cylinder with distillate or the water-alcohol solution, to carry out the measurements of specific gravity at 20 °C and the ASV, until reaching the mark (70 ml);
- place the cylinder on the instrument;
- immerse the thermo probe into the liquid;
- immerse carefully the floater clean and dry;
- agitate (for a few second) floater into the liquid slowly, with circular movement, in order to uniform the temperature (being careful not to create air bubbles adhering to the floater or the wire of floater which might falsify the measurement);
- carry out the reading:
 - a. push the cylinder against the appropriate centre-cylinders ;
 - b. wait that display of the unit shows the “**zero condition**” or the “**zero**” of the hydrostatic balance; if this should not occur (+000 or -000) carry out a manual operation of zeroing:
 - press the key [T] ;
 - on the display appears “zero” therefore it is returned to the precedent view
 - c. hang the floater in suspension and verify that it not touch the internal walls of the cylinder or the thermo probe. Correct the position of the centre-cylinders, if it is necessary, loosening the knobs that blocks it (figure 3) and put in place it again in an appropriate mode;

- 1 – Suspension
2 – Hook to replace
the floater not in use

Fig. 5

- d. read the ASV value on display
- The ASV and the specific gravity at 20 °C appear subsequently on the display to the read specific gravity; the figures(numbers) change until when the data not reach the stability.


Fig. 6

- a short acoustic signal and the ASV value shown with white fonts on black background highlight the achievement of this condition. The value remains on the display until the removal of the floater from the suspension;
- the specific gravity shown is then related to 20°C in accordance with the tables required by the official method OIV. The correction is reliable only for hydroalcoholic mixture and in the liquid temperature range 15 - 25°C. Outside this range, the results can undergo variations;

N.B. At the end of the determination, **always** place the floater in use on the appropriate hook

5.1.1 Relative density of a liquid

When the specific gravity of the liquid in analysis does not fall within the range of values that correspond to an Alcoholic Strength by Volume, the display shows only the relative density of the liquid and the temperature

5.1.2 Correction of the floater or compensation

Through the procedure of “Floater compensation” both the mass and the volume of the floater in use **are not modified**.

The so-called “compensation factor” of the floater is modified instead, so that all the floaters used will give the same response in the same standard conditions. These conditions are the determination of the relative density of distilled water and the determination of the ASV or alcoholic strength of a hydroalcoholic mixture, of known ASV and ACCREDIA RMP certified, provided on request by our company.

To access the function keep the MENU button for more than 3 s and select Correction Factor with arrows. Confirm with OK key.

The display shows the current value of the compensation factor that can be increased or decreased by five units at a time by pressing [+] and [-]; the compensation factor remains stored even if the instrument is turned off.

It is advisable to check regularly this factor.

5.1.2.1 Compensation of the floater with the Gibertini certified hydroalcoholic mixture ACCREDIA RMP

It is done as follows:

- fill the cylinder with 70 ml of water-alcohol solution with a known titre in a temperature comprises between 15 and 25°C, immerse the thermo probe, the floater and agitate;
- keep MENU button for more than 3 s, select Correction Factor and confirm with OK button
- wait the “zero” of the balance and hang the floater on the suspension;
- the values of the factor of current correction, specific gravity of the water-alcohol solution, the specific gravity of the same one at 20°C and the alcoholic strength appear on the display;
- if the value of alcoholic strength read on display correspond(coincede) to that declared on the label of water-alcohol solution used as reference, then the instrument is in condition of ideal reading; in the opposite case it is necessary to change(vary) the factor of correction of floater. Press the keys + and – taking in consideration that key + increase the read alcoholic strength and the key – decrease it. Taking into consideration that between the variations of compensation factor and those of alcoholic strength there does not exist a precise correspondence, and that such variations come very slowly: therefore it is necessary to proceed waiting a few seconds before read the results and stopping at the value nearest this reference.
- At the end of the determination, always replace the floater in use on the appropriate hook.


Fig. 7a: Compensation of the floater

5.1.2.2 Compensation of the floater with distilled water

It is done as follows:

- fill the cylinder with 70 ml of distilled water in a temperature comprises between 15 and 25°C, immerse the thermo probe, the floater and agitate.
- keep MENU button for more than 3 s, select Correction Factor and confirm with OK button
- Wait the “zero” of the balance and hang the floater on the suspension;
- the values of the factor of current correction, the specific gravity that the water should have at the temperature measured by thermo probe and the specific gravity of the water measured effectively by the instrument appear on the display;
- if these two values are equal then the instrument is in condition of ideal reading; in the opposite case it is necessary to change the correction factor of floater by means of ± keys


Fig. 7b: Compensation of the floater

5.1.3 Printing the value of ASV and of relative density

To have a printout of the value of the ASV or of the relative density of the liquid under consideration, just press [PRINT] key for more than 3 s. The printout is so composed, in the two possible cases (printout may vary according to printer model in use):

Date: 06-03-2017	Hour : 16:07	Pr. 012
Alcoholic strength: [%ASV]	:	13.24
Temp. of distillate	:	23.7 °C
Spec. gravity at 20 °C	:	0.98270
Spec. gravity of distillate	:	0.98170

Date: 06-03-2017	Hour : 15:20	Pr. 001
Temp. of liquid	:	24.5 °C
Spec. gravity of liquid	:	1.10510

On request is available a dedicated software to set a printout headset with end user name or other data.

5.1.4 Zeroing of the balance

Although the instrument automatically maintains the zero, the operator may carry out a manual zeroing operation, if it is necessary:

- make sure that the suspension is free and the floater is hung on the hook
- press the key Tare [T]; on the display appears "zero";
- if the suspension is engaged, on the display appears "Measure in progress" and zeroing it is not performed

5.2 TOTAL DRY EXTRACT

The total dry extract includes all matter that is non-volatile under specified physical conditions. The official European method is densimetric

The total dry extract is calculated indirectly basing upon the value of the specific gravity of the alcohol-free must/wine, i.e. of the must/wine from which the alcohol was eliminated, and which was brought back to the initial volume by adding distilled water.

The total dry extract is expressed by the concentration in g/l of a sucrose solution with the same specific gravity of the alcohol-free must/wine.

The utilization of the instrument allows one to determine the value of the total dry extract of a wine or a must for specific gravity including values from 0.990 and 1.160, in a temperature range of 15 °C and 25 °C. Under 15 °C and over 25 °C the results are not completely reliable.

Total Dry Extract determination it is possible only with dedicated software (not included) and for covering determination procedure please refer to software operating manual.

5.3 EVALUATION OF MUSTS

The degrees Baumé, Babo, Brix, Oechsle and Plato are in relation with the relative density $d_{20/20}$ of a must.

The use of the instrument allows to determine automatically the density and the value of the degrees Baumé, Babo, Brix, Oechsle and Plato of a must, for density $d_{20/20}$ between 1.04 and 1,16, in the range of temperature 15°C - 25 °C. Below 15°C and above 25°C the results are not completely reliable.

Evaluation of musts it is possible only with dedicated software (not included) and for covering determination procedure please refer to software operating manual.

6. SETTINGS

To modify some settings, the options menu presents different possibilities that are described in the following.

6.1 Date/Time setting

Select Date/Time Setting and follow indications on display to modify Date/Time of the unit

6.6 INTERNAL/EXTERNAL CALIBRATION

The instruments provides automatically, after settled time interval, to start the calibration procedure.

If the operator senses it necessary, he can activate a supplementary calibration selecting covering function on menu and following indications on display.

Performing an External Calibration it is compulsory to use a 50g Class E2 calibrated mass to ensure accuracy of the instrument and we suggest the use of the accessory Metrological Pan (see Par. 5)

7. SETUP MENU

To access at SETUP menu keep SETUP PARAMETERS key for at least 3 s.

Numerous functions for optimizing and configuring the scale operating parameters and communication via the RS232 serial port are available within the menu.

The instrument is configured during the final testing phase in an optimal way for all operating applications and does not require any modification of the relative parameters by the operator. If however it is necessary to configure one or more parameters differently, we recommend contacting Gibertini technical service.

The SETUP menu includes the following functions:

SERIAL PORT:	allows the configuration of the communication parameters of the RS232 port
PRINTER TYPE:	allows the selection of the type of serial printer connected to the RS232
CALIBRATION:	enables / excludes external / internal calibration operations
AUTOZERO:	scale configuration parameter
FILTRATION:	scale configuration parameter
IDLE MODE:	scale configuration parameter
LCD CONTRAST	allows the contrast (backlight) of the display to be adjusted
TIME FORMAT	allows the configuration of the TIME / DATE format
LANGUAGE	allows the selection of the display language ITALIAN / ENGLISH
STORE	enables saving of any changes made

To modify the desired parameter, enter the respective function and follow the instructions on the display for any modification. At the end of the configuration, enable the saving of the change by pressing the STORE button.

8. USE AND MAINTENANCE SUGGESTIONS

A dedicated microprocessor provides, besides controlling the most important functions of the instrument, to diagnose internal problems of the balance, which may have effects on the precision, suggesting a programmed service after 10000 hours of functioning with qualified personnel. The displayed message is **“The instrument requires maintenance press a key to continue”** which, however **does not** compromise the use of the balance. To continue press any key.

9. ACCREDIA LAT CALIBRATION

Gibertini Elettronica is an accredited laboratory Accredia LAT CENTER n. 94 for the characterization of electronic balances with capacity till 30 kg and microvolumes and Accredia RMP (Reference Material Producer) for the production of Certified hydroalcoholic solution. For details please consult accreditation documents at www.accredia.it

The certifications issued have an official value: they are documents that, according to the law and where applicable, can be used for all the procedures of certification, homologations and accreditation of products, services and systems of business quality.

Gibertini, therefore, can emit the calibration certificate that guarantees:

- the maintenance of the referral of apparatus used by the Centre with national samples of unity of the International System of Units (S.I.);
- the metrological correctness of the measuring procedure adopted by the Centre.

UPON REQUEST: ACCREDIA RMP Certified hydroalcoholic solutions from 5%ASV to 60%ASV.

10. QUICK RESTORE OF DEFAULT PARAMETERS

If you have mistakenly changed the SETUP parameters with an invalid configuration, you can restore the default parameters (factory settings) through a simple method.

- Remove the power supply connector of the scale located on the back of the same.
- Press and hold the first multifunction button on the left.
- Reinsert the power connector
- Wait for the word "SET FACTORY PARAMETER" to appear on the display.
- Now release the previously pressed key.
- In this phase the scale automatically restores the factory parameters

DISPOSAL – INFORMATION FOR USERS


According to the 2012/19/UE, 2002/96/CE and 2003/108/CE Directives, concerning the reduction in the use of hazardous substances in electrical and electronic apparatus, as well as the disposal of waste materials.

The symbol of a crossed box applied on the apparatus or on the packaging indicates that the product must be collected separately from other waste materials at the end of its useful life.

The separate waste collection of the apparatus which has reached the end of its useful life is organised and managed by the producer. The user who desires to get rid of present apparatus must therefore contact the seller and follow the given instructions.

Suitable separate waste collection for future sending of the disused apparatus for recycling, treatment and environmentally friendly disposal, contributes towards preventing any possible negative effects on the environment and on health and encourages the reuse and recycling of the materials the apparatus is made of.

Unauthorised disposal of the product by the user will lead to payment of the administrative sanctions in force in the country where it is put on the market.

customer care

info@gibertini.com

Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434

Fax (+39) 02 3541438

www.gibertini.com

sales@gibertini.com