

MicroDens

QUICK GUIDE



GIBERTINI

INTRODUCTION

Thank you for choosing **MicroDens**, the portable digital densimeter from Gibertini Elettronica.

MicroDens, is designed to accurately measure the density of liquids in seconds.

Compact, precise and intuitive, it represents the evolution of density measurement technology and is the ideal solution for laboratories, quality control, production and field applications.

Designed to combine laboratory precision and operational convenience, **MicroDens** integrates the performance of a benchtop density meter into a handheld format.

Precise

Accuracy $\pm 0.0001 \text{ g/cm}^3$

Laptop

Rechargeable lithium battery
with over 10 hours of battery life

Versatile

Ideal for all Food&Beverage
liquids analysis.

Intuitive

Color touch display and user-friendly interface

Traceable

Stores up to 2000 measurements with timestamp,
and data export via Bluetooth and USB-C.

PACKAGE CONTENTS

Upon receiving the MicroDens, check that the delivery is complete before using it for the first time.

The two PTFE tubes included in the package ensure a spare part is always available.

MicroDens

Portable Densimeter

Li-Ion battery

Rechargeable battery

PTFE tube

Filling tube ($\times 2$)

USB-C cable

For charging

Hard case

For transport

Documentation

Quick Guide



WARNING: in case of supply incomplete contact Gibertini Elettronica.

INITIAL ASSEMBLY

Just a few steps separate the instrument from the first measurement: plugging the battery, assembling the tube and a full charge.

Battery Installation

The connector is directional and prevents insertion with incorrect polarity: if you feel resistance, do not force it and check the orientation.

1. Place the instrument with the display facing downwards
2. Open the battery door by pressing the locking clip
3. Insert the battery in the housing
4. Close the door until you hear the click



Recharge

It is recommended to recharge the instrument during periods of inactivity, so that it is always ready for use.

- Charging time: about 3-4 hours
- Battery life: ~10 hours when used continuously
- Connector: USB-C (above display)

Filling Tube Assembly



1. Take one of the supplied PTFE tubes
2. Screw the tube into the appropriate fitting on the instrument body
3. Tighten firmly (without overtightening)

☐ For laboratory use, leave the filling tube attached; for outdoor transport, place the instrument and tube in the case.



TIP: : always keep a spare tube in the case.

INITIAL SETUP

Once assembly is complete, it only takes a few minutes to set up the instrument and prepare it for the first measurement.

Turning On

Press the PWR button for 3-4 seconds to turn on the instrument.

Upon startup, MicroDens automatically performs:

- System initialization
- Piston zero calibration
- Sensor reading

❑ These automatic checks ensure that the sensor starts from a correct reference condition. Wait for the sequence to complete before operating.

Initial Configuration

Language Setting

Menu > Settings > Device Settings > Language

Available languages: IT, EN, DE, FR, ES

Date and Time




Menu > Settings > Device Settings > Date & Time

Use the arrow keys to change the values, SET to confirm.

Unit of Measurement

Menu > Settings > Global Units

Density: g/cm³, kg/m³, SG
Temperature: °C, °F



❑ **Default settings** are already enabled and optimize battery life. You can customize them by going to Menu > Settings > Display.

INTERFACE AND NAVIGATION

Physical Keys

KEY	FUNCTION
PWR	On/Off (3 sec)
UP ▲	Scroll up / Increase value
DOWN ▼	Scroll down / Decrease
SET	Confirm selection / Access submenu
PUSH	Activate filling/rinsing pump

MicroDens features a touchscreen display with integrated touch controls



The physical 'PUSH' button is dedicated to the sample aspiration and expulsion functions and cannot be used for menu navigation.

Touch Display

The 2.8" touchscreen display makes it easy to :

Direct selection of menu options

Navigating between pages

Confirm/cancel operations

Saving results

Main Menu

Measurement Profile

Select measurement method and sample ID

Settings

Device Setup

Methods

Creating, Editing, or Deleting Methods

Calibration

Calibration procedure and history

History

Measurement archive (2000) and BLE data sending / Printer

Service

O-ring cleaning

TAKE A MEASUREMENT

Preparation : Startup Requirements

Check battery level

Icon at the top right with charge

Measuring cell clean

Free from residues of previous samples

Sample Quantity


keep at least 10 mL in the container ; the instrument aspirates ~1 mL

Waste container

Make sure you have a suitable container for disposing of the tested sample.

Measurement

1. Select the METHOD and set the Sample ID
2. Immerse the filling tube in the sample
3. Press PUSH to aspirate the sample (~1 mL)
4. Wait for stabilization (background turns GREEN)
5. Read the result on the display
6. Press SAVE to store data (optional)
7. Press PUSH to eject the sample
8. Repeat the procedure until it comes view the icon stable measurement : “ **Detected** ” (bottom left).



If the cell remains full for more than ~5 min, the instrument stops and asks for emptying.

Measurement

PRECISE

result in 5-20 sec.,
maximum accuracy

FAST

result in 1-5 sec,
routine checks

MANUAL

Operator decides
when to save

The GREEN background indicates a stable measure or in Limits.

The RED background indicates the measure is Out of Limits of the method settings .

CLEANING PROCEDURE

Thorough cleaning after each work session is essential to maintain the product over time, consistently accurate results and ensure proper operation of the instrument.

Standard Procedure:

1. Immerse the tube in the solution detergent indicated in the instruction manual
2. Press PUSH to fill with liquid
3. Wait 5-10 seconds
4. Press PUSH to rinse
5. Repeat 3-5 times until completely clean.

Solutions detergents and solvents recommended

Always select the solvent based on the nature of the sample.
For stubborn residue, rinse with acetone followed by isopropyl alcohol (IPA) and water.

SAMPLE TYPE	SOLVENT
Aqueous / Sugars	Deionized water
Spirits / Wines	Deionized Water or Water + Solution Ethanol
Oils / Fats	Isopropanol (IPA)
Difficult residues	Acetone followed by rinsing with IPA and water



CAUTION: Do not use concentrated acetone or aromatic solvents (benzene, toluene).

O-Ring Maintenance

O-ring maintenance is simple and infrequent: the tool guides you through each step and requires only a minimal amount of silicone grease.

Frequency and maintenance recommended : every 9-12 months lubricate the piston O -ring by going to Menu > Service > Clean Pump



TROUBLESHOOTING

Most problems can be resolved with simple checks.
If the problem persists, contact Gibertini technical support.

The instrument does not turn on

- Check battery integrity and connection
- Press and hold PWR for at least 4 seconds

Unstable measurement

- Check that there are no air bubbles in the sample
- Use PRECISE mode
- Wait for the temperature to stabilise

Inaccurate results

- Perform calibration using ISO 3696 water
- Clean the measuring cell thoroughly
- Check sample temperature (0-60°C)

Technical Support

GIBERTINI ELETTRONICA Srl

Via Vincenzo Bellini 37
20026 Novate Milanese (MI) - Italy

Tel: +39 02 3541434
info@gibertini.com
www.gibertini.com

© 2026 Gibertini Elettronica Srl - All rights reserved
MicroDens Quick Guide - Rev. 1.0

Thanks

Gibertini Elettronica would like to thank everyone who has contributed to the development of MicroDens with their valuable feedback and suggestions.

Our technical support team is always available to help you use the tool and ensure you have the best possible experience.

Enjoy using MicroDens!

The Gibertini Elettronica Team

