INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА

**INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА**

INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА
INSTALLAZIONE ED USO
INSTALLATION AND USE
УСТАНОВКА И ПОЛЬЗА

ScaleComm

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**Communication Software for
Gibertini balances**

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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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*All rights reserved.
Forbidden **any form** of print, duplication, reproduction and publication of this manual, or part of it,
without the written agreement of GIBERTINI ELETTRONICA Srl*

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Two modules constitute the application: a DLL library for the data exchange with the balance and an example module for Excel®.


To install insert the CDROM in the driver and wait for the automatic startup of the installation program or execute **SETUP.EXE** and then follow the instructions.

After installation the program modules will be available in the directory “**ScaleCom**” (usually C:\Programmi\ScaleCom or C:\Program Files\ScaleCom).

Two Excel® modules are available: **Demo1.xlt** and **Demo2.xlt**, the first one with the minimum software requirement for interfacing with the DLL, the second with the addition of a more complex acquisition task and a small graphic.

The two modules may be copied in the Microsoft Office model directory for further development of applications: so the models will be displayed selecting the **New** command in the **File** menu.

Open Excel® and then **Demo1** or **Demo2**, enter the menu **Tools**, select the option **References** and enable the communication module as described below (this is required only the first time the application is activated).

Example of Italian window

Exit the **Tool** function and run the application....

ScaleComm installation will generate a new class **Scale_RS232** with the following properties:

Communication Port:

ComPort Default value 1

Baudrate:

Baud Default value "9600"

Parity:

Parity Default value "E"

Data Bits:

DataBit Default value "7"

Stop Bits:

StopBit Default value "1"

Length of answer string (if NumCar=0 the string is accepted until an end-of-string CR character is received):

NumCar Default value 0

End-of-string character (usually CR):

StopChar Default value Chr(13)

Command for data acquisition:

ReadCommand Default value "B" + Chr(13)

Timeout (if the balance does not respond before timeout seconds the error string "NO-RS232" is returned):

MaxWaitTime Default value 0.5

If **MaxWaitTime** < 0 ScaleComm send command but it does not attend the answer.

Answer from instrument (read-only property):

ReadData Returns the received string or "NO RS232" on error.

**Example of a simple application program written in VBA (Visual Basic) using the communication dll
(other examples may be found in the Excel forms)**

```
...  
REM    MyScale definition as a Scale_RS232 class  
Private MyScale As New Scale_RS232  
REM    Temporary variable definition  
Dim TempString As String
```

```
...  
REM    Set communication port number to COM1  
MyScale.ComPort = 1  
REM    Set BaudRate=9600  
MyScale.Baud = "9600"  
REM    Read Data in TempString  
TempString = MyScale.ReadData  
END
```

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