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1 GENERAL SPECIFICATIONS

SO₂ GRAVIMETER is an automatic instrument for determination of the total sulphur dioxide (SO₂) (on the distillate obtained by steam distillation), the free sulphur dioxide for white, rosé and red wines (after treatment with BUBBLE for red wines)

Repeatability

Total sulphur dioxide: 6 mg/l
Free sulphur dioxide: 2 mg/l

Execution time (distillation excluded)

(approximate values depending on the sample type)

Total sulphur dioxide: 4 - 5 minutes
Free sulphur dioxide: 2 - 3 minutes

Autonomous working time

(approximate values obtained on the base of the maximum limits stated by the CEE legislation for the wines)

Total sulphur dioxide (max. limit for white and rosé wines 210 mg/l) about 30 analyses

Equipment supplied

- Containers for the reagents with driers to protect the solutions from the air
- 3 300 ml conical flasks
- 2 magnetic bar
- Parallel output for the printer
- Instruction manual
- 10 tips
- 100g Mass
- Reagents

General specifications

- Mains voltage: 230 V -15% +10%
- Environmental working conditions: +15 to +35 °C
- Power consumption: 120VA
- Fuse: T500 mA
- Dimensions: 42x40x45 cm (w x l x h) 12Kg

2 OPERATIVE INSTRUCTIONS

2.1 Preparation of reagents

Sulphuric acid H ₂ SO ₄ 35% ca	dilute 1:3 starting from concentrated acid
Indicator starch paste in water 1%	already prepared
Iodine (I ₂) N/50	from a solution N/10 dilute 1:5

2.2 Installation and preliminary operations

Connect the instrument (without conical flask) to the mains and switch on with the switch situated on the rear panel. The instrument is now ready for use.

Perform the filling operations pressing the "FILL" key

Make sure that the tubes are filled completely and that the reagents flow out freely from the tips. If necessary, replace the closed tips

2.3 Precautions

Perform a calibration at power on

Press and hold the MENU key for about 1 second.

The display shows the message SET UP PARAMETERS

Press the MENU key again until EXTERNAL CALIBRATION is displayed

Confirm by pressing OK

Pressing the MENU key displays the different options.

SET UP PARAMETERS, DATE AND TIME , EXTERNAL CALIBRATION

Proceed as follows

Select EXTERNAL CALIBRATION

When this function is activated, the calibration process comes immediately executed

All menu functions are automatically exited and the message appears on the display "LOAD WEIGHT"

In this phase the scale "waits" on the weighing plate for a calibration mass

Whose weight must be 100g

Once the "LOAD" phase is completed and the calibration value is acquired, the phase it is necessary to remove the calibration weight from the plate and the display shows "DOWNLOADED WEIGHT".

If all the phases are respected and the weights are correct the calibration is accepted and the display shows "CALIBRATED INSTRUMENT" for 3 seconds.

If, on the other hand, errors have occurred, the writing will appear

"CALIBRATION ERROR" and the previous calibration value will be maintained

Carefully clean the external surface of the flask, wiping away any drops before placing it on the instrument.

ATTENTION : DO NOT TOUCH THE CONICAL FLASK WHILE IS DOSING IODINE

3 DETERMINATION OF THE TOTAL SULPHUR DIOXIDE **For determination of total SO₂ on sample different from wine contact GIBERTINI SERVICE.**

Using the electronic distillator SUPER DEE, or DEE with the steam generator VADE, distil in steam current 50 ml of wine and 2 ml of H₂SO₄ diluted 1:10. Collect 250 ml of distillate in a conical flask containing 10 ml of NaOH 1 N. Add 3ml of H₂SO₄ and 5ml of starch paste, introduce a magnetic bar into the flask and place it in its seat. Press "START" and then "OK"
The analysis requires the addition of iodine solution until a blue colour is obtained.

The result (in mg/l) is shown on the display .

WARNING !

After one or more distillation, with SUPER DEE or DEE distillators, to determine the total sulphur dioxide, it is necessary to clean the column by performing at least 4 – 5 distillation of distilled water. The column is cleaned when the obtained distillate (250 ml) can be titrated with 0,1 ml of NaOH 0,1 N in presence of phenolphthaleine.

4 DETERMINATION OF THE FREE SULPHUR DIOXIDE

4.1 For white wine

Put 50 ml of wine in a conical flask, add 3ml of H_2SO_4 and 5ml of starch paste , introduce a magnetic bar make up to 250 ml with distilled water and place it in its seat.

Press "START" and then "OK"

The analysis requires the addition of iodine solution until a blue colour is obtained.

The result (in mg/l) is shown on the display .

4.2 Using the "BUBBLE"

The use of the BUBBLE allows to extract the free sulphur dioxide from a wine or a must and to fix it in a NaOH solution. The colourless solution can be titrated with SO_2 GRAVIMETER. The preliminary treatment with BUBBLE is necessary to determine the free sulphur dioxide in red and very red wines, or in musts and wines added with ascorbic acid.

- Put 200 ml of distilled water and 20 ml of sodium hydroxide 0,1M in a conical flask.
- Suspend the conical flask on the seat of the BUBBLE and then introduce the cylinder of the round flask.
- Measure 50 ml of wine into a volumetric flask at 20 °C and bring to about 10 °C.
Transfer the cooled sample to the flask of the BUBBLE and add 10 ml of sulphuric acid diluted 1:3. Close carefully and press BUBBLE key .
The bubbling will stop automatically after 10 minutes.
- Remove the conical flask add 3ml of H_2SO_4 and 5ml of starch paste , introduce a magnetic and place it in its seat.
Press "START" and then "OK"
The analysis requires the addition of iodine solution until a blue colour is obtained.

The result (in mg/l) is shown on the display .

5 ERRORS AND TECHNICAL FAULTS

Error	Cause	Solution
ANOMALOUS RESULT		
	Reagents altered	Replace reagents
	Spillage of reagents from tips	Verify the tip and/or joint and eventually replace
LACK OF TONING	Lacked distribution of Iodine	Verify aspiration and distribution. Eventually replace the tube with the joints.

customer care

info@gibertini.com

Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434

Fax (+39) 02 3541438

www.gibertini.com

sales@gibertini.com