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*All the information given in this manual is what available at the time of printing.  
The manufacturer reserves the right to make changes to the product at any time without notice.  
It is advisable to verify if there are any updates.*

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## DISPOSAL – INFORMATION FOR USERS

**According to the 2011/65/CE e 2003/108/CE Directives, concerning the reduction in the use of hazardous substances in electrical and electronic apparatus, as well as the disposal of waste materials.**

The symbol of a crossed box applied on the apparatus or on the packaging indicates that the product must be collected separately from other waste materials at the end of its useful life.

The separate waste collection of the apparatus which has reached the end of its useful life is organised and managed by the producer. The user who desires to get rid of present apparatus must therefore contact the seller and follow the given instructions.

Suitable separate waste collection for future sending of the disused apparatus for recycling, treatment and environmentally friendly disposal, contributes towards preventing any possible negative effects on the environment and on health and encourages the reuse and recycling of the materials the apparatus is made of.

Unauthorised disposal of the product by the user will lead to payment of the administrative sanctions in force in the country where it is put on the market.

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CONNECT TO SUPERDEE UPPER CONDENSER INLET

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CONNECT TO SUPERDEE LOWER CONDENSER INLET

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CONNECT TO THERMOCHILLER **RETURN** NIPPLE

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*All rights reserved.*

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or  
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

# 1 UNPACKING

If the unit has been shipped by carrier, carefully unpack the cartons and stand the distillation unit on a bench. Identify the box containing the reagents, bottles and various connectors.

A visual check should be made of all items for evidence of any damage that, if found, must be reported to the carrier and the suppliers within 3 days.

## 1.1 STANDARD ACCESSORY KIT

### Glassware

3	100 Volumetric flasks with plug
1	200 Volumetric flasks with plug
2	300 ml Conical flasks with plug
1	500 ml Dispenser
1	100 ml Dropper bottle
3	Dropper plastic pipettes

### Reagents

1	500 ml Calcium oxide CaO suspension 12%
1	100 ml Concentrated antifoam solution
1	500 ml Solution for liqueurs

### Fittings

1	Water inlet tubing and clips
1	Water drainage tubing and clip

### Spare parts

1	Cartridge for water (column of deionisation)
1	Tube for phial
1	Tube for peristaltic pump
1	Tube for steam

#### NOTE

The flasks are calibrated and interchangeable in mass, not in volume.

Remember that our flasks explicitly do not indicate ml, but only 100 and 200.

We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary, but it is imperative that the volume of wine and distillate is the same at the same temperature.

## 1.2 ELECTRICAL REQUIREMENTS

**Power surge:**

220V (-10/+15%) - 50 Hz

**Power consumption:**

3300 W max

2600 W normal

ON/OFF SWITCH


Figure 1

## 1.3 WATER REQUIREMENTS

**WATER**

Minimum pressure : 2 bar

Max temperature : 20 °C

Flow rate : 5 litres/min.

**Warning** – Verify the cleanness in relation to the water “hardness”.


Figure 2

Connect the large bore corrugated drainage tube to the stainless steel water discharge pipe. Clamp the tube with a hose-clip. (Fig. 2). Be careful to extract the tubes completely from the pipe before fixing drainage corrugated tube

Connect the mains water supply to the nipple on the Tee piece using the clear PVC tubing and clamp with a hose-clip. (Fig. 2)

## 1.4 DIMENSIONS AND WEIGHT

Dimensions: L30 x D48 x H102 cm

Net weight: 41 kg

## 2 SPECIFICATIONS

### 2.1 GENERAL DESCRIPTION

An automatic distillation unit designed for distillation of wines, spirits and beers for determination of alcohol content and volatile acidity.

The measured quantity of sample, (100 or 200 ml for alcohol content, 20 ml for volatile acidity), is placed in the distillation chamber, together with the reagents. Reagents are added to neutralise the sample and to make a conductive solution.

The distillate is collected in a volumetric flask mounted on a precision balance arm. The balance will stop the distillation when a preset volume has been collected.

For volatile acidity, the distillate is collected in a 300 ml conical flask.

The SUPER DEE can be used for determination of other analytical parameters of the wines (e.g. sorbic acid, total sulphur dioxide). For further informations see the book "Analysis methods of wines and alcoholic beverages" Gibertini – current Edition.

The SUPER DEE has an ecologic device for the waste water: the water passes through the cooling coil only during the distillation, and it cleans the drainage basin.

The SUPER DEE complies with the directives 2006/42 EEC and 85/375/EEC for the safety of the analysts. The instruments are also in accordance with the directive 2004/108 EEC (electromagnetic compatibility).

### 2.2 EXECUTION TIMES

Wines	100 ml	3 - 4 minutes
Wines	200 ml	6 - 7 minutes
Volatile acidity	250 ml	6 - 7 minutes
Spirits		4 - 8 minutes

### 2.3 COMPARISON WITH THE OFFICIAL METHODS

The instruments are in conformity with OIV methods as prescribed by Regulation n. 606/2009 EEC, both for the Alcoholic Strength by Volume (ASV) and for volatile acidity.

#### 2.3.1 Alcoholic strength by volume (ASV)

For the wines, the current regulations permit the use of any model of apparatus for distillation or extraction with a steam current on condition that it satisfies the following test:

- Distill 5 times in succession a 10 % vol water-alcohol mixture. After the last distillation, the distillate should present an Alcoholic Strength by Volume of almost 9,9 % vol. In other words, during each distillation, the instrument does not lose a quantity of alcohol higher than 0,02 % vol.

For spirits, the following test must be satisfied:

- Distil 200 ml of water-alcohol solution, with well-known tritre near 50 % vol. The instrument does not lose a quantity of alcohol higher than 0,1 % vol.

#### 2.3.2 Volatile acidity

The ECC legislation does not describe a unique type of distillation apparatus, but it states that any apparatus or technique, which satisfies the following 3 requirements, constitutes an official international apparatus or technique:

- 1) The steam of water produced by the generator must be free from CO<sub>2</sub>. Distill 20 ml water collecting 250 ml of distillate, adding 0,1 ml di NaOH 0,1 N and 2 drops phenolphthalein 1 % in alcohol: the colour must remain stable pink for at least 10 seconds.
  - 2) 99,5 % of acetic acid contained in a 0,1 M solution must be found in the distillate. Distil 20 ml acetic acid 0,1 M. Collect 250 ml distillate and titrate with at least 19,9 ml NaOH 0,1 N.
  - 3) A maximum of 0,5 % of a 1 M solution of lactic acid may distil. Distil 20 ml of a M solution of lactic acid, collect 250 ml of distillate and tritrate with maximum 1,0 ml NaOH 0,1 N.
-

### 3. OPERATIVE MENU

#### 3.1 LANGUAGE SELECTION MENU

To select a language different than the one currently set, operate as follows:

- Switch off the distillator with the ON/OFF SWITCH.
- Switch on again pressing simultaneously the **SHIFT** key on the keyboard and keep pressed until the end of the initial test.
- The display shows the current language and the shifting arrows in correspondence of **A** and **C** keys. Select the language with these two keys.
- Confirm with the **MODE** key.

#### 3.2 OPERATIVE MENU

*See also Chapter 6 (ANALYSIS – Example of distillation)*

The instrument display is composed by two lines of text. The first shows the operative mode and the second describes the function corresponding to **A**, **B** and **C** keys.

<b>Mode: Alcohol</b> <b>Alcohol Volat.</b>
---

Switch On the instrument. The display shows the menu with the analysis type selection (alcohol or volatile)

By selecting the Alcohol Mode, the instrument it is ready for a distillation with electrodes. By selecting Volatile Mode, the instrument is ready for a distillation with steam and it starts to warm the boilers.

For a correct distillation with steam, it is recommended to predispose the instrument in volatile mode some minutes before the distillation and if needed to perform a distillation of water only.

By pressing the **MODE** key, it is possible to proceed to the next menu,

<b>Set : 80</b> <b>80            160            240</b>
--

The software has 3 definite set points : **80**, **160** and **240** grams, selectable by pressing the **A**, **B** or **C** key respectively

The precision balance shows the gram.

<b>Set : 87</b> <b>-                    +</b>
--

It is possible to modify these values by pressing the **SHIFT** key with the **A** (-) or **B** (+) key

By pressing the **MODE** key, it is possible to proceed to the distillation start menu,

<b>Standby</b> <b>Start</b>
--------------------------------

It is possible to start the distillation by pressing the **A** key, correspondent to the **START** function. The distillation does not start without a volumetric flask or a conical flask on the balance.

Once the distillation has started, the display shows the weight value collected and the total volume to be collected.

<b>15            -&gt; 80</b> <b>Stop</b>
--

When the distillate reaches the selected quantity (80), the instrument does not accept a new start if the flask is not removed.

It is possible, at any moment, to stop the distillation by pressing the **B** key (**Stop**).

If the distillation is running, only the **STOP** key (**B**) is active. If the instrument is not active, it is possible to go to the following menu by the **MODE** key, to select the quantity of distillate to collect.

## 4 DETERMINATION OF VOLATILE ACIDITY

Pour into the distillation chamber 20 ml wine, free of CO<sub>2</sub>. It is advisable to add 1 ml tartaric acid 50 % (for example if the wine has a fixed acidity less than 6 g/l).

Lower the plexiglas shield, close the distillation chamber and insert a conical flask on the balance, otherwise the distillation DOES NOT START (see figure 3), and press START push button.

The distillation stops automatically as soon as about 250 ml of distilled liquid have been collected. The time of distillation is 6-7 minutes.

At the end of the distillation, add to the distillate a few drops of phenolphthaleine 1 % and titrate with NaOH 0,1 N until the distillate appears of a stable pink colour for at least 10 seconds.

Be **n** the millilitres used.

$$\text{Gross volatile acidity} = n \times 0,3$$

*In order to obtain the correct value of volatile acidity, it is necessary to determine the amount of sulphur dioxide present in the distillate, whose quantity must be subtracted from the gross volatile acidity.*

**Figure 3**

### 4.1 SUBTRACTION OF SULPHUR DIOXIDE SO<sub>2</sub>

Add to the already titrated distillate a few drops of sulphuric acid solution diluted 1:4 and about 5 ml of starch paste 1%. The solution will become colourless.

Titrate the free sulphur dioxide with I<sub>2</sub> N/100 until the solution turns blue. Be **n'** the millilitres used.

Now make the solution alkaline by adding 20 ml of saturated sodium borate solution.

The solution will turn pink .

Titrate with I<sub>2</sub> N/100 until the appearance of the blue colour. **n''** is the millilitres used.

Subtract from ml of NaOH 0,1 N (**n**) a quantity equivalent to:  $\frac{n'}{10} + \frac{n''}{20} = n'''$

Correcy volatile acidity = (**n** - **n'''**) x 0,3 (g/l acetic acid).

The automatic enological analyser GIBERTINI QUICK allows the titration of the distillate for the determination of volatile acidity, either gross or correct.

## 5 WINES AND SPIRITS DISTILLATION

### 5.1 THE DISTILLATION PROCEDURE FOR WINES

May the young fermenting wines or sparkling wines have any CO<sub>2</sub> removed by agitation or filtering through filter paper, or by short immersion in an ultrasonic bath.

Distillation procedure	Observations
1. The wine is poured up to the mark of the special flask <sup>5</sup>	The temperature of the sample should be 20 °C. One can also operate at room temperature: from the sampling and the bringing back to volume after distillation a maximum deviation of $\pm 2$ °C is allowed.
2. The content is poured into the distillation chamber of the SUPER DEE	
3. The flask is washed 3 time with a few ml of distilled water that is also poured into the chamber	
4. The sample is fully alkalized by the addition <sup>6</sup> of a suspension of milk of lime (2M), containing 120 g/litre CaO	The use of phenolphthalein indicator paper can be useful. However, the coulour change in the wine is obtained, visually, because of the alkalinity of the wine. This procedure is necessary to avoid the distillation of the other volatile substances (acetic acid, sulphurous anhydride, aldehydes, etc.), which will interfere with the density. In fact, the desired goal is to obtain a distillate similar to a hydroalcolohic solution, that is to say a pure alcohol-water mixture.
5. A few drops of antifoam is added	This is an aqueous suspension of oxyethylenic polymers. Foaming in the bulb slows the distillation and disrupts the passage of ethanol into the distillate. The presence of sugars in the sample enhances the foam formation.
6. Add approximately 2 ml distilled water to the same flask the sample was measured in. Collect approximately 75 ml of distillate therein.	The addition of water serves to avoid the loss of alcohol. The distillation is stopped automatically by a balance mechanism when the correct quantity has been distilled.
7. Add distilled water to the flask nearly up to volume. It is stoppered and brought to volume once temperature of 20 °C or the temperature of the sample before distillation is reached.	The use of a thermostatic bath is suggested
8. The flask is stoppered and mixed with care to avoid bubbles	The aim is to ensure a homogeneous distilate because, the determination of gravity will proceed from this.

<sup>5</sup> The flasks are calibrated and interchangeable in mass, not in volume.

Remember that our flask explicitly DO NOT indicate ml, but only 100-200.

We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary, but it is imperative that the volume of wine and distillate is the same at the same temperature.

<sup>6</sup> The exact quantity depends on the sample.

## 5.2 THE DISTILLATION PROCEDURE FOR SPIRITS

Spirits present different characteristics from those of wines:

- Values of ASV up to 70%
- Variable composition

In fact spirits with very different composition exist: from cream whisky to fruit creams, from very sweet liqueurs such as vodka with fruit or limoncello to liqueurs with extracts of aromatic plants or roots, or even drinks with egg or coffee.

In addition, the higher alcoholic strength results in:

- a very slow distillation in order to make sure of the complete extraction of alcohol;
- collection of a greater volume of distillate (around 95 ml).

It is difficult to find a common distillation method for all the different kinds of spirits because the various ingredients interfere with the distillation in different ways, sometimes blocking the passage of the alcohol.

The basic procedure for spirits is listed below. This procedure allows carrying out a good distillation in a short time, without excessive foam, and resulting in a clear distillate; the specific methods applicable to different alcohol products are indicated in the subsequent paragraph.

The distillation procedure	Observations
1. The sample is poured up to the mark of the special flask <sup>7</sup>	It is recommended that the temperature of the sample is 20 °C. One can also operate at room temperature: a maximum deviation of $\pm 2$ °C is allowed between the sampling and the bringing back to volume after distillation.
2. Pour the contents into the bulb of the distillatory	
3. Wash the flask 3 times with 10 ml of distilled water and add washings to distillatory bulb	
4. Add the “ <b>solution for liqueurs</b> ” <sup>8</sup> : solution 1:1 of NaCl 1% and [KAl(SO <sub>4</sub> ) <sub>2</sub> ] 12%	The aim is to increase the conductivity of the sample and to suppress the aromatic substances.
5. Add a few ml of antifoam	Its addition avoids the formation of foam which would be an obstacle to the normal distillation of the alcohol.
6. Add about 2 ml cold distilled water to the same flask used for measurement of the sample. Collect $\approx$ 95 ml of distillate therein.	
7. Add distilled water nearly up to volume, stopper and bring to volume when the temperature reaches 20 °C or the temperature of sample before distillation	USE OF A THERMOSTATIC BATH IS ADVISED
8. Stopper the flask and mix with care to avoid bubbles	

Regarding brandy, etc. the method is the same as the one indicated for wines, with the exception of the addition of a small quantity of lime milk: because of the low acidity of such products, 1-2 ml is sufficient. In case of alcoholic strength higher than 50 % vol, it is advisable to dilute to 1:2

<sup>7</sup> See note 1

<sup>8</sup> For the quantity see Table 1

### 5.3 COMPARISON BETWEEN THE SUPER DEE AND THE OIV GLASS DISTILLING UNIT

Tests of distillation with spirits with extremely varied compositions were carried out. Each liqueur was distilled :

- with the SUPER DEE distilling unit following the above procedure;
- with a glass distilling unit using the official procedures of OIV.

This latter distilling unit satisfies the official OIV condition.

The results obtained with the two apparatus are shown in the **TABLE 1**.

A specific method was developed to be used with the SUPER DEE distilling unit for each different sample type. Generally it is recommended:

- very sweet liqueurs: use small quantities of distilled water to wash the flask (less than 10 ml); in many cases this solution blocks the formation of emulsions which can result in an incorrect distillation;
- liqueurs containing milk or cream: the solution for liqueurs should not be added (it may curdle the milk creating a clot that is difficult to clean out of the flask).

The results obtained with the use of the SUPER DEE distilling unit are without doubt comparable with those obtained with the OIV glass distilling unit, with the advantage of much shorter distillation times. The differences found from the experimental data on some samples compared to the values indicated on the label fall within the tolerance expected by the current legislation (0,3 % vol on indicated ASV).

For further information see the book “Analysis methods of wines and alcoholic beverages”, Gibertini – current Edition.

**TABLE 1**

<b>SAMPLE</b>	<b>INGREDIENTS</b>	<b>DEE PROCEDURE</b>	<b>%vol DEE</b>	<b>Declared</b>
<b>Liqueur type Strega</b>	Alcohol water, sugar, natural flavours, saffron	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	40.37	40
<b>Lemon cream</b>	Water, sugar, lemon infusion, milk, cream, alcohol, glucose syrup	100 ml +12 ml antifoam	22.05	22
<b>Rakee of raspberries</b>	rakee, raspberries juice, sugar, natural flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	32.02	32
<b>Limoncello</b>	Alcohol, sugar, infusion of lemon rind, natural flavours	100 ml + 3 ml antifoam +10 ml milk of lime	30.29	30
<b>Whisky cream</b>	Dairy cream, sugar, Whisky, water, natural flavours	100 ml +15 ml antifoam	14.71	15
<b>Rakee</b>	Rakee	100 ml +5 drop antifoam +3 ml solution for liqueurs	40.12	40
<b>Pina Colada</b>	sugar, alcohol, cream, rum, milk, flavours	100 ml + 9 ml antifoam	16.80	17
<b>China Martini</b>	Sugar, alcohol, flavours, colorants: caramel; acidifier: citric acid	100 ml +15 ml antifoam + 10 ml milk of lime	31.00	31

*(continue)*

<b>SAMPLE</b>	<b>INGREDIENTS</b>	<b>DEE PROCEDURE</b>	<b>%vol DEE</b>	<b>Declared</b>
<b>White Vermuth</b>	Wine, sugar, alcohol, flavours	100 ml +1.5 ml antifoam + 5 ml milk of lime	16.06	16
<b>Strong anisette</b>	Alcohol, sugar, natural flavours	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	40.02	40
<b>Liqueur with coffee</b>	sugar, alcohol, coffee, natural flavours	100 ml + 12 ml antifoam	32.89	33
<b>Cherry</b>	Alcohol, sugar, egriot juice, natural flavours	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	30.03	30
<b>Tequila</b>	Distillate of agave techilana	100 ml +1.5 ml antifoam +2 ml solution for liqueurs	38.14	38
<b>Pastis</b>	Alcohol, sugar, 2 <sup>0</sup> / <sub>00</sub> anetolo, infusion of vegetable substances, natural flavours,	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	44.76	45
<b>Banana Vodka</b>	vodka, sugar, banana juice, flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	25.17	25
<b>Brandy</b>	Alcohol hydrate, Brandy, natural flavours, colorant: caramel	100 ml + 3 ml antifoam +2 ml solution for liqueurs	36.07	36
<b>Liqueur with chocolate and mint</b>	Milk, sugar, cream, glucose syrup, alcohol, milk protein, infusion of cacao, flavours, colorant: caramel	100 ml + 12 ml antifoam	16.06	16
<b>Cachaca</b>	Distillate of sugar-cane	100 ml + 3 ml antifoam +3 ml solution for liqueurs	43.66	43
<b>Amaretto (alcoholized for exclusive use in confectionery)</b>	Alcohol, caramel, flavours	25 ml + 75 ml distilled water (collect 100 ml) + 3 ml antifoam +3 ml solution for liqueurs	% vol value × 4 = 70.00	70
<b>Sake</b>	Distillate of rice	100 ml + 1 ml antifoam +10 ml solut. for liqueurs	15.85	16
<b>Sambuca</b>	Alcohol, sugars, anicette, natural flavours	100 ml + 5 ml antifoam +8 ml solution for liqueurs	39.86	40
<b>Cointreau</b>	Alcohol, sugars, orange flavours	100 ml + 5 ml antifoam +7 ml solution for liqueurs	39,81	40

**NOTE**

The samples are only indicative, so they must not be considered exhaustive.

## 6 ANALYSIS – Example of distillation

### 6.1 ANALYSIS TYPE SELECTION

Figure 4

80 / 160 for alcohol distillation  
 240 for volatile

Push **SHIFT** to modify the data.  
 The display shows + 0 -

#### TO SELECT THE KEY FOR USE :

80	=	100 ml	key <b>A</b>	} ALCOHOL
160	=	200 ml	key <b>B</b>	
240	=	240 ml	key <b>C</b>	} VOLATILE

### 6.2 SELECTION MODE - ALCOHOL / VOLATILE

To select ALCOHOL mode push the **A** key.

Confirm with the **MODE** key.


Figure 5

To select VOLATILE mode push the **B** key.

Confirm with the **MODE** key.


Figure 6

### 6.3 START DISTILLATION

Start distillation pushing the **A** key.

Figure 7

#### WARNING

Before starting the distillation check that the plug is indeed closed (turn clockwise the knob), the flask is installed on the balance and the protection is completely lowered. If there aren't these conditions the distillation don't start and an acoustic signal informs of the anomaly.

Example of distillation:  
-plug and drainage closed  
-flask on the balance  
-protection down.

**Figure 8**

**6.4 DISTILLATION CHAMBER CLEANING**

After the distillation open the drainage and clean the distillation chamber with the shower supplied with SUPER DEE

**Figure 9**

When the SUPER DEE distilling unit is not in use leave the drainage open.

**Figure 10**

## 7 ERROR MESSAGES

There are indications for the main anomalies. The display shows the messages with the operative situations and the lights **WARNING** and **ALARM**.

Remember that, for a correct operation, the system must have a capacity of 5 litres/min. and that the water must have a temperature of 20°C max and a pressure of at least 2 bar.

- 1) In case of NO water or low water pressure:

**Error !  
Water press. low**

*THE INSTRUMENT DOES NOT START (Fig. 11)*

- 2) After 10 seconds the display shows the message

**No water**

*THE INSTRUMENT DOES NOT START*

If the error situation continues the distillation will be cancelled and the program will wait to confirm the reading of the message by operator.

**No water  
Ok**

*CONFIRM THE READING OF THE MESSAGE.*


- 3) In case of cooling water too hot the yellow light flashes. Push the START key to start a distillation. The instrument will wait 10 seconds and if the error situation remains the operator will select if to continue or cancel the distillation.

**Water too hot  
Go on                      Cancel**


*CONFIRM OR CANCEL THE DISTILLATION  
(Fig. 12)*

### ATTENTION!

**The temperature of the cooling water upper 20 °C compromises the result of the distillation whether in alcohol or Volatile acidity mode.**

*Figure 11*

*Figure 12*

Many situations of error have the indication of the error message until the restored normal situation.

- 4) In case of drainage open the following message will appear:

**Attention!  
Drainage open**

*THE DISTILLATION DON'T START.  
THE INSTRUMENT WAITS FOR THE CLOSURE  
OF THE DRAINAGE (Fig. 13)*


- 5) In case of open plexiglas the following message will appear:

**Attention!  
Plexiglass open**

*THE DISTILLATION DOES NOT START.  
THE INSTRUMENT WAIT FOR THE CLOSURE  
OF THE PLEXIGLASS (Fig. 14)*

*Figure 13*

*Figure 14*

- 6) If the flask is not positioned correctly, errors in the quantity of distillate collected may occur. Position the flask as shown in the figure on the side, so as to avoid contact between the flask and the glass tube for the collection of the distillate.

## 8 MAINTENANCE

### 8.1 DISTILLATION CHAMBER

After each distillation the distillation chamber must be washed with water to prevent contamination of the successive sample. Periodically clean the distillation chamber introducing approximately 500 ml of distilled or demineralised water and approximately 20 ml of sodium chloride NaCl 1% solution. Start the distillation in ALCOHOL MODE. To remove deposits or incrustations that the samples can have left on the electrodes and/or on the glass can be done in the two following ways:

1- Pour into the distillation chamber a dilute solution of HCl 0,1N. Await some hour and eventually use a tube-brush to clean the chamber.

2- Introduce into the distillation chamber a solution 50% NaOH 1M + sodium hypochlorite 2,5 ÷ 4%, in sufficient quantity to cover the incrustations. Also in this case to await some hour and use eventually a tube-brush.

In any case rinse always accurately the distillation chamber before reuse the instrument.

**WARNING! Don't press the START button when into the distillation chamber there is a cleaning solution. For more safety remove the power socket from the main.**

### 8.2 ELECTRODES

The electrodes will gradually diminish in size and will have to be replaced at some stage. It is difficult to predict the frequency of replacement because it will depend on total usage and type of sample (approx. 2000 tests). **For the replacement call the SAT (Service) of GIBERTINI or your local supplier.**

### 8.3 RECTIFICATION COLUMN

During the normal use of the SUPER DEE, the column doesn't require a particular cleanness. However, if the column is contaminated by acid residues, it is necessary to clean it following indication reported on procedure 70.30.136 Rev.1 available as appendix of this user manual

### 8.4 EMPTYING FUNCTION

- Switch off the DISTILLATOR with the ON/OFF SWITCH.
- Switch on pressing simultaneously the **B** key on the keyboard.
- The display shows the message "**Take off water tube Wait for tank drain**"
- Empty the cooling coil using the shower.
- Wait for the complete emptying of the distilled water tank by means of the steam pump (the water finishes to go out by the nozzle and the pump becomes very noisy). Switch off the instrument.
- Incline the instrument backwards to help the exit of water.

## 8.5 DEIONISATION COLUMN (REFILL R3)


Figure 15

Replace the column when it is completely red (after about 400 volatiles). It is advisable however to replace the cartridge when it is black only in the last 3-5 cm.

Replacement of the cartridge:

- unscrew the nuts that keep the cartridge in position;
- incline the cartridge leaving the bottom end in position. Close the water inlet hole (on top) with a finger to prevent water from spilling over the components underneath;
- substitute the cartridge paying attention that the inlet and outlet gaskets of the cartridge are in the correct position;
- screw in the fixing nuts.

Verify that the water pipes do not have bents that may prevent water to flow freely.

## 8.6 BOILERS

If the display shows **CHECK BOILERS**, push the 2 red keys (Fig. 15).

### NOTE

Arrange a possible periodic maintenance with the GIBERTINI SAT (Service) - tel. (+39) 02 3541434

## 9 OPTIONAL CONNECTION WITH “THERMOCHILLER”

In some operating environments, it may happen that the network water, used as a heat exchanger in the refrigerant, is available at temperatures that are not low enough to allow adequate cooling, reducing the temperature difference necessary to ensure effective distillation.

In these cases the use of the Gibertini liquid circulation cooler mod. Thermochiller (or other suitable cooler) is strongly recommended.

### 9.1 DISTILLER INTENDED FOR USE WITH “THERMOCHILLER” (FACTORY SET)

If the distiller is purchased with the Thermochiller as accessory, the instrument is set at the factory with the appropriate electronic configuration and complete with connections, pipes and accessories necessary for use with the cooler.

Compared to the standard version, the electronic display, every time a distillation is started, displays the warning message for a few seconds for the preparation and ignition of the relative Thermochiller and then the distillation itself is performed.

**ATTENTION: in this configuration, if Thermochiller is not activated, no liquid flow through the refrigerant cooling circuit and the distillation will not be successful, with possible overheating of the system.**

### 9.2 DISTILLER INTENDED FOR USE WITH TAP WATER

If the distiller is subsequently equipped with the Thermochiller accessory, the instrument will not have the appropriate electronic setting for use with the cooler.

To configure the instrument for use with ThermoChiller proceed as follows:

- turn off the instrument and restart it later by holding down the "C" key
- the chiller use selection screen will be displayed
- confirm by pressing the YES key
- wait for the normal restart of the instrument

After that the distiller will operate in the "Chiller" configuration and must be used only with the relative cooler in operation.

To return to the mains water cooling configuration, carry out the above procedure, selecting "NO" in the chiller use selection screen.

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CONNECT TO THERMOCHILLER **OUTLET** NIPPLE

APPENDIX1: Procedure 70.30.136 How to clean the distillation chamber and rectification column

	<b>CLEANING MODE</b>	<b>WHEN TO DO IT</b>	<b>PROCEDURE</b>
1-	<b>Standard cleaning</b>	Periodically, after a series of VOLATILE distillations.	<ul style="list-style-type: none"> <li>a. Pour in the distillation chamber about 20ml of distilled water, and perform distillation in VOLATILE mode.</li> <li>b. Repeat once or twice.</li> </ul>
2-	<b>Cleaning of rectification column 1<sup>st</sup> mode</b>	When the column appears to be polluted by acid residuals (i.e. when the distillate of water requires more than 0.1ml of NaOH 0.1N to be titrated).	<ul style="list-style-type: none"> <li>a. Fill the distillation chamber with about 500 ml of distilled water + about 20 ml of NaCl 1% and perform a distillation in ALCOHOL mode.</li> </ul> <p><b>N.B.-</b> <i>This type of washing performs also a mechanical action on the column, since the water goes up from the distillation chamber into the column.</i></p>
3-	<b>Cleaning of rectification column 2<sup>nd</sup> mode</b>	When the column filling is made black by the products of the distillations, or when acidity persists during volatile distillations (the distillate of water requires more than 0.1ml of NaOH 0.1N to be titrated).	<ul style="list-style-type: none"> <li>a. Pour in the distillation chamber about 1ml of a bleach (sodium hypochlorite), add about 500 ml of distilled water and perform a distillation in ALCOHOL mode.</li> <li>b. Repeat cleaning as indicated in point 2 (1<sup>st</sup> mode).</li> <li>c. Pour in the distillation chamber about 20 ml of distilled water, and perform distillation in VOLATILE mode.</li> <li>d. Repeat distillation of water until the water distillate requires less than 0.1ml of NaOH 0.1N to be titrated.</li> </ul>
4-	<b>Cleaning of distillation chamber</b>	When the distillation chamber appears dirty or there are residual inside.	<ul style="list-style-type: none"> <li>a. Pour in the distillation chamber about 100 ml of bleach, and let it rest for about 30 minutes.</li> </ul> <p><b>DO NOT DISTILL !!!</b></p> <ul style="list-style-type: none"> <li>b. Drain and wash.</li> <li>c. Pour in the distillation chamber about 20 ml of distilled water, and perform distillation in VOLATILE mode.</li> <li>d. Repeat distillation of water until the water distillate requires less than 0.1ml of NaOH 0.1N to be titrated.</li> </ul>

customer care

info@gibertini.com

Tel. (+39) 023541434

ISO 9001:2015

LAT N° 094

---

**GIBERTINI**  
**ELETTRONICA s.r.l.**

Via Bellini 37  
20026 Novate Mil.se  
Milano (Italy)

Tel. (+39) 02 3541434  
Fax (+39) 02 3541438  
www.gibertini.com  
sales@gibertini.com

---