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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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CONNECT TO SUPERDEE UPPER CONDENSER INLET

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CONNECT TO SUPERDEE LOWER CONDENSER INLET

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CONNECT TO THERMOCHILLER **RETURN** NIPPLE

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1 UNPACKING

If the unit has been shipped by carrier, carefully unpack the cartons and stand the distillation unit on a bench. Identify the box containing the reagents, bottles and various connectors.

A visual check should be made of all items for evidence of any damage that, if found, must be reported to the carrier and the suppliers within 3 days.

1.1 STANDARD ACCESSORY KIT

Glassware

- 3 100 Volumetric flasks with plug
- 1 200 Volumetric flasks with plug
- 1 500 ml Dispenser
- 1 100 ml Dropper bottle
- 2 Dropper plastic pipettes

Reagents

- 100 ml Concentrated antifoam solution
(code 0700034)
- 500 ml Calcium oxide CaO suspension 12%
(code 0700033)
- 500 ml Solution for liqueurs
(code 0700096)

Fittings

- 1 Water inlet tubing and clip
- 1 Water drainage tubing and clip

NOTE

The flasks are calibrated and interchangeable in mass, not in volume.

Remember that our flasks explicitly do not indicate ml, but only 100 and 200.

We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary, but it is imperative that the volume of wine and distillate is the same at the same temperature.

1.2 ELECTRICAL REQUIREMENTS

Power source:
230V (-10/+15%) - 50 Hz

Power consumption:
3300 W max
2600 W normal


Figure 1

1.3 WATER REQUIREMENTS

WATER

Minimum pressure : 2 bar
Max temperature : 20 °C
Flow rate : 5 litres/min.

Warning – Verify the cleanness in relation to the water “hardness”.


Figure 2

Connect the large bore corrugated drainage tube to the stainless steel water discharge pipe. Clamp the tube with a hose-clip. (Fig. 2). Be careful to extract the tubes completely from the pipe before fixing drainage corrugated tube

Connect the mains water supply to the nipple on the Tee piece using the clear PVC tubing and clamp with a hose-clip. (Fig. 2)

1.4 DIMENSIONS AND WEIGHT

Dimensions: L30 x D48 x H102 cm

Net weight: 41 kg

2 SPECIFICATIONS

2.1 GENERAL DESCRIPTION

An automatic distillation unit designed for distillation of wines, spirits and beers for determination of alcohol content.

The measured quantity of sample, (100 or 200 ml for alcohol content), is placed in the distillation chamber, together with the reagents. Reagents are added to neutralise the sample and to make a conductive solution.

The distillate is collected in a volumetric flask mounted on a precision balance arm. The balance will stop the distillation when a preset volume has been collected.

The SUPER DEE can be used for determination of other analytical parameters of the wines (e.g. sorbic acid, total sulphur dioxide). For further informations see the book *"Analysis methods of wines and alcoholic beverages"* Gibertini – current Edition.

The SUPER DEE has an ecologic device for the waste water : the water passes through the cooling coil only during the distillation, and it cleans the drainage basin.

2.2 EXECUTION TIMES

Wines	100 ml	3 - 4 minutes
Wines	200 ml	6 - 7 minutes
Spirits		4 - 8 minutes

2.3 COMPARISON WITH THE OFFICIAL METHODS

The instruments are in conformity with OIV methods as prescribed by Regulation n. 606/2009 EEC, both for the Alcoholic Strength by Volume (ASV) and for volatile acidity.

2.3.1 Alcoholic strength by volume (ASV)

For the wines, the current regulations permit the use of any model of apparatus for distillation or extraction with a steam current on condition that it satisfies the following test:

- Distil 5 times in succession a 10 % vol water-alcohol mixture. After the last distillation, the distillate should present an Alcoholic Strength by Volume of almost 9,9 % vol. In other words, during each distillation, the instrument does not lose a quantity of alcohol higher than 0,02 % vol.

For spirits, the following test must be satisfied:

- Distil 200 ml of water-alcohol solution, with well-known titre near 50 % vol. The instrument does not lose a quantity of alcohol higher than 0,1 % vol.

3. OPERATIVE MENU

3.1 LANGUAGE SELECTION MENU

To select a language different than the one currently set, operate as follows:

- Switch off the distillator with the ON/OFF SWITCH.
- Switch on again pressing simultaneously the **SHIFT** key on the keyboard and keep pressed until the end of the initial test.
- The display shows the current language and the shifting arrows in correspondence of **A** and **C** keys. Select the language with these two keys.
- Confirm with the **MODE** key.

3.2 OPERATIVE MENU

See also Chapter 5 (ANALYSIS – Example of distillation)

The instrument display is composed by two lines of text. The first shows the operative mode and the second describes the function corresponding to **A**, **B** and **C** keys.

Switching on the instrument, the display shows the standard setpoints of the distillate to collect (see Figure 3).

By pressing the **MODE** key, it is possible to select the quantity:

Set : 80
80 160 240

The software has 3 definite setpoints : **80**, **160** and **240** grams, selectable by pressing the **A**, **B** or **C** key respectively

The precision balance shows the gram.

Set : 87
- +

It is possible to modify these values by pressing the **SHIFT** key with the **A** (-) or **B** (+) key

By pressing the **MODE** key, it is possible to proceed to the distillation start menu,

Standby
Start

It is possible to start the distillation by pressing the **A** key, correspondent to the **START** function. The distillation does not start without a volumetric flask or a conical flask on the balance.

Once the distillation has started, the display shows the weight value collected and the total volume to be collected.

15 -> 80
Stop

When the distillate reaches the selected quantity (80), the instrument does not accept a new start if the flask is not removed.

It is possible, at any moment, to stop the distillation by pressing the **B** key (**Stop**). If the distillation is running, only the **STOP** key (**B**) is active.

4 WINES AND SPIRITS DISTILLATION

4.1 THE DISTILLATION PROCEDURE FOR WINES

May the young fermenting wines or sparkling wines have any CO₂ removed by agitation or filtering through filter paper, or by short immersion in an ultrasonic bath.

Distillation procedure	Observations
1. The wine is poured up to the mark of the special flask ⁵	The temperature of the sample should be 20 °C. One can also operate at room temperature: from the sampling and the bringing back to volume after distillation a maximum deviation of ± 2 °C is allowed.
2. The content is poured into the distillation chamber of the SUPER DEE	
3. The flask is washed 3 time with a few ml of distilled water that is also poured into the chamber	
4. The sample is fully alkalized by the addition ⁶ of a suspension of milk of lime (2M), containing 120 g/litre CaO	The use of phenolphthalein indicator paper can be useful. However, the coulour change in the wine is obtained, visually, because of the alkalinity of the wine. This procedure is necessary to avoid the distillation of the other volatile substances (acetic acid, sulphurous anhydride, aldehydes, etc.), which will interfere with the density. In fact, the desired goal is to obtain a distillate similar to a hydroalcolohic solution, that is to say a pure alcohol-water mixture.
5. A few ml of diluted antifoam is added	This is an aqueous suspension of silicones or of oxyethylene polymers. Foaming in the bulb slows the distillation and disrupts the passage of ethanol into the distillate. The presence of sugars in the sample enhances the foam formation.
6. Add approximately 2 ml distilled water to the same flask the sample was measured in. Collect approximately 75 ml of distillate therein.	The addition of water serves to avoid the loss of alcohol. The distillation is stopped automatically by a balance mechanism when the correct quantity has been distilled.
7. Add distilled water to the flask nearly up to volume. It is stoppered and brought to volume once temperature of 20 °C or the temperature of the sample before distillation is reached.	The use of a thermostatic bath is suggested
8. The flask is stoppered and mixed with care to avoid bubbles	The aim is to ensure a homogeneous distilate because, the determination of gravity will proceed from this.

⁵ The flasks are calibrated and interchangeable in mass, not in volume. Remember that our flask explicitly DO NOT indicate ml, but only 100-200. We DO NOT withdraw 100 ml, but a repeated quantity. So accuracy in ml is NOT necessary, but it is imperative that the volume of wine and distillate is the same at the same temperature.

⁶ The exact quantity depends on the sample.

4.2 THE DISTILLATION PROCEDURE FOR SPIRITS

Spirits present different characteristics from those of wines:

- Values of ASV up to 70%
- Variable composition

In fact spirits with very different composition exist: from cream whisky to fruit creams, from very sweet liqueurs such as vodka with fruit or limoncello to liqueurs with extracts of aromatic plants or roots, or even drinks with egg or coffee.

In addition, the higher alcoholic strength results in:

- a very slow distillation in order to make sure of the complete extraction of alcohol;
- collection of a greater volume of distillate (around 95 ml).

It is difficult to find a common distillation method for all the different kinds of spirits because the various ingredients interfere with the distillation in different ways, sometimes blocking the passage of the alcohol. The basic procedure for spirits is listed below. This procedure allows carrying out a good distillation in a short time, without excessive foam, and resulting in a clear distillate; the specific methods applicable to different alcohol products are indicated in the subsequent paragraph.

The distillation procedure	Observations
1. The sample is poured up to the mark of the special flask ⁷	It is recommended that the temperature of the sample is 20 °C. One can also operate at room temperature: a maximum deviation of ± 2 °C is allowed between the sampling and the bringing back to volume after distillation.
2. Pour the contents into the bulb of the distillatory	
3. Wash the flask 3 times with 10 ml of distilled water and add washings to distillatory bulb	
4. Add the “ solution for liqueurs ” ⁸ : solution 1:1 of NaCl 1% and [KAl(SO ₄) ₂] 12%	The aim is to increase the conductivity of the sample and to suppress the aromatic substances.
5. Add a few ml of diluted antifoam	Its addition avoids the formation of foam which would be an obstacle to the normal distillation of the alcohol.
6. Add about 2 ml cold distilled water to the same flask used for measurement of the sample. Collect \approx 95 ml of distillate therein.	
7. Add distilled water nearly up to volume, stopper and bring to volume when the temperature reaches 20 °C or the temperature of sample before distillation	USE OF A THERMOSTATIC BATH IS ADVISED
8. Stopper the flask and mix with care to avoid bubbles	

Regarding brandy, etc. the method is the same as the one indicated for wines, with the exception of the addition of a small quantity of lime milk: because of the low acidity of such products, 1-2 ml is sufficient. In case of alcoholic strength higher than 50 % vol, it is advisable to dilute to 1:2

⁷ See note 1

⁸ For the quantity see Table 1

4.3 COMPARISON BETWEEN THE SUPER DEE SV AND THE OIV GLASS DISTILLING UNIT

Tests of distillation with spirits with extremely varied compositions were carried out. Each liqueur was distilled :

- with the SUPER DEE SV distilling unit following the above procedure;
- with a glass distilling unit using the official procedures of OIV.

This latter distilling unit satisfies the official OIV condition.

The results obtained with the two apparatus are shown in the **TABLE 1**.

A specific method was developed to be used with the SUPER DEE SV distilling unit for each different sample type. Generally it is recommended:

- very sweet liqueurs: use small quantities of distilled water to wash the flask (less than 10 ml); in many cases this solution blocks the formation of emulsions which can result in an incorrect distillation;
- liqueurs containing milk or cream: the solution for liqueurs should not be added (it may curdle the milk creating a clot that is difficult to clean out of the flask).

The results obtained with the use of the SUPER DEE SV distilling unit are without doubt comparable with those obtained with the OIV glass distilling unit, with the advantage of much shorter distillation times. The differences found from the experimental data on some samples compared to the values indicated on the label fall within the tolerance expected by the current legislation (0,3 % vol on indicated ASV).

For further information see the book “*Analysis methods of wines and alcoholic beverages*”, Gibertini – current Edition.

TABLE 1

SAMPLE	INGREDIENTS	DEE PROCEDURE	%vol DEE	Declared
Liqueur type Strega	Alcohol water, sugar, natural flavours, saffron	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	40.37	40
Lemon cream	Water, sugar, lemon infusion, milk, cream, alcohol, glucose syrup	100 ml +12 ml antifoam	22.05	22
Rakee of raspberries	rakee, raspberries juice, sugar, natural flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	32.02	32
Limoncello	Alcohol, sugar, infusion of lemon rind, natural flavours	100 ml + 3 ml antifoam +10 ml milk of lime	30.29	30
Whisky cream	Dairy cream, sugar, Whisky, water, natural flavours	100 ml +15 ml antifoam	14.71	15
Rakee	Rakee	100 ml +5 drop antifoam +3 ml solution for liqueurs	40.12	40
Pina Colada	sugar, alcohol, cream, rum, milk, flavours	100 ml + 9 ml antifoam	16.80	17
China Martini	Sugar, alcohol, flavours, colorants: caramel; acidifier: citric acid	100 ml +15 ml antifoam + 10 ml milk of lime	31.00	31

(continue)

SAMPLE	INGREDIENTS	DEE PROCEDURE	%vol DEE	Declared
White Vermuth	Wine, sugar, alcohol, flavours	100 ml +1.5 ml antifoam + 5 ml milk of lime	16.06	16
Strong anisette	Alcohol, sugar, natural flavours	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	40.02	40
Liqueur with coffee	sugar, alcohol, coffee, natural flavours	100 ml + 12 ml antifoam	32.89	33
Cherry	Alcohol, sugar, egriot juice, natural flavours	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	30.03	30
Tequila	Distillate of agave techilana	100 ml +1.5 ml antifoam +2 ml solution for liqueurs	38.14	38
Pastis	Alcohol, sugar, 2 ⁰ / ₀₀ anetolo, infusion of vegetable substances, natural flavours,	100 ml +1.5 ml antifoam +5 ml solution for liqueurs	44.76	45
Banana Vodka	vodka, sugar, banana juice, flavours	100 ml +1.5 ml antifoam +3 ml solution for liqueurs	25.17	25
Brandy	Alcohol hydrate, Brandy, natural flavours, colorant: caramel	100 ml + 3 ml antifoam +2 ml solution for liqueurs	36.07	36
Liqueur with chocolate and mint	Milk, sugar, cream, glucose syrup, alcohol, milk protein, infusion of cacao, flavours, colorant: caramel	100 ml + 12 ml antifoam	16.06	16
Cachaca	Distillate of sugar-cane	100 ml + 3 ml antifoam +3 ml solution for liqueurs	43.66	43
Amaretto (alcoholized for exclusive use in confectionery)	Alcohol, caramel, flavours	25 ml + 75 ml distilled water (collect 100 ml) + 3 ml antifoam +3 ml solution for liqueurs	% vol value × 4 = 70.00	70
Sake	Distillate of rice	100 ml + 1 ml antifoam +10 ml solut. for liqueurs	15.85	16
Sambuca	Alcohol, sugars, anicette, natural flavours	100 ml + 5 ml antifoam +8 ml solution for liqueurs	39.86	40
Cointreau	Alcohol, sugars, orange flavours	100 ml + 5 ml antifoam +7 ml solution for liqueurs	39,81	40

NOTE

The samples are only indicative, so they must not be considered exhaustive.

5 ANALYSIS – Example of distillation

5.1 ANALYSIS SELECTION

Put wine or liqueur sample into the distillation chamber and establish the quantity of the distillate by pushing the “**MODE**” key.

Select one of the 3 set points (80, 160, 240) through the “**A**”, “**B**” or “**C**” keys, respectively.

It is possible to modify the selected value by pushing together “**SHIFT**” and “**A**” (+) or “**B**” (-) once or more times as to increase or decrease one or more units.

Figure 3

5.2 DISTILLATION START

Press “**MODE**” and the “**A**” key to begin the distillation.


Figure 4

WARNING

Before starting the distillation check that the plug is indeed closed (turn clockwise the knob), the flask is installed on the balance and the protection is completely lowered. If there aren't these conditions the distillation don't start and an acoustic signal informs of the anomaly.

Example of distillation:

- plug and drainage closed
- flask on the balance
- protection down.

Figure 5

5.3 DISTILLATION CHAMBER CLEANING

After the distillation open the drainage and clean the distillation chamber with the shower supplied with SUPER DEE SV.

Figure 6

When the SUPER DEE SV distilling unit is not in use leave the drainage open.

Figure 7

6 ERROR MESSAGES

There are indications for the main anomalies. The display shows the messages with the operative situations and the lights **WARNING** and **ALARM**.

Remember that, for a correct operation, the system must have a capacity of 5 litres/min. and that the water must have a temperature of 20°C max and a pressure of at least 2 bar.

- 1) In case of NO water or low water pressure:

Error !
Water press. low

THE INSTRUMENT DOES NOT START (Fig. 8)

- 2) After 10 seconds the display shows the message

No water

THE INSTRUMENT DOES NOT START

If the error situation continues the distillation will be cancelled and the program will wait to confirm the reading of the message by operator.

No water
Ok

CONFIRM THE READING OF THE MESSAGE.

- 3) In case of cooling water too hot the yellow light flashes. Push the START key to start a distillation. The instrument will wait 10 seconds and if the error situation remains the operator will select if to continue or cancel the distillation.

Water too hot
Go on Cancel

CONFIRM OR CANCEL THE DISTILLATION (Fig. 9)

ATTENTION!
The temperature of the cooling water upper 20 °C compromises the result of the distillation


Figure 8


Figure 9

Many situations of error have the indication of the error message until the restored normal situation.

- 4) In case of drainage open the following message will appear:

Attention!
Drainage open

*THE DISTILLATION DON'T START.
THE INSTRUMENT WAITS FOR THE CLOSURE
OF THE DRAINAGE (Fig. 10)*

- 5) In case of open plexiglass the following message will appear:

Attention!
Plexiglass open

*THE DISTILLATION DOES NOT START.
THE INSTRUMENT WAIT FOR THE CLOSURE
OF THE PLEXIGLASS (Fig. 11)*


Figure 10


Figure 11

- 6) If the flask is not positioned correctly, errors in the quantity of distillate collected may occur. Position the flask as shown in the figure on the side, so as to avoid contact between the flask and the glass tube for the collection of the distillate.

7 MAINTENANCE

7.1 DISTILLATION CHAMBER

After each distillation the distillation chamber must be washed with water to prevent contamination of the successive sample.

Periodically clean the distillation chamber introducing approximately 500 ml of distilled or demineralised water and approximately 20 ml of sodium chloride NaCl 1% solution. Start the distillation.

To remove deposits or incrustations that the samples can have left on the electrodes and/or on the glass can be done in the two following ways:

- Pour into the distillation chamber a dilute solution of HCl 0,1N. Await some hour and eventually use a tube-brush to clean the chamber.

or

- Introduce into the distillation chamber the solution B, supplied with the instrument, in sufficient quantity to cover the incrustations. Also in this case to await some hour and use eventually a tube-brush.

In any case rinse always accurately the distillation chamber before reuse the instrument.

WARNING! Don't press the START button when into the distillation chamber there is a cleaning solution. For more safety remove the power socket from the main.

7.2 ELECTRODES

The electrodes will gradually diminish in size and will have to be replaced at some stage. It is difficult to predict the frequency of replacement because it will depend on total usage and type of sample (approx. 2000 tests).

For the replacement call the SAT (Service) of GIBERTINI - Phone (+39) 023541434 or your supplier.

7.3 RECTIFICATION COLUMN

During the normal use of the SUPER DEE, the column doesn't require a particular cleanness. However, if the column is contaminated by acid residues, it is necessary to clean it following indication reported on procedure 70.30.136 Rev.1 available as appendix of this user manual

7.4 EMPTYING FUNCTION

- Switch off the DISTILLATOR with the ON/OFF SWITCH.
 - Switch on pressing simultaneously the **B** key on the keyboard.
 - The display shows the message "**Take off water tube Wait for tank drain**"
 - Empty the cooling coil using the shower.
 - Wait for the complete emptying of the distilled water tank by means of the steam pump (the water finishes to go out by the nozzle and the pump becomes very noisy). Switch off the instrument.
 - Incline the instrument backwards to help the exit of water.
-

8 OPTIONAL CONNECTION WITH “THERMOCHILLER”

In some operating environments, it may happen that the network water, used as a heat exchanger in the refrigerant, is available at temperatures that are not low enough to allow adequate cooling, reducing the temperature difference necessary to ensure effective distillation.

In these cases the use of the Gibertini liquid circulation cooler mod. Thermochiller (or other suitable cooler) is strongly recommended.

8.1 DISTILLER INTENDED FOR USE WITH “THERMOCHILLER” (FACTORY SET)

If the distiller is purchased with the Thermochiller as accessory, the instrument is set at the factory with the appropriate electronic configuration and complete with connections, pipes and accessories necessary for use with the cooler.

Compared to the standard version, the electronic display, every time a distillation is started, displays the warning message for a few seconds for the preparation and ignition of the relative Thermochiller and then the distillation itself is performed.

ATTENTION: in this configuration, if Thermochiller is not activated, no liquid flow through the refrigerant cooling circuit and the distillation will not be successful, with possible overheating of the system.

8.2 DISTILLER INTENDED FOR USE WITH TAP WATER

If the distiller is subsequently equipped with the Thermochiller accessory, the instrument will not have the appropriate electronic setting for use with the cooler.

To configure the instrument for use with ThermoChiller proceed as follows:

- turn off the instrument and restart it later by holding down the "C" key
- the chiller use selection screen will be displayed
- confirm by pressing the YES key
- wait for the normal restart of the instrument

After that the distiller will operate in the "Chiller" configuration and must be used only with the relative cooler in operation.

To return to the mains water cooling configuration, carry out the above procedure, selecting "NO" in the chiller use selection screen.

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CONNECT TO THERMOCHILLER **OUTLET** NIPPLE

APPENDIX1: Procedure 70.30.136 How to clean the distillation chamber and rectification column

	CLEANING MODE	WHEN TO DO IT	PROCEDURE
1-	Standard cleaning	Periodically, after a series of VOLATILE distillations.	<ul style="list-style-type: none"> a. Pour in the distillation chamber about 20ml of distilled water, and perform distillation in VOLATILE mode. b. Repeat once or twice.
2-	Cleaning of rectification column 1st mode	When the column appears to be polluted by acid residuals (i.e. when the distillate of water requires more than 0.1ml of NaOH 0.1N to be titrated).	<ul style="list-style-type: none"> a. Fill the distillation chamber with about 500 ml of distilled water + about 20 ml of NaCl 1% and perform a distillation in ALCOHOL mode. <p>N.B.- This type of washing performs also a mechanical action on the column, since the water goes up from the distillation chamber into the column.</p>
3-	Cleaning of rectification column 2nd mode	When the column filling is made black by the products of the distillations, or when acidity persists during volatile distillations (the distillate of water requires more than 0.1ml of NaOH 0.1N to be titrated).	<ul style="list-style-type: none"> a. Pour in the distillation chamber about 1ml of a bleach (sodium hypochlorite), add about 500 ml of distilled water and perform a distillation in ALCOHOL mode. b. Repeat cleaning as indicated in point 2 (1st mode). c. Pour in the distillation chamber about 20 ml of distilled water, and perform distillation in VOLATILE mode. d. Repeat distillation of water until the water distillate requires less than 0.1ml of NaOH 0.1N to be titrated.
4-	Cleaning of distillation chamber	When the distillation chamber appears dirty or there are residual inside.	<ul style="list-style-type: none"> a. Pour in the distillation chamber about 100 ml of bleach, and let it rest for about 30 minutes. <p>DO NOT DISTILL !!!</p> <ul style="list-style-type: none"> b. Drain and wash. c. Pour in the distillation chamber about 20 ml of distilled water, and perform distillation in VOLATILE mode. d. Repeat distillation of water until the water distillate requires less than 0.1ml of NaOH 0.1N to be titrated.

customer care

info@gibertini.com

Tel. (++39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI

ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (++39) 02 3541434

Fax (++39) 02 3541438

www.gibertini.com

sales@gibertini.com