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This manual was written and produced with the best care. However, errors cannot be fully excluded. GIBERTINI ELETTRONICA Srl does not take any responsibility and accepts no liabilities of any kind that may occur because of errors in the manual.

GIBERTINI ELETTRONICA Srl reserves the right to change the contents and/or the informations contained in this manual without advice, at any time without notice.

All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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1.	DESCRIPTION	20
1.1	TECHNICAL FEATURES	20
1.2	PRINCIPLE OF OPERATION	21
2.	INSTALLATION	21
3.	USE	23
3.1	MAIN MENU	23
3.2	INPUT A STRING	23
3.3	ENTERING A NUMBER	23
3.4	FLOW CELL ASPIRATION	24
3.5	FLOW CELL WASHING	24
3.6	SELECT METHOD	24
4.	ABS READING	25
4.1	READING BY CUVETTE	25
4.2	READING BY FLOW CELL	26
5.	MODIFY METHOD.....	26
6.	EXECUTE METHOD.....	28
6.1	FLOW CELL REQUEST	28
6.2	TEMPERATURE CHECK	28
6.3	FLOW CELL WASHING REQUEST	28
6.4	LAMP CHECK AND WATER READING	28
6.5	CHECK REAGENTS	29
6.6	STANDARD READING	29
7.	END-POINT METHOD	30
8.	KINETIC METHOD	30
9.	FIXED-TIME METHOD	30
10.	SETTING	31
10.1	LANGUAGE	31
10.2	FLOW CELL SETTING	31
10.2.1	ASPIRATION	31
10.2.2	WASHING	32
10.2.3	CALIBRATION PUMP	32
11.	PRINTER	33
12.	STORED RESULTS	34
13.	QUALITY CONTROL	35
13.1	QUALITY CONTROL EXECUTION	35
13.2	Q.C. STATISTIC DATA	35
14.	PERFORMANCES OF REAGENTS FOR WINEMATIC.....	37

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1. DESCRIPTION

1.1 TECHNICAL FEATURES


WINEMATIC is a photometer used for biochemical tests.

Light source	Digitally controlled Halogen lamp (12V - 20Watt)
Wavelength	8 standard interference filters: 340nm, 420nm, 492nm, 505nm, 520nm, 546nm, 578nm, 620nm + led 280nm.
Photometric Range	0 - 2.5 O.D (led 280nm 0 - 1.6 O.D)
Flow cell system	Flow cell with 10mm light path, interchangeable with disposable macro, semi-micro, or special optical glass cuvettes.
Delay time	4 sec.
Reaction volume	900 µL per test min.
Temperature control	Peltier elements, between 15°C - 45°C.
Incubator	8 positions.
Display	Graphic, 240x128 pixel.
Printer	Thermal graphic, 24 characters/line.
Communications	USB port.
Dimensions	28 x 36 x 22 cm (L x W x H)
Power supply	100÷240 VAC 50/60 Hz With special power cable EMC compatible.

1.2 PRINCIPLE OF OPERATION

WINEMATIC is a photometer that measures the concentration of a substance in a liquid via optical principles. Through the implementation of methods designed for the particular substance it is possible to reproduce chemical reactions that turn into a solute capable of absorbing light of a specific wavelength. The physical law that determines the concentration according to absorption of light is known as the Lambert-Beer law:

$$C = k * \log\left(\frac{I_0}{I_1}\right)$$

Where:

- I_0 è the intensity of light incident.
- I_1 è the intensity of refracted light.
- k is a constant that depends on physical parameters such as temperature T , the optical path l , the wavelength of incident light λ .
- C is the concentration of solute.




The quantity $\log\left(\frac{I_0}{I_1}\right)$ is called as *absorbance (ABS)* or *optical density (OD)*.

It is proportional to the concentration by k . WINEMATIC directly measure the absorbance and calculate the factor k using reagents and calibrators. This makes it possible to determine the concentration C .

2. INSTALLATION

The installation of the instrument must be performed by a qualified technician. In any case, the following precautions should always be performed:


- Place the instrument into a dry and clean site.
- Adjust the instrument so that it is slightly tilted towards the operator
- Keeping clean the instrument, especially the front suction tube to avoid contamination between samples and between examinations.
- **Use only the power cord supplied.** It should be noted that the instrument has internal self-fuse.
- Make sure the power supply is connected to an electrical ground .
- Connect the waste connector on rear side to a bottle using the tube.
- Check O.D. linearity by using periodically standard calibrators.
- If the measured absorbance is greater than 2.500 ABS, repeat the reading or dilute the sample
- Check periodically the regular calibration of the peristaltic pump.
- Keep clean and dry the surfaces of the flow cell.
- Wait at least 25 minutes after switching the instrument to make sure it has reached thermal equilibrium.
- Do not expose the instrument to heat and keep it away from direct sunlight.

3. USE

3.1 MAIN MENU

The main screen has six items:

With the keyboard it is possible to select following functions, pushing the number on the left:

1. **EXECUTE METHOD:** Allow to execute a method. Push this section to start immediately with the execution (see 3.6).
2. **MODIFY METHOD:** Allow to modify a method or create a new method. . Push this section to start immediately with the modifications (see 3.6).
3. **STORED RESULTS:** Allows to see the archive of the last 15 stored results.
4. **READ ABS:** Allow to perform an absorbance reading manually. This can be accomplished in two ways: pressing the key **4** is accessed for reading through disposable cuvette, pressing the **GO** button accesses the readings through flow cell.
5. **SETTINGS:** Allow to access on the setting menu.
6. **FLOW CELL WASH:** Allow to perform a flow cell washing.

3.2 INPUT A STRING

The input of a string is done by inserting a character by time.

- Using the number keys to each key has an associated set of characters. Pressing the same button repeatedly will display these characters in sequence. Pressing a button other than the one preceding it automatically switches to the next character. To confirm the entry value push **ENTER** press once or wait one second.
- To cancel the entry value push **DEL** and to finish without entering a string, press **ESC**.


3.3 ENTERING A NUMBER

The release of a number is via the numeric keys. To confirm using **ENTER**, use **ESC** to delete and to exit without confirming.

3.4 FLOW CELL ASPIRATION

The aspiration of the liquid through the flow cell follow these steps:

- **LIQUID REQUEST:** Tool buzzer alert the operator to immerse the suction tube into the liquid. The open end of the tube should always draw the liquid into the lower part of the container. The figure shows an example of a cuvette.
- Once the tube immersed in the liquid, press the **GO** button to start the aspiration
- **LIQUID REMOVAL:** Over the aspiration of the liquid (depending on the settings) the instrument alerts the operator via buzzer and wait three seconds to allow the removal of the tube from the liquid.
- **AIR GAP CREATION:** Allows, when enabled, creating a gap of air in through a suction tube into the air.

3.5 FLOW CELL WASHING

The flow cell washing follow these steps:

- **WASHING SOLUTION REQUEST:** The tool alerts via display to dip the tube in the washing solution. Apply the requirements used for aspiration. Once the tube immersed in the liquid, press the **GO** button to start the aspiration.
- **WASHING SOLUTION REMOVAL:** Over the aspiration of the liquid (depending on the settings) the instrument alerts the operator via buzzer and wait three seconds to allow the removal of the tube from the liquid.
- **PAUSE:** Where washing is started manually from the main menu tool performs a break (depending on settings). During the execution of a method of washing is performed without pause.
- **WASTE:** The instrument empty the liquid from the tube.

3.6 SELECT METHOD

Both the function 1-EXECUTE METHOD that the function 2-MODIFY METHOD main menu requires the selection of a method to execute or edit. It is via the menu SELECT METHOD:

1-EXECUTE METHOD

SELECT METHOD:

1-SLOT0

2-SLOT1

3-SLOT2

4-SLOT3

5-SLOT4

F1 – SEARCH ID:

NEXT →

2-MODIFY METHOD

SELECT METHOD:

1-SLOT0

2-SLOT1

3-SLOT2

4-SLOT3

5-SLOT4

F1 – SEARCH ID:

PRESS +/- TO ADD METHODS

NEXT →

The instrument has available 100 memory slots each capable of containing a method and numbered from 1 to 100. The list of methods is organized by five method by page, through ← and → button you go from page to page.

The selection of the method may occur in two ways:

- **QUICK SELECTION:** The selection is done via keyboard by pressing the highlighted number (from **1** to **5**) to the left of each item.
- **SELECTION BY METHOD ID:** Selection by ID is done by pressing the **F1** key directly from the keyboard and then entering the slot number where you stored the method then press **ENTER**. Pressing the **PRINT** button you can print the list of the methods previously stored. Is recommended to print the list to quickly locate to the slot where you stored the selected method. If you want to enter a new method is necessary to find a free memory slot and select it. To return on the main menu press **ESC**.

Pressing the **PRINT** button you can print the list of the methods previously stored. Is recommended to print the list to quickly locate to the slot where you stored the selected method. If you want to enter a new method is necessary to find a free memory slot and select it. To return on the main menu press **ESC**.

4. ABS READING

ATTENTION: to work with IPT method or with 280nm led remove the flow cell and use quartz or disposable cuvettes, use only disposable cuvette which do not absorb at 280 nm.

4.1 READING BY CUVETTE

Allows to measure the absorbance of a liquid contained in the cuvette from the main menu by selecting the function 4-READ ABS then select the **F2** key. In this mode, you can not use the cell to flow.

The screen is:

ABSORBANCE TEST	
1 – FILTER:	340nm
2 – CELL TEMPERATURE	37.0°C
3 – TIMER: 0 sec.	
ABS: 0.000	
F1 - ZERO	F2 - ABS


The options are:

1. **FILTER:** Repeatedly pressing the **1** button to select the wavelength at which to read. **Each time you change the filter you need to read again the zero reference absorbance.**
2. **CELL TEMPERATURE:** Pressing **2** button you can change the temperature of the cell by entering a new reading. On the same line you can read the current temperature of the cell.
3. **TIMER:** The tool allows the operator setting a timer in seconds. Pressing **3** button is likely to enter a number of seconds and then immediately placing the countdown. At the end of the counting instrument alerts the operator via buzzer. Pressing any button turns off the buzzer. The use of other features stops the countdown.
4. **F1-ZERO:** Pressing the **F1** key the instrument reads the absorbance, the set as zero reference and shows **0.000**.
5. **F2-ABS:** Pressing the **F2** key tool reads the absorbance and displays. Before reading the absorbance of a sample is always necessary to reset the instrument. If this was not done the tool alerts you both on screen and with the buzzer. If you change the wavelength zero is cancelled and the instrument requires re-reading.

To print the reading press **PRINT**. To return on the main menu press **ESC**.

4.2 READING BY FLOW CELL

Allows to measure the absorbance of a liquid aspirated by peristaltic pump. It is accessible directly from the main menu by pressing the **GO** button or selecting from the main menu function 4-READ ABS then selecting the **F1** key. The screen is:

The options are:

1. **FILTER:** Repeatedly pressing the **1** button to select the wavelength at which to read. **Each time you change the filter you need to read again the zero reference absorbance.**
2. **CELL TEMPERATURE:** Pressing **2** button you can change the temperature of the cell by entering a new reading. On the same line you can read the current temperature of the cell.
3. **TIMER:** The tool allows the operator setting a timer in seconds. Pressing **3** button is likely to enter a number of seconds and then immediately placing the countdown. At the end of the counting instrument alerts the operator via buzzer. Pressing any button turns off the buzzer. The use of other features stops the countdown.
4. **WASH FLOW CELL:** Pressing **4** button to wash the flow cell.
5. **F1-ZERO:** Pressing the **F1** key the instrument set as zero reference the last and shows **0.000**.
6. **GO:** Pressing the **GO** key the instrument start the sample aspiration. At the end of the aspiration appears on the screen the reading absorbance value. Once finished the flow cell is being empty.

To print the reading press **PRINT**. To return on the main menu press **ESC**.

5. MODIFY METHOD

To access the modify or the inclusion of a new method choose, from the main menu, the function 2-MODIFY METHOD. To modify a method type the ID or selecting an empty slot to create a new one. In both cases, the input parameters procedure is the same. For this reason we will refer only to the release of a new method.

The input of the parameters is done by pages. Each page contains a number of parameters and the sequence of pages is fixed.

To move from page to page using the buttons **←** and **→**. To select a parameter and change the value press the key repeatedly referred to the numerical parameter or to the release of an alphanumeric string (eg. method name, sample volume, reagent volume, ... qc) follow the procedure described in section 3.2. When the operator has completed the changes, press **ESC**. It appears a screen that allows you to save your changes. If you decide to save press **F1** key, otherwise, if you decide not to save you press **F2** key and in this case all changes are lost.

Any screen you can press the **DEL** key to delete the method.

Depending on the method set some options and some parameters may be present or not. Follows the table of adjustable parameters for each type of programmable methods:

PARAMETER	ENDPOINT	FIXED TIME	KINETIC	NOTE
ID	YES	YES	YES	ID linked to the method
NAME	YES	YES	YES	Name linked to the method
TYPE	END POINT	FIXED TIME	KINETIC	Type of the method
ZERO	WATER REAGENT SAMPLE	WATER REAGENT	WATER REAGENTE	Typo of zero
UNITS	See table	See table	See table	Measure unit
TEMPER. [°C]	YES	YES	YES	Execution temperature method
CALIB.	FACTOR STANDARD	FACTOR STANDARD	FACTOR STANDARD	Type of factor (k) calculation
DECIMALS	0 – 3	0 – 3	0 – 3	Number of decimals
SAMPLE[μL]	YES	YES	YES	Sample volume
REAGENT(1)[μL]	YES	YES	YES	Reagent volume
REAGENT(2)[μL]	NO	NO	YES	Second reagent volume
FILTER(1)	YES	YES	YES	Primary wavelength
FILTER(2)	YES	NO	NO	Secondary wavelength (bichromatic method)
WEIGHT FILTER 2	0 – 1	NO	NO	Secondary wavelength
K-FACTOR	YES	YES	YES	Factor (k)
OFFSET	YES	NO	NO	Offset
STD NUMBER	1 – 8	1	1	Number of standard
STD#n	YES	YES only STD#1	YES only STD#1	"n" value standard concentration
FITTING	LINEAR REGR LINEAR P-W	NO	NO	Type of approximation for multistandard methods
READ. TIME	NO	YES	NO	Time reading for Fixed-Time
DIRECTION	NO	INCR DECR	INCR DECR	Expected direction of the reaction
START dABS	NO	YES	YES	Expected direction edge of initial absorbance
INCUBATION(1)[sec]	NO	YES	YES	Incubation time
INCUBATION(2)[sec]	NO	NO	YES	Incubation time of second reagent (only if REAGENT2≠0)
REAG. BIAS	NO	YES	NO	Select YES to include in calculating the variation of the reagent.
READINGS #	NO	NO	YES	Number of readings to be performed
DELTA TIME [sec]	NO	NO	YES	Time between reading and the subsequent
NORMAL MAX	YES	YES	YES	Maximum expected concentration
NORMAL MIN	YES	YES	YES	Minimum expected concentration
LINEARITY MAX	YES	YES	YES	Concentration above which the method is no linear anymore
REAG. MAX ABS	YES	YES	YES	Maximum expected absorbance of the reagent
REAG. MIN ABS	YES	YES	YES	Minimum expected absorbance of the reagent
FIT MIN	NO	NO	0 – 1	Minimum FIT for interpolation curve

6. EXECUTE METHOD

To execute a method choose, from the main menu, the function 1-EXECUTE METHOD and select the required method (see paragraph 3.6). The execution of a method is in the form of assisted procedure: the instrument through a series of screens guide the operator during the execution with appropriate requests for information. Initially run a series of preliminary steps common to every type of method (described below). According to the print settings selected, the instrument may or may not print the method before executing.

6.1 FLOW CELL REQUEST

ATTENTION: to work with IPT method or with 280nm led remove the flow cell and use quartz or disposable cuvettes, use only disposable cuvette which do not absorb at 280 nm.

The first request concerns the flow cell. If the operator intends to work with the flow cell must respond YES by pressing the **F1** key or the **GO** button, or whether it will work with cuvette must answer NO pressing **F2** key. Pressing **ESC** button is aborting the execution.

If the operator has chosen to work with flow cell the **GO** button is activated and can continue execution of the method even in subsequent requests. If the operator has chosen not to use the flow cell the **GO** button is disabled. Depending on whether or not you use the cell flow, the instrument will run or not the request of the aspirations and washing the sample and reagents.

USE FLOW CELL ?

F1/GO-YES F2-NO

6.2 TEMPERATURE CHECK

The instrument checks that the temperature is that one required by the method waiting until the cell reaches the proper temperature reading. Alternatively, you can skip checking by ignoring the waiting, the temperature will come to the setting value but the operator doesn't have the check.

METODO: NAME
ID: ID N°

PLEASE WAIT TEMPERATURE
TEMPERATURE
37.0 °C

F1/GO-IGNORE ESC-ABORT

6.3 FLOW CELL WASHING REQUEST

The instrument asks whether or not to wash the flow cell. If the operator wishes to wash the flow cell must respond YES by pressing the **F1** key, otherwise, if you are not washing the flow cell must answer NO pressing **F2** or the **GO** button.

6.4 LAMP CHECK AND WATER READING

Initially, the instrument requires the reading of water. This reading is used to check the lamp and, as appropriate, such as zero reference.

6.5 CHECK REAGENTS

According to the settings of the method, the instrument can ask whether the reading of absorbance of the reagents.

If in the running method they were specified as zero reference then the instrument performs the reading of the absorption, otherwise the reading is done only if you specified the parameters of minimum or maximum absorbance of the same. Control absorbance is read and displayed on the screen where the absorbance does not fall within the range defined in the method is reported the error.

REAGENT ABS R2: 0.130

F1 – REPEAT

F2/GO - CONTINUE

According to the settings of the method, the instrument will store the values of absorbance of the reagents, asking whether to use the values saved or re-reading.

6.6 STANDARD READING

If the method require a calibration and calibration has never run for the selected method, the instrument requires reading in sequence starting from the standard number 1. Once read you can accept it (**F1/GO**), repeat (**F2**) or terminate (**ESC**) the execution of the method.

METHOD: NAME

ABS STD#n: x.xxx


F1/GO-CONTINUE F2-REPEAT ESC-ABORT

During the reading of standard the instrument will print all the readed absorbances by each concentrations specified by the method. In case of a multi-standard method, after the reading shows a graph of the readed absorbances and interpolate the curve that can be printed by pressing **PRINT**. In addition to the chart is printed a list of standards and the fit calculated by the instrument.

In the case that the instrument has already made a calibration of the method selected, the operator will be asked whether to repeat the calibration or hold the last calibration performed.


7. END-POINT METHOD

An END-POINT method is used to determine the concentration of a substrate. It consists of a reaction between the substrate that is to be determined and a substance known as chromogen, whose products are substances that are relatively stable over time. Usually this reaction is catalyzed by one or more enzymes. In this type of reaction the enzyme is always in excess. After the addition of sample to the reagents, reaction comes to completion after a fixed period of time. At this point the reaction products are stable and it is possible to measure the absorbance. The concentration value obtained by multiplying the absorbance value measured by the factor (k) previously calculated based on data provided by methods and / or obtained by calibration using a standard.

8. KINETIC METHOD

A kinetic method is used to determine the activity of an enzyme. It consists of a reaction where the substrate concentration is constant on the reagent and the speed of reaction is proportional to the concentration of the enzyme in the sample. By measuring the speed of reaction (i.e. speed with which changes the concentration of substrate and / or products over time) it goes back to the activity of the enzyme in solution. The instrument measures the absorbance change at regular intervals of time, then calculate the speed to obtain the concentration of the enzyme.


9. FIXED-TIME METHOD

A fixed-time method is a non-linear end-point method, with only two meaningful readings. The difference between these two measurements versus time allows to calculate the concentrations of interest and thus the result of the analysis.

10. SETTING

5-SETTING menu allows you to set the operating parameters of the instrument:

The available functions are:

1. **LANGUAGE:** This feature allows you to select your language ITALIAN or ENGLISH
2. **FLOW CELL SETTINGS:** This feature allows to access to the configuration menu of the flow cell to adjust the suction volume (aspiration setting), washing volume (washing setting) and the calibration (pump calibration) of the peristaltic pump.
3. **CLOCK ADJUST:** This feature allows to adjust the date and time.
4. **PRINTER:** This feature allows to modify the printer settings.
5. **LAMP CALIBRATION:** This feature allows to calibrate the lamp energy.
6. **ADVANCED:** This feature allows to perform maintenance operations such as calibration lamp or thermostats. It is password protected and should be used only by qualified service personnel.

10.1 LANGUAGE

Access to select the English or Italian language.

10.2 FLOW CELL SETTING

Through the features offered by this menu you can perform the calibration of the peristaltic pump and the configuration of the volumes of aspiration.

10.2.1 ASPIRATION

The function 1-ASPIRATION SETTINGS configures the suction volume and air-gap option. The suction procedure allows to perform two consecutive aspirations. The first concerns the suction of the liquid to be analysed while the second (optional) concern in the air-gap and allows the creation of a volume of separation between two consecutive aspirations which contributes to avoid contamination between successive samples.

ASPIRATION SETTINGS

1-VOLUME [uL]: 900
2-CREATE AIR GAP: ON
3-AIR GAP VOLUME [uL]: 250
DEFAULT: 900/ON/250
450/ON/150

The possible options are:

1. **VOLUME:** Volume of sample to aspirate (uL micro-liter).
2. **CREATE AIR GAP:** If it is ON abilitate the second aspiration, if it is OFF disabilitate.
3. **AIR GAP VOLUME:** Volume of second aspiration (uL micro-liter).

It is advisable to avoid using suction volumes lower than the recommended default.

10.2.2 WASHING

The function 2-WASHING SETTINGS concern to configure the volume of washing liquid inlet and the waiting time, that means, the time during which the washing liquid remains within the hydraulic system. This pause occurs only if the wash is called from menu: during the execution of a method is never executed.

WASHING SETTINGS

1-VOLUME [uL]: 1000
2-TIME [SEC]: 10
DEFAULT: 1000uL/10sec.

The possible options are:

1. **VOLUME:** Volume of liquid to aspirate (uL micro-liter).
2. **TIME:** Waiting time (seconds).

10.2.3 CALIBRATION PUMP

The function 3-PUMP CALIBRATION allows to calibrate the peristaltic pump by linking the rotation of the pump with the volume aspirated. It is necessary to calibrate the pump at the moment of installation of the instrument. A regularly calibration of the pump is essential because it allows to recover the errors that arise over time in the event of excessive wear of the tube. It is recommended to perform regular calibration in case of intensive use of the instrument.

PUMP CALIBRATION

1-ASPIRATED VOLUME [uL]: 1000
2-STEPS COUNTED: 6000
3-SAVE CALIBRATION
CALCULATE uL/STEPS: 0.167
F1-TEST ASPIRATION
F2-TEST WASHING
PRESS "GO" BUTTON TO GENERATE

The peristaltic pump is equipped with a stepping motor. Each step of the engine is a volume suction of the pump. The instrument calculates the number of steps that the engine must run to aspirate the required volume and activates the motor according to this calculation. Each time the operator press the **GO** button the instrument activates the pump and gives to the video a number of step generated so that the operator can know the number of steps needed to aspirate a certain volume of liquid. The calibration procedure is as follows:

1. Prepare a tube with a known volume of distilled water to start the calibration. For example 3000 uL.
2. Insert the chosen volume in the instrument by option 1.
3. Insert the aspiration tube into the tube with water and hold the **GO** button until the instrument aspirate all the volume of water contained into the tube, avoiding to aspirate over this period.
4. Read on screen the number of steps generated under STEPS COUNTED and place it into the instrument through option 2.
5. If the operator are sure to follow the above operations correctly, save the calibration using option 3 otherwise repeat steps from 1 to 4.
6. Test the suction and washing through options **F1** and **F2** to verify the accuracy of the calibration.

The volume of aspiration to be set must be chosen so that the execution of the analysis is not compromised by small changes in the fluid aspirated. The default settings displayed on the screen will work properly and is not recommended by the operator to change these volumes.

11. PRINTER

The configuration of the printer allows the operator to activate it or not both generally and for printing methods and / or results.

PRINTER	
1-PRINTER ENABLE:	ON
2-AUTOPRINT METHOD:	OFF
3-AUTOPRINT RESULT:	ON

The available options permit to:

1. **PRINTER ENABLE:** enable or disable the printer.
2. **AUTOPRINT METHOD:** autoprint of the method parameters at the beginning of each method execution.
3. **AUTOPRINT RESULTS:** autoprint the results after each sample results.

12. STORED RESULTS

The analyser is equipped with a memory able to store the values of the last fifteen tests. Memory management is automatic and will always show the results from most recent to oldest, which will be the first result to be discarded if it exceeds the maximum number of results stored (fifteen). To access the stored results from the main menu select **3**:

SELECT RESULT

1	ID0070	21/01/10
2	ID0069	21/01/10
3	ID0068	20/01/10
4	ID0067	19/01/10
5	ID0066	17/01/10

F1-SEARCH ID:

NEXT →

Pressing the **F1** key you can directly indicate the ID associated with the progressive results to display.

ACETIC ACID
ID 070 SAMPLE: _____
RESULT:
1.728 ug/dL
FACTOR: 1.000

ESC-EXIT

To delete a result press two times the **Del** key during the result display.

13. QUALITY CONTROL

For each method is possible to program until three quality control levels (Control1, Control2, Control3). Each statistic can store until 31 points for each control level.

13.1 QUALITY CONTROL EXECUTION

Quality control is performed by selection of option **1** on the acceptance screen before each sample execution.

METHOD	
TYPE: END POINT	
1-Q.C. ?:	Control 1
2-ID:	72
3-SET SAMPLE NAME	

F1-CONTINUE	ESC-ABORT

Pressing **1** is possible to select the QC level (Control1, Control2, Control3, None).

SAVE RESULT IN Q.C. STATS?	
F1-YES	F2-NO

Once obtained the result, pressing **F1** key is possible to save the same on the memory.


13.2 Q.C. STATISTIC DATA

To see the statistic of a method controls, is necessary to enter on the method parameters by the function 2-MODIFY METHOD from main menu (see Chapter 5).

EDIT METHOD	
1- Q.C.REF. :	Control 1
2-Min:	0.000
3-Max:	1.000
4-VIEW Q.C. STATS	
5-RESET Q.C. STATS	
← PREV	FIRST →


When the instrument contain more than 2 Q.C. results is possible to view or reset the statistical data by the options **4** and **5**.

Pressing **4** button the statistic graphic appear on the screen:

A solid line represents the mean value, a thick dashed line represents the +/- 2s values, and a dotted line to middle part represents the +/- 3s values.

Pressing the **Enter** key is possible to verify the values for each point of the series by moving with the **←** and **→** key. The screen shows any alerts (warns).

Warns:

- Range: the point is out of the programmed range of the control.
- 2s: the point is above the average for a value greater than 2s.
- 3s: the point is above the average for a value greater than 3s.

14. PERFORMANCES OF REAGENTS FOR WINEMATIC

CODE	REAGENT	LINEARITY (Note 2)	RANGE (Note 2)	LINEARITY (Note 3)		CV% (Note 1)
EN001L	ACETIC ACID	1,2	0,05 - 1,2	24	g/L	0,5-1,5
EN003	CITRIC ACID	0,8	0,05 - 0,8	16	g/L	0,5-1,5
EN004	D-GLUCONIC ACID	2	0,05 - 2	40	g/L	0,5-1,5
EN005	L-LACTIC ACID	0,6	0,05 - 0,6	12	g/L	0,5-1,5
EN006	D-LACTIC ACID	0,6	0,05 - 0,6	12	g/L	0,5-1,5
EN007	L-MALIC ACID	1,2	0,05 - 1,2	24	g/L	0,5-1,5
EN008	AMMONIACAL NITROGEN	60	5 - 60	1200	mg/L	0,5-1,5
EN009	ACETALDEHYDE	120	10 - 120	2400	mg/L	0,5-1,5
EN010	ALPHA-AMMINIC NITROGEN	180	5 - 180	3600	mg/L	0,5-1,5
EN011	ANTHOCYANS	800	5 - 800	16000	mg/L	0,5-1,5
EN012	CHLORIDE	1,5	0,8 - 1,5	30	mg/L	0,5-1,5
EN013	COPPER (no red wine)	4	0,2 - 4	80	mg/L	0,5-1,5
EN014	GLUCOSE-FRUCTOSE	2,5	0,05 - 2,5	50	g/L	0,5-1,5
EN140	GLUCOSE	4	0,05 - 4	80	g/L	0,5-1,5
EN015	GLYCEROL	13	0,5 - 13	260	g/L	0,5-1,5
EN016	IRON	40	0,5 - 40	800	mg/L	0,5-1,5
EN0180	TOTAL POLYPHENOLS	4	0,1 - 4	80	g/L	0,5-1,5
EN019	CATECHIN	60	0,5 - 60	1200	mg/L	0,5-1,5
EN020	CALCIUM	160	2 - 160	3200	mg/L	0,5-1,5
EN021W	TARTARIC ACID	5	1 - 5	100	g/L	0,5-1,5
EN022	SO2 FREE	150	1,5 - 150	3000	mg/L	0,5-1,5
EN023	SO2 TOTAL (no red wine)	200	10 - 200	4000	mg/L	0,5-1,5
EN025	MAGNESIUM	40	1,5 - 40	800	mg/L	0,5-1,5
EN026	POTASSIUM (semi-quantitative)	3000	125 - 3000	60000	mg/L	0,5-1,5
EN055	UREA	50	125 - 3000	1000	mg/L	0,5-1,5

Note 1: calculation made by taking an average of 20 samples

Note 2: linearity calculated without predilution of the sample.

Note 3: linearity calculated with a 1:20 predilution of the sample.

info@gibertini.com
Tel. (++39) 023541434

ISO 9001:2015

LAT N° 094

GIBERTINI
ELETTRONICA s.r.l.

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (++39) 02 3541434
Fax (++39) 02 3541438
www.gibertini.com
sales@gibertini.com