PRECISION SCALES

CENT-2 CE



Simple Weighing

Eng / Deu / Fra

RS-232 / USB

CE Approved for legal metrology

Percentage

GIBERTINI

Via Bellini, 37
20026 Novate Milanese
(Milano) - italy

Tel. (+39) 023541434
info@gibertini.com
sales@gibertini.com

Features

Vacuum fluorescent display (VFD)
Waterproof and acid-resistant membrane keyboard, easy and intuitive to use, with TARE, ON/OFF, PRINT, and MODE keys
Indicator of stable weight
Functions available from the menu: display in grams (g), pounds (lb), ounces (oz), carats (ct), pieces (pcs), and percentage (%).
Access for weighing from below the scale platform.
Selectable “fast/slow” response time
Language : English / Deutsch / French

Technical Specification

Serial Output	RS232 I/O configurable - USB
Temperature	10° - 30°C (optimal 18° - 28°)
Power Tension	90 – 264 VAC, 47 - 60 Hz
Power Absorption	1,5A (18W)
Dimensions (L x D x H)	228 x 381 x 90 mm
Net Weight	5,6 Kg

Model	Capacity (g)	Resolution (g)	Verification Division (e)	Linearity (g)	Repeatability (g)	Response Time	Pan mm	Calibration
LS 6200C CE	6200	0.01	0,1	± 0,015	± 0,01	2/3 sec.	200x200	Internal

Version with CE marking (legal metrology) in accordance with EC Directive 2014/31/EU

Subject to technical changes without notice.
Instruments intended for scientific purpose and/or for internal control only. They must not be used in environments with danger of explosions. Use of our
scales is forbidden in the cases provided by art. 1, point 2, letter a) of directive 2014/31/UE except for those models approved with CE mark (legal
metrology). Electronic balances are sensible to the variation of the acceleration of gravity. They must be calibrated on the working site (2014/31/UE).
All our instruments comply with directive 2014/30/UE (electromagnetic compatibility)

For further information visit **www.gibertini.com**